



TOPSPIN

Installation Guide

for Red Hat

Enterprise Linux WS 3

Copyright (C) 2004 by Bruker BioSpin GmbH

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form, or by any means without the prior consent of the publisher.

Part-No: H9556SL1/02/01

revision state: 12/17/2004

Product names used are trademarks or registered trademarks of their respective holders.

Bruker software support is available via phone, fax, e-mail, Internet, or ISDN.
Please contact your local office, or directly:

Address: Bruker BioSpin GmbH
Service & Support Department
Silberstreifen
76287 Rheinstetten
Germany

Phone: +49 (721) 5161 455
Fax: +49 (721) 5161 943
E-mail: nmr-software-support@bruker-biospin.de
FTP: <ftp.bruker.de> / <ftp.bruker.com>
WWW: www.bruker-biospin.de / www.bruker-biospin.com
ISDN: on request

Contents

Chapter 1	Introduction	9
1.1	About this manual	9
1.2	Conventions	10
1.3	Hardware and software requirements for TopSpin	11
1.3.1	Operating system requirements	11
1.3.2	Special PC Hardware requirements	11
1.3.3	Spectrometer requirements	12
1.3.3.1	Avance systems with AQS rack	12
1.3.3.2	Avance systems with AQX/AQR rack	14
1.3.4	Software requirements	17
1.4	Important Red Hat Linux Enterprise WS 3 configurations	18
1.4.1	KDE - desktop manager; recommended by Bruker BioSpin	18
Chapter 2	Express Installation and Configuration of TOPSPIN	19
2.1	About this chapter	19
2.2	Installing TOPSPIN	19
2.2.1	Installation in parallel	20
2.2.2	What do you want to do with TopSpin?	24
2.3	Configuring TOPSPIN	24
2.3.1	If you want to use a default configuration on a datastation	25
2.3.2	If your PC controls a spectrometer	25
2.3.3	If you want to configure a processing-PC like your spectrometer-PC:	28

Part I

TopSpin Installation and Configuration

Chapter 3	General Installation Information	31
3.1	Installation order	31
Chapter 4	Installing TOPSPIN	33
4.1	NMR SUPERUSER and NMR ADMINISTRATION PASSWORD	33
4.1.1	How to change the NMR ADMINISTRATION PASSWORD	34
4.1.2	How to define a new NMR ADMINISTRATION PASSWORD	35
4.2	Packages on the TOPSPIN DVD	35
4.3	Installing TOPSPIN	36
4.3.1	Installation in parallel	36
4.3.2	What do you want to do with TopSpin?	44
4.4	Perl and GCC	44
4.5	Installing TopSpin patches	45
4.6	The usage of the NMR-GUIDE	46

4.6.1	Requirements for NMR Guide 4.0	46
4.6.2	Java Requirements for NMR-GUIDE 4.0	47
4.6.3	Working with the NMR-GUIDE	47
4.6.3.1	Start up of the NMR-GUIDE	47
4.6.3.2	Functionality test of NMR-GUIDE	48
4.6.3.3	Using the NMR-GUIDE search engine	48
Chapter 5	The TOPSPIN license	51
5.1	Questions and answers about the TOPSPIN license	51
5.2	Ordering and installing the TOPSPIN license	55
5.2.1	Ordering a license	55
5.2.2	How do I determine the correct hostid	56
5.2.3	Installing a Demo license	56
5.2.4	Installing a full license on a license server	56
5.2.5	Installing a full license on a license client	57
5.2.6	Modifying the license file	57
Chapter 6	Configuration of TOPSPIN	59
6.1	Take over the configuration of XWIN-NMR/TOPSPIN into TOPSPIN 1.3	59
6.1.1	You have a XWIN-NMR configuration that you like to use for TopSpin (on the same PC)	60
6.1.2	You have a XWIN-NMR configuration that you like to use for TOPSPIN (on another PC)	61
6.2	Configuring TOPSPIN	61
6.2.1	If your PC controls a spectrometer	62
6.2.2	If you want to configure a processing-PC like your spectrometer-PC:	71
6.2.3	If you want to use a default configuration on a datastation	75
6.2.4	If you do not know the configuration of your spectrometer	75
6.2.5	If you added a new hardware component to your spectrometer	77
6.3	Firmware update of the spectrometer hardware	77
6.4	Configuring the BSMS-Keybaord	79
6.5	A log of the entire system configuration	79
6.6	Backup of VIFs (=Very Important Files) with nmr_save	80
Chapter 7	Plotting and Printing with TopSpin	85
7.1	Printing / plotting data	85
7.1.1	Plotting and printing	85
7.1.2	Print/Plot from the Menu	85
7.1.3	Plot data from the Processing guide	87
7.1.4	Plot data with the TOPSPIN PLOT EDITOR	87
7.1.5	Store (Export) a Data Window as Graphics File	87
7.1.6	TOPSPIN PLOT EDITOR	88
Chapter 8	User Management	89
8.1	Adding a new user account	89
8.1.1	Change group membership of an existing user	91

Part II

RED-HAT Enterprise Linux WS 3 Installation-Tips and Configuration

Chapter 9	RED HAT Enterprise Linux WS 3	95
9.1	Hints for installation of Red Hat	95
9.2	Scratch installation of Desktop computer	95
9.2.1	Necessary additional packages	98
9.2.1.1	Additional packages for spectrometer control - bootparamd	98
9.2.1.2	Additional packages for documentation view - Acrobat Reader	99
9.3	Installation of additional Red Hat packages on already installed systems	100
9.3.0.3	uucp, nfs-utils and redhat-config-nfs	101
9.3.1	Additional packages for the compilation of AU programs	102
9.3.1.1	glibc-kernheaders	102
9.3.1.2	glibc-devel	103
9.3.2	Additional package for sending data by email	104
9.3.2.1	Mutt	104
9.4	Check if all necessary components are installed and active	104
9.5	Dual Boot computer	105
Chapter 10	Network configuration	107
10.1	Configuring the laboratory network	107
10.2	Some useful commands	108
10.2.1	Configuring the laboratory network by DHCP	108
10.2.2	Configuring the laboratory network by using a fixed TCP/IP address	111
10.2.3	Spectrometer PC that should not be connected to the laboratory network	115
10.3	Configuring the spectrometer network	118
10.3.1	Configuration of the spectrometer network	118
10.3.2	Changing the hostname	121
10.4	Email client	121
10.5	Network diagnostic commands	122
10.6	Accessing remote files and directories	122
10.6.1	Sharing directories from a Linux host using Samba	122
10.6.2	Mounting Windows 2000/XP partitions on a Linux host using Samba	123
10.6.3	Mounting Windows 2000/XP partitions from a Linux host using NFS	125
10.6.4	FTP and Telnet	125
Chapter 11	Installing a Printer	127
11.1	Introduction	127
11.2	Print system , CUPS'	128
11.2.1	Installation of a Local printer	128
11.2.2	Installation of a Unix network printer	133
11.2.3	Installation of a Windows network printer	138
11.2.4	Installation of a Jet Direct printer	143

Part III

Tips, Tricks and Troubleshooting

Chapter 12	Troubleshooting TopSpin	151
12.1	TopSpin does not start	151
12.2	TopSpin troubleshooting	152
12.2.1	Troubleshooting while no acquisition is running	152
12.2.2	Troubleshooting while an acquisition is running	153
12.2.3	Send debug information to Bruker BioSpin	155
12.3	Communication problems between PC and CCU	155
12.3.1	Communication problems but telnet spect works	155
12.3.2	Communication problems and telnet spect does not work	157
12.3.2.1	Software related problems	157
12.3.2.2	Hardware related problems	162
12.3.3	Files which are needed for the spectrometer network	162
12.3.3.1	File entries made for the spectrometer network	162
12.3.3.2	Files created for the spectrometer network	164
12.3.4	View the output of the CCU console	165
12.3.4.1	View/change the prom parameters	166
12.3.4.2	Viewing the boot process	166
12.3.4.3	Running CCU diagnostic tests	170
12.4	General steps that you can try to locate or solve problems	170
12.4.1	Check the WWW for FAQ's, known bugs or known problems	170
12.4.2	RED HAT help	171
12.4.3	Checking the history of TopSpin	171
12.4.4	Rebooting the PC	171
12.4.5	Rebooting the PC and the CCU with 'reviveccu'	172
12.5	Various problems in Linux	172
12.5.1	The entire screen is frozen	172
12.5.2	Programs like shmm, reviveccu and touser are not found	173
12.5.3	Xserver does not start after reboot / not enough disk space available	173
12.6	Trouble under Linux	174
12.6.1	You cannot log in as normal user	174
12.7	Useful configuration hints	174
12.8	Spectrometer hardware tests	176
12.8.1	Spectrometer internal boards	176
Chapter 13	Important Linux features	177
13.1	The Boot Floppy	177
13.1.1	Creating a boot floppy	177
13.2	Important Commands	178
13.2.1	rpm	178
13.2.2	cd	178
13.2.3	ls	179
13.2.4	chmod	179
13.2.5	chown / chgrp	180
13.2.6	grep	181
13.2.7	find	182
13.2.8	ps	183

13.2.9	kill	184
13.2.10	compress, gzip	184
13.2.11	tar	186
13.3	Important Files	187
13.3.1	/etc/shadow	187
13.3.2	/etc/passwd	188
13.3.3	/etc/group	190
13.3.4	.bashrc	191
13.3.5	umask	192
13.4	Linux shells	193
13.4.1	How to use the bash	194
13.5	Useful tips for Linux	196
13.5.1	Configure the screen resolution with Xconfigurator	196
13.6	How to avoid much typing in Linux	197
13.6.1	General	197
13.6.2	Linux shell scripts	198
13.6.3	cron - crontab	200
Chapter 14	File and directory handling	201
14.1	Permissions and Type of Access	201
14.1.1	Permissions of a file	201
14.1.2	Permissions of a directory	202
14.1.3	Types of Access	202
14.1.4	Example	202
Chapter 15	Useful tools for Red Hat Enterprise Linux WS 3	205
15.1	How to use a memory stick under Linux	205
15.1.1	PC with SCSI drives	205
15.1.2	PC with IDE drives	208
15.2	Firewall configuration	209
15.2.1	Firewall configuration	210
15.3	Tools for internal communication	212
15.3.1	Installation description of an additional package: xinetd	213
15.3.2	Check if all necessary components are installed	213
15.4	Setting up permissions for rsh, rlogin and rcp	214
15.5	Creating a screenshot of the desktop	215
15.6	Creating an icon on the desktop	215
15.6.1	Create an icon for opening a pdf file with Acrobat Reader	215
15.7	Red Hat Network: Hotfixes to keep the system secure	217
15.7.1	Activate your Linux license at Red Hat Network	217
15.7.2	Register the PC to the Red Hat Network	219
15.7.3	Activate the Red Hat Alert Notification Tool	224
15.8	The Emergency Recovery CD / Image CD	226
15.8.1	Creating a boot floppy	227
15.8.2	Creating the Rescue Disk	227
15.8.3	Creating the Emergency Recovery CD	227
15.8.4	Using the Emergency Recovery CD	228
15.9	Using the CD-Writer	229
15.10	Add a new harddisk	231

Chapter 16	Appendix	233
16.1	Some notes for the helium handling	233
16.1.1	Automatic notification for helium refill	233
16.1.1.1	Set minimum limit of helium level	233
16.1.1.2	Create shell script 'heliumwarn' for email notification	234
16.2	Script for saving important files	236
16.3	Where Do I Find What?	237
Chapter 17	Books and Web addresses on Linux	239
Chapter 18	Bruker addresses	241

Chapter 1

Introduction

1.1 About this manual

This manual describes the installation of the TopSpin 1.3 under RED HAT ENTERPRISE LINUX WS 3.

It is available:

- as a hard copy 'Installation Guide for TOPSPIN: RED HAT ENTERPRISE LINUX WS 3'
- as pdf file on the toplevel of the TOPSPIN DVD
- as TOPSPIN Online Help: click **Help** → **Installation Guide**
- as the up to date version on the Bruker web-server:

<http://www.bruker-biospin.de/NMR/nmr-softw/passwd/docu/index.html>

- Note that this manual does not contain the TOPSPIN Release letter. That is a separate document which is not available as a hard copy but delivered as pdf and text files on the TOPSPIN DVD and also on the Bruker web-server (you can find it at the same URL as the Installation Guide)

<http://www.bruker-biospin.de/NMR/nmr-softw/passwd/docu/index.html>

This manual is subdivided into an introduction and three main parts. It depends on your local demands which of these are necessary for you at the moment.

The first part shows how to install and configure TOPSPIN and some essential additional software packages. It is recommended to follow the instructions step by step.

Part two describes the scratch installation of RED HAT ENTERPRISE LINUX WS 3 and some configuration steps like setup for network and printers. If you want to install Bruker BioSpin TOPSPIN on a PC which is already in use, you can skip this part.

In the last part, you can find some additional information for working with RED HAT ENTERPRISE LINUX WS 3 and TOPSPIN and also some tips for troubleshooting.

Note that all text-messages and graphics shown in this book are taken from the current version of the respective software products. Small differences to other versions are possible but generally speaking they should be very similar.

If you have any comments or suggestions, or if you find any errors in this manual, please do not hesitate to contact us at:

`nmr-software-support@bruker.de`

1.2 Conventions

The following conventions will be used throughout this manual:

Bold-Courier: Linux commands to be entered from the keyboard

Bold-Times: Linux commands/applications to be clicked with the mouse

Bold-Courier-Italics: NMR Suite commands to be entered from the keyboard

Bold-Times-Italic: NMR Suite commands to be clicked with the mouse

Courier: a file name or directory name

Courier Small Italics: responds of the system

Courier Small: the contents of a file

Times italics: any name which is not a file name e.g. host names, user names

etc.

< >: place holder

1.3 Hardware and software requirements for TopSpin

1.3.1 Operating system requirements

The following operating systems are tested and supported for TOPSPIN 1.3.

Linux

Red Hat Enterprise Linux WS 3

Windows

Windows 2000 / Windows XP

Any other operating systems might also work but because TOPSPIN 1.3 was not tested there, usage of not supported operating system is at one's own risk.

1.3.2 Special PC Hardware requirements

To run TOPSPIN 1.3 we recommend the following computer hardware:

- PC¹ with CPU Clock \geq 1 GHz
- Memory \geq 512 MB
- Video (graphics)-card memory: 64 MB ²
- 2 Ethernet cards (the second card is only required if your PC controls a spectrometer)
- DVD device³
- Mouse with 3-button possibilities

If your hardware does not meet these requirements, TOPSPIN might still run but with a much lower performance.

Furthermore, we recommend setting the display resolution to 1280*1024 or higher.

Installation was tested on a HP XW4100 PC running under Red Hat Enterprise Linux WS 3 and TopSpin 1.3.

-
1. TOPSPIN only runs on x86-based systems. It does not support 64bit architecture
 2. It is recommend not to use shared memory graphics, because in some cases memory problems could be observed.
 3. TopSpin is delivered on DVD's with "-R" format

1.3.3 Spectrometer requirements

TOPSPIN 1.3 is released for datastations and Avance-AQS, Avance-AQX/AQR.

For upgrade possibilities of your spectrometer with TOPSPIN 1.3, please check with your local Bruker office.

Hardware components that are TOPSPIN version or operating system dependent are listed below. Components that are not listed here do not need an upgrade for TOPSPIN 1.3.

Please note that AQS and AQX/AQR components cannot be interchanged.

If your Avance spectrometer does not meet the above specifications, please get in touch with your local Bruker representative for your upgrade options.

1.3.3.1 Avance systems with AQS rack

- CCU:

LINUX, independent of your TOPSPIN version	
CCU	required minimum EC level
CCU 10 (part # H9503, EC level 00 to 19)	06
CCU 10 (part # H9503, EC level > 19)	22
CCU 11 ^a (part # H9503V1)	00 (recommended 01)

Windows, independent of your TOPSPIN version	
CCU	required minimum EC level
CCU 10 (part # H9503)	any
CCU 11 ^a (part # H9503V1)	00 (recommended 01)

a. Note that CCU 11 also exists for AQX/AQR spectrometers.

- **FCU 3:**

FCU 3	required minimum EC Level
part # H5822 (64k memory)	04
part # H9598 (256k memory)	02

- **RCU:** (Note: RCUs of the AQX/AQR rack can be used also in a AQS rack)

RCU	EC level
RCU 1/2 (part # Z003678)	any
RCU 1/3 (part # Z022488)	any

- **Multiple Receiver RCU:**

Multiple Receiver RCU	EC Level
RCU 1/4 (part # Z052488)	any
RCU 2/4 (part # Z062488)	any
RCU 3/4 (part # Z072488)	any
RCU 4/4 (part # Z082488)	any
RCU 1/5 (part # Z003206)	any
RCU 2/5 (part # Z003207)	any
RCU 3/5 (part # Z003208)	any
RCU 4/5 (part # Z003209)	any

- **TCU 3:**

TCU 3	required minimum EC Level
part # H5813, EC level 00 to 19	09
part # H5813, EC level 20 to 29	28
part # H5813V1	any
part # H5813F2	any
part # H5813F3	any

1.3.3.2 Avance systems with AQX/AQR rack

RCU	required minimum EC level
RCU 1/1 (part # Z002488) 64 kB SRAM & 1MB DRAM	06
RCU 1/2 (part # Z012488) 256 kB SRAM & 2MB DRAM	06
RCU 1/3 (part # Z022488) 1 MB SRAM & 2MB DRAM	06

WINDOWS	
CCU	required minimum EC level
CCU 04 (part # H2552)	04
CCU 05 (part # H2570, EC 00 to 19)	05
CCU 06 (part # H2570, EC 20 to 29)	24
CCU 08 (part # H2570, EC 30 to 39)	35
CCU 08 (part # H2570, EC 40 to 49)	43
CCU 09 (part # H2570)	any
CCU 11 ^a (part # H2570V1)	any

a. Note that CCU 10 does not exist for AQX/AQR spectrometers. Further note that CCU 11 also exists for AQS spectrometers.

LINUX	
CCU	required minimum EC level
CCU 04 (part # H2552)	05
CCU 05 (part # H2570, EC 00 to 19)	06
CCU 06 (part # H2570, EC 20 to 29)	26
CCU 08 (part # H2570, EC 30 to 39)	38
CCU 08 (part # H2570, EC 40 to 49)	46
CCU 09 (part # H2570, EC 50 to 59)	55
CCU 09 (part # H2570, EC 60 to 69)	62
CCU 11 (part # H2570V1)	any

TCU	required minimum EC level
TCU 1 main board (part # H5811 or 5812)	03
TCU 0 main board (part # H2558)	07
TCU extention board (part # H2562)	03
TCU extention board (part # H2562) in DMX/DRX with > 5 channels	20

Note: 'normal' RCUs and multiple receiver RCUs should not be mixed!

Multiple Receiver RCU	EC level
RCU 1/4 (part # Z052488)	any
RCU 2/4 (part # Z062488)	any
RCU 3/4 (part # Z072488)	any
RCU 4/4 (part # Z082488)	any
RCU 1/5 (part # Z003206)	any
RCU 2/5 (part # Z003207)	any
RCU 3/5 (part # Z003208)	any
RCU 4/5 (part # Z003209)	any

FCU0	Required minimum EC level
64k memory: part # H2556 64k memory: part # H2554	
FCU0: for use without cortab	04
FCU0: for DSX without use of cortab	05
FCU0: for use with cortab with EC level between 00 and 19	07
FCU0: for use with cortab with EC level ≥ 20	21

Tomo FCU	required minimum EC level
part # T5565 (64k memory) for use without cortab	00
part # T5565 (64k memory) for use with cortab	02
part # T6765 (256k memory) for use with or without cortab EC level 00 to 19	00
part # T6765 (256k memory) for use with or without cortab with EC level 20 to 29	21

1.3.4 Software requirements

! Because of the variety of Linux systems and their possible problems with some PC-hardware components, Bruker only supports RED HAT ENTERPRISE LINUX WS 3 and the hardware components described in section 1.3.2. This software/hardware combination was successfully tested with the TopSpin 1.3. Any other Linux systems and PC hardware components might work, but are not tested and supported by Bruker.

Requirements for TOPSPIN:

- Red Hat Enterprise Linux WS 3
- Adobe Acrobat Reader (*WWW*)

Requirements for NMR-GUIDE:

- Web Browser: Mozilla (*CD* or *WWW*)

1.4 Important Red Hat Linux Enterprise WS 3 configurations

1.4.1 KDE - desktop manager: recommended by Bruker BioSpin



All information described in this manual is based on the desktop manager KDE. If you use GNOME, there will be differences in names, directories and files.

In the login window, type in:

- your login name
- your password
- choose as *Session* **KDE**¹
- press **Enter**

1. If you choose KDE as login session the first time, you will be asked if KDE should become your default session type. If you agree KDE will be used automatically in the future.

Chapter 2

Express Installation and Configuration of TOPSPIN

2.1 About this chapter

This chapter gives a brief description of how to install Bruker BioSpin's TOPSPIN.

This chapter is directed to users who would like to install TOPSPIN 1.3 on a PC, on which XWIN-NMR or a version of TOPSPIN is already installed and configured correctly. Therefore this chapter provides an overview and provides not information on all aspects of configuring TOPSPIN. Further details can be found in the following chapters.

For more details on the installation see chapter 4

For more details on the *TOPSPIN* commands like *cf*, *expinstall* etc. please refer to the description of these commands in the Acquisition Reference manual.

2.2 Installing TOPSPIN

The main decision you have to make during installation is where do you want to install the new TopSpin version. If you have previous NMR SUITE or TOPSPIN

installation on this PC available, you can either:

- install the new version in parallel to the previous one
- or
- install the new version in the same directory as the previous one

The installation in parallel has the advantage, that the previous version remains unchanged. After a parallel installation you can start working with the new version and if you want to go back to the previous version for whatever reasons you can just start and reuse it.

The installation of TOPSPIN 1.3 offers an automatic import of your individual configuration of a previous version, so that there is no advantage to install the new version in the same directory as the old one.

2.2.1 Installation in parallel

1. Close ALL WINDOWS on the desktop.
2. Insert the DVD 'TOPSPIN 1.3'
3. Start of the installation
 - If the automatic start of the DVD is enabled, the shell described in the next step appears automatically, go to step 4
 - If the shell of the next step does not appear automatically you have click on the DVD icon on your desktop and click on the file `install` ¹
4. If you are not logged in as root, a Linux shell will appear that asks you for the root password:

```
Please enter root password if prompted for it:  
Password:
```

 - Enter root password and press enter
5. The window *TOPSPIN Installation Manager* will appear
6. The window *Welcome!* will appear
 - Click **Next**

1. Please note that there is also a file called 'install.cmd' on the DVD. This file is the installation script for windows, under Linux you have to click on the file 'install'

7. The window *TOPSPIN Release Letter* will appear:

Please read the NMR Release letter! You can find therein all information about news and bug fixes of the current version.

 - When you have finished reading the Release Letter click **Next**
8. The window *Please select a target directory for TOPSPIN* will appear:

Enter an installation path which does not exist so far, e. g. */opt/topspin1.3*

 - Click **Next**
9. If the installation directory does not exist so far, you will be asked if the installation process should create it
 - Click **Yes**
10. A window *Please select the setup type* will appear:
 - a) Select the setup type. Do one of the following:
 - Select *Data Processing only* if you want to use TOPSPIN for data processing.
 - Select *Data processing, acquisition and automation* if you want to use TOPSPIN for spectrometer control. This will e. g. install **Diskless** automatically.
 - *Customized setup* (allows components selection) if you want to install selected programs only.
 - b) Click **Next**
- 10.1. If *Customized* was selected in step 10 ...
 - a) ... a window *Please select products to install* will appear:
 - Select the packages which you want to install (note that DISKLESS is only required if your PC controls a spectrometer or you want to configure it like this (see chapter 6.2.2)).
 - Click **Next**
11. The window *Select a spectrometer configuration to import* will appear. You can either import no spectrometer configuration, or browse for a spectrometer configuration or select one of the configurations you get offered within this window.¹
12. If a previous topspin installation is available in the global environment variable Path, you will be asked if you want to define the current installation path

instead. It is recommended to set the current installation into the PATH variable, so enable:¹

Yes, add the new TOPSPIN installation directory to the global PATH

- Click **Next**

13. A window appear that asks for the installation path of FlexLM. It is strongly recommended to accept the default path²!

- Click **Next**

14. A window appear that displays the program selection of step 10a

- Click **Next**

15. A window 'Set NMR Super User' will appear:

a) Do one of the following:

- Accept the default NMR SUPERUSER (nmrsu) or
- Enter a special user for this purpose, e.g. named *nmr-admin* or
- Enter any other user

If the user does not exists, the installation will create it. For further details see chapter 4.1 or the TOPSPIN Users Guide.

b) Click **Next**

16. A window 'Password Input' will appear and asks for the NMR ADMINISTRATION PASSWORD³ (for more details about this feature see chapter 4.1)

- Type in the NMR ADMINISTRATION PASSWORD
- Re-enter the NMR ADMINISTRATION PASSWORD
- Click **Next**

17. If *Data Processing, acquisition and automation* or *Customized* with the Disk-

-
- The installation process will offer all configurations which are available in the TOPSPIN / XWIN-NMR installation that are defined in the /etc/profile.d/xwinnmr(csh/.sh) resp. top-spin(csh/.sh) directory.
 - In case more than one TopSpin version is installed in parallel it is your decision which one should be defined here
 - It is strongly recommended to accept the default path for the FlexLM installation: /usr/local/flexlm ! If you choose a user defined path, it may cause some problems with other applications using the FlexLM
 - Certain commands like *cf* and *expinstall* ask for the NMR ADMINISTRATION PASSWORD.

less package was selected in step 8 a window will appear to inform that the spectrometer has to be shut down.

- a) If the PC is not yet connected to a spectrometer or if the CCU is not booted¹ you can promptly
 - Click **Next**
- b) If your PC is currently connected to a spectrometer and the CCU is booted:
 - a) Open a Linux shell
 - b) enter the command:
telnet spect
 - c) login as *root*
 - d) enter the command:
init 5
→ this will shut down the CCU
 - e) Close the Linux shell
 - f) Click **Next** in the window that asked you about shutdown the CCU

Note: The CCU is now either down or switched off and should be turned on again at the end of this procedure.

→ The installation of TOPSPIN will start now.

18. After installation is finished a window appears that informs about all installed programs
 - Click **Seen**
19. Depending if you have installed DISKLESS one of the following windows will appear:
 - If you have installed DISKLESS a window will appear that asks you to reboot the PC. Make sure that you can reboot the PC now (check for any unsaved windows (e. g. an editor window) and check if no other users are logged in into this PC) then click on **Yes**
 - If you have not installed DISKLESS a window will appear that informs you to log out and log back in before you start TOPSPIN, click **Seen**

1. Note that during a first time installation of the TOPSPIN the CCU can not be booted.

20. The window wherein you entered the root password informs you that you have to press *Enter* to close this window
 - Press **Enter**
21. If you have started the installation manually by clicking the file 'install' in the Konqueror window, you have to close this Konqueror window.
22. To eject the DVD click right on the DVD icon on the desktop and choose **eject**

2.2.2 What do you want to do with TOPSPIN?

- it should control a spectrometer or you want to use a special spectrometer configuration
 - you have to execute *cf* and *expinstall*
- it should be used for data processing only?
 - you have to execute *expinstall*

For more details please refer to the next chapter

2.3 Configuring TOPSPIN

The idea of this chapter is to give you some guidelines to configure TOPSPIN for using it on a workstation for data processing, or for operating a spectrometer. We do not discuss the commands in detail here. For a full description refer to the TOPSPIN manuals.

Start TOPSPIN:

- Click the icon TOPSPIN on the desktop¹, or
- Click <redhat> → **Other** → **TOPSPIN<vers.>**

After the first startup of TOPSPIN you have the possibility to configure the software depending upon the hardware of your spectrometer. If you do this, you are probably in one of the following five situations:

1. You want to configure a processing-only-PC with a default configuration on a datastation → go to chapter 2.3.1 (for details see chapter 6.2.3)
2. You have installed the new software on the same PC as a previous XWIN-NMR/TOPSPIN installation, so the old spectrometer configuration is available on

1. How to create a icon? See chapter 15.6.1

the same PC → go to chapter 2.3.2 (for details see chapter 6.2.1)

3. You want to configure a processing-only-PC like your spectrometer-PC, so you have the spectrometer configuration on another computer → go to chapter 2.3.3 (for details see chapter 6.2.2)
4. You have to configure the software for controlling a spectrometer on a new hard disk without a backup of the spectrometer configuration, so you have no configuration files → go to chapter 6.2.4
5. You have to configure the software for controlling a spectrometer with a new spectrometer hardware component → go to chapter 6.2.5

2.3.1 If you want to use a default configuration on a datastation

A default configuration can be used on a PC that is used for processing only. It is automatically performed (no **cf** required) during the installation of TOPSPIN ...

- in a new TOPSPIN installation directory
- if no TOPSPIN/XWIN-NMR was installed on this computer so far

The default configuration name is *Bruker_default_av500* and corresponds to an Avance 500 MHz spectrometer.

For data processing only you have to execute **expinstall** and select *Installation for Datastation (Default)*.

For further details please refer to the respective chapters (see chapter 6.2.3).

2.3.2 If your PC controls a spectrometer:

1. Be sure that the spectrometer is booted (see chapter 9.3.0.3 and 12.3.4.2)
2. Open a dataset
 - Type **new** and create a new dataset.
 - You can also select a dataset (e.g. 'exam1d_1H') acquired with another TOPSPIN or XWIN-NMR version
3. Enter **cf** (The command **cf** now configures the software according to the hardware of your spectrometer.)
 - a) A window will appear:
Password request. Please enter the NMR ADMINISTRATION PASSWORD

Type in the NMR ADMINISTRATION PASSWORD. Click **OK**.

- b) A window will appear:

Spectrometer configuration: ...

Select your previous configuration and press **Edit** to confirm this configuration.

- c) A window will appear:

Edit configuration. ...

The selected configuration is displayed. Choose whether to configure it as a spectrometer or datastation. Choose the type of your spectrometer e.g. **AVANCE** → **ENTER**. Check or type in the 1H frequency of the magnet e. g. **500.13** → **ENTER**. Use the debug mode only in case of troubleshooting.

- d) A window will appear:

Specify the channel to which external devices are connected. The RS232/485 channel for external devices-Table appears.

Check if the tty-numbers of the hardware components are set correctly.

- e) A window will appear:

Additional configuration. Security check. Enable peak power check (POWCHK). Enable puls power check (PULCHK).

Choose the options if you want to use this security options. *Powchk* and *Pulchk* require special probes and *cortab*. If in doubt contact your local Bruker representative, please.

- f) Only if you have a sample changer!

Should the Sample Changer control the Lift? (yes/no)

The answer depends on how the air flow is connected, normally you will answer: **yes** → **ENTER**

If you have a sample changer, this question appears:

Delay between SX and next command [sec]?

Type in the length of the delay e. g. **10** → **ENTER**

- g) The *nuclei-table* appear

Check if the frequency of the nuclei are set correctly. If in doubt click **Restore** then click **Save**

- h) The results of **cf** appears in a text window. Check this list carefully and compare the result with an list from an earlier configuration (if available) click **Print** and store the paper output with your other spectrometer documentation, then click **OK**. These data are saved in:

```
<TOPSPINHOME>/conf/instr/<spect>/uxnmr.info
```



If during 'cf' an error message occurs that invites you to do a firmware update see chapter 6.3

After finishing **cf** the window 'config' will appear showing a list of configuration steps. You will be asked for the NMR ADMINISTRATION PASSWORD several times:

- Continue with **expinstall**
 - a) Type in the NMR ADMINISTRATION PASSWORD and **ENTER**
 - b) Click **Next** if you either have not modified Bruker BioSpin Parameter files, AU-programs and Pulse-programs or if you have archived them in a different directory. Your individual Parameter files, AU-programs and Pulse-programs etc. that are named different to the original Bruker BioSpin files will not be destroyed during the installation.
 - c) A window appears: Check *Installation for Spectrometer* → Click **Next**
 - d) In the appearing window: Check *the required system* → Click **Next**
 - e) In the appearing dialog box: Select the spectrometer configuration name → Click **Next**
 - f) In the appearing dialog box: Check the items you want to install or, accept the default selection → Click **Next**
 - g) In the appearing dialog box: Select the default printer and plotter and the desired paper format → Click **Next**
 - h) In the appearing window: Enter the desired spectrometer frequency, acquisition mode and pre-scan-delay or, accept the default values → Click **Next**
 - i) In the appearing information box: Check the configuration selection and, if it is correct → Click **Finish**
and wait ...
- Continue with **edsolv**

Check the list of lock solvents → click **Save** if it is correct → Click **Close**
- Continue with **edhead**

Select the current probe → click **Define as current probe** → Click **Exit**

- Continue with **edprosol**
Check the parameters for your observed nuclei, and if necessary change and/or recalculate them → Click **Save** to store to disk → Click **Exit**
- Continue with **edlock**
Check the lock table values (field, lock phase, ...), change them and click **Save** if necessary, otherwise click **Abort**
- Continue with **edscon**
Check the spectrometer constants, change them and click **Save** if necessary, otherwise click **Cancel**

After **cf**, **expinstall**,... have finished insert a sample and make sure that you can lock on the lock solvent. Type **ii** and read with **rpar** a standard Bruker BioSpin parameter file (e. g. '**rpar PROTON all**'), type **getprosol**, **rga**, **wobb** and **zg** and collect a normal ^1H NMR spectrum.

2.3.3 If you want to configure a processing-PC like your spectrometer-PC:

If the PC does not control a spectrometer but is used in connection with a particular spectrometer, e.g. for setting up experiments, processing and or plotting, you should configure it as that spectrometer (see FAQ 210 on www.bruker-biospin.de) (see chapter 6.2.2).

Part I

TopSpin Installation and Configuration

Chapter 3

General Installation Information

3.1 Installation order

A complete new installation on a PC connected to a spectrometer involves the following steps:

1. Installing RED HAT ENTERPRISE LINUX WS 3 (see chapter 9.1)
2. Installing additional packages (see chapter 9.2)

Step 1 to 2 are normally already done by the manufacturer.

3. Configuring the network (see chapter Chapter 10)
4. Installing TOPSPIN (see chapter 4.3)
5. Installing TOPSPIN license (see chapter Chapter 5)
6. Configuring of TOPSPIN (see chapter 6.1)
7. Installing a printer (see chapter Chapter 11)

If your PC has already been used for TOPSPIN/XWIN-NMR, steps 1 and 2 are probably already done. In that case, you can start with step 3.

Chapter 4

Installing TOPSPIN

4.1 NMR SUPERUSER and NMR ADMINISTRATION PASSWORD

TOPSPIN offers a complete new password check for the configuration tools of TOPSPIN. During the TOPSPIN installation, you are prompted to define:

1. The username for the so called NMR SUPERUSER. The NMR SUPERUSER can be *nmr_{su}* (this is the default) or any other user (if the user does not exist so far, it will automatically be created during installation process).
 - After the installation, the NMR SUPERUSER is the owner of all TOPSPIN program files. Log in as this user allows you to remove these files, change file permissions etc. The name of the NMR SUPERUSER will be stored in the readable text file: `<TOPSPINHOME>/conf/nmrsuperuser`
2. the NMR ADMINISTRATION PASSWORD to be used for TOPSPIN configuration commands. This password can be freely chosen and is not connected to any user. TOPSPIN asks for this password by commands like *cf, expinstall* etc. The encrypted NMR ADMINISTRATION PASSWORD is stored in the file:

`<TOPSPINHOME>/conf/nmradminpassword`



Note: The NMR SUPERUSER login password and the NMR ADMINISTRATION PASSWORD have different purposes and are totally independent. Changing one of them does not affect the other.

With the information about the NMR ADMINISTRATION PASSWORD and NMR SUPERUSER TOPSPIN is ready to be configured. The NMR ADMINISTRATION PASSWORD will be asked for every time a configuration type command is entered. If this password is typed in correctly, the command will be executed regardless of which user is logged in the operating system. As a result all newly created configuration files will be owned by the user that was logged in at the time the configuration command was executed. Any other user can modify the configuration during his login session, just by knowing the correct NMR ADMINISTRATION PASSWORD. This is possible because all configuration files are read- and writeable for everybody, but TOPSPIN will not write them without the correct NMR ADMINISTRATION PASSWORD!

TOPSPIN also offers a so-called *Security Mode*. With this feature all configuration files are read- and writeable only for the NMR SUPERUSER.

4.1.1 How to change the NMR ADMINISTRATION PASSWORD

The *installnmr* script can be used to change the NMR ADMINISTRATION PASSWORD and the NMR SUPERUSER. If it is called with the option *'-secure'*, it changes the permissions in a way, that only the NMR SUPERUSER can execute the configuration commands and the resulting files are readable and writable only for him, all other users have read-only permissions.

To change the NMR ADMINISTRATION PASSWORD, you have to:

1. Login as NMR SUPERUSER
2. Open a Linux shell.
3. Go to the TOPSPIN installation directory, in a default installation this is

```
/opt/topspin
```

4. Enter the command:

```
prog/bin/installnmr <TOPSPINHOME> <NMRSUPERUSER>
```

If you installed TOPSPIN in the default directory and you want to use the user 'nmrsu' as NMR SUPERUSER, you have to type in

```
prog/bin/installnmr /opt/topspin nmrsu
```

5. Enter the old password and new password as requested.

4.1.2 How to define a new NMR ADMINISTRATION PASSWORD

To define a new NMR ADMINISTRATION PASSWORD, for example because you have forgotten it, you have to:

1. Open a Linux shell
2. Login as root with the command
`su -`
3. Delete the file
`<TOPSPIN_HOME>/conf/nmradminpassword`
4. Enter the command
`prog/bin/installnmr <TOPSPINHOME> <NMRSUPERUSER>`
5. Enter new password as requested.

4.2 Packages on the TOPSPIN DVD

1. TOPSPIN: Acquisition and processing software
2. TOPSPIN PLOT EDITOR: Object oriented WYSIWYG plot editor (automatically installed as part of the TOPSPIN package)
3. ICON-NMR: Icon-driven interface for Routine Spectroscopy and Automation
4. NMR-GUIDE: Web-browser based teaching and training program
5. DISKLESS: Spectrometer CCU operating system
6. NMR-SIM: Program for numerical simulation of NMR experiments
7. AMIX-Viewer: new multiple object viewer (requires a separate license)
8. AURELIA: Analysis of 2D/3D/4D NMR data (requires separate license)
9. GLP: Program for 'Good Laboratory Practice' tests (requires separate license)
10. FLEXLM: License manager required for starting TOPSPIN programs.
11. GNU TOOLS: System tools required by TOPSPIN.
(automatically installed as part of the TOPSPIN package)
12. PERL: language interpreter required for certain TOPSPIN commands.
(automatically installed as part of the TOPSPIN package)

4.3 Installing TOPSPIN

This chapter describes the installation of the TOPSPIN.

The main decision you have to make during installation is where do you want to install the new TopSpin version. If you have a previous NMR SUITE or TOPSPIN installation on this PC available, you can either:

- install the new version in parallel to the previous one
- or
- install the new version in the same directory as the previous one

The installation in parallel has the advantage, that the previous version remains unchanged. After a parallel installation you can start working with the new version and if you want to go back to the previous version for whatever reasons you can just start and reuse it.

The installation of TOPSPIN 1.3 offers an automatic import of your individual configuration of a previous version, so that there is no advantage to install the new version in the same directory as the old one.

4.3.1 Installation in parallel

1. Close ALL WINDOWS on the desktop
2. Insert the DVD 'TOPSPIN 1.3'
3. Start of the installation
 - If the automatic start of the DVD is enabled, the shell described in the next step appears automatically, go to step 4
 - If the shell of the next step does not appear automatically you have click on the DVD icon on your desktop and click on the file `install`¹
4. If you are not logged in as root, a Linux shell will appear that asks you for the root password:



Please enter root password if prompted for it:

Password:

- Enter root password and press enter


1. Please note that there is also a file called 'install.cmd' on the DVD. This file is the installation script for windows, under Linux you have to click on the file 'install'

5. The windows *TOPSPIN Installation Manager* and *Welcome to the Bruker Software Installation Manager* will appear :




- Click **Next**

6. A window *Release Letter TOPSPIN* will appear:



- Please read the NMR Release letter! You can find therein all information about news and bug fixes of the current version.
- When you have finished reading the Release Letter click **Next**


7. A window *Please select a target directory for TOPSPIN* will appear:



Enter an installation path which does not exist so far, e. g. */opt/topspin1.3*


- Click **Next**

8. If the installation directory does not exist you will be asked if the installation process should create it automatically
 - Click **Yes**
9. A window *Please select the setup type* will appear:



- a) Select the setup type. Do one of the following:
 - Select *Data Processing only* if you want to use TOPSPIN for data processing.
 - Select *Data processing, acquisition and automation* if you want to use TOPSPIN for spectrometer control. This will e. g. install **Diskless** automatically.
 - *Customized setup* (allows components selection) if you want to install selected programs only.
 - b) Click **Next**
- 9.1.** If *Customized* was selected in step 9 ...

... a window *Please select products to install* will appear:



- Select the packages which you want to install (note that DISKLESS is only required if your PC controls a spectrometer or you want to configure it like this (see chapter 6.2.2)).
 - **TOPSPIN**: Acquisition and processing software
 - TOPSPIN PLOT EDITOR: Object oriented WYSIWYG plot editor
 - ICON-NMR: Icon-driven interface for Routine Spectroscopy and Automation
 - DISKLESS: Spectrometer CCU operating system
 - NMR-SIM: Program for numerical simulation of NMR experiments
 - NMR-GUIDE: web-browser based teaching and training program
 - NMR-Data: spectroscopic example data
 - AURELIA-AMIX-Viewer: Viewer of 2D/3D/4D NMR data (requires separate license)
 - NMR-GLP: Program for 'Good Laboratory Practice' tests (requires separate license)
- Click **Next**

10. The window *Select a spectrometer configuration to import* will appear. You can either import no spectrometer configuration, or browse for a spectrometer configuration or select one of the configurations you get offered within this win-

dow.¹




The configurations that are shown in this picture are just an example. On your PC this list could look different.

11. If a previous topspin installation is available in the global environment variable Path, you will be asked if you want to define the current installation path instead.² It is recommended² to set the current installation into the PATH variable, so enable:

Yes, add the new TOPSPIN installation directory to the global PATH

- Click **Next**


12. A window appear that asks for the installation path of the FlexLM. It is strongly recommended to accept the default path!³



- Click **Next**


-
1. The installation process will offer all configurations which are available in the TOPSPIN / XWIN-NMR installation that are defined in the /etc/profile.d/xwinnmr(.csh/.sh) resp. topspin(.csh/.sh) directory.
 2. In case more than one TopSpin version is installed in parallel it is your decision which one should be defined here
 3. It is strongly recommended to accept the default path for the FlexLM installation: /usr/local/flexlm ! If you choose a user defined path, it may cause some problems with other applications using the FlexLM

13. A window appear that displays the programs which will be installed



- Click **Next**

14. A window will appear that prompted for a definition of an NMR SUPERUSER



a) Do one of the following:


- Accept the default NMR SUPERUSER (*nmrsu*) or
- Enter a special user for this purpose, e.g. named *nmr-admin* or
- Enter any other user

If the user does not exist, the installation will create it. For further details see chapter 4.1 or the TOPSPIN Users Guide.


b) Click **Next**

15. A window 'Password Input' will appear and asks for the NMR ADMINISTRATION

TION PASSWORD¹ (for more details about this feature see chapter 4.1)



- Type in the NMR ADMINISTRATION PASSWORD
 - Re-enter in the NMR ADMINISTRATION PASSWORD
- c) Click **Next** → The installation of TOPSPIN will start now.
16. If *Data Processing, acquisition and automation* or *Customized* with the Diskless package was selected in step 8 a window will appear to inform that the spectrometer has to be shut down



- a) If the PC is not yet connected to a spectrometer or if the CCU is not booted² you can promptly
- Click **Next**
- b) If your PC is currently connected to a spectrometer and the CCU is booted (see chapter 9.3.0.3 and 12.3.4.2)
- a) Open a Linux shell
- b) enter the command:
- telnet spect**

1. Certain commands like **cf** and **expinstall** ask for the NMR ADMINISTRATION PASSWORD.


2. Note that during a first time installation of the TOPSPIN the CCU can not be booted.

- c) login as *root*
- d) enter the command:
init 5
→ this will shut down the CCU
- e) Close the Linux shell
- f) Click **Next** in the window that asked you about shutdown the CCU

Note: The CCU is now either down or switched off and should be turned on again at the end of this procedure.

→ The installation of TOPSPIN will start now.


17. After installation is finished the window 'Setup is complete' will appear:




- Click **Seen**

18. Depending if you have installed DISKLESS one of the following windows will appear:

- If you have installed DISKLESS a window will appear that asks you to reboot the PC. Make sure that you can reboot the PC now (check for any unsaved windows (e. g. an editor window) and check if no other users are logged in into this PC) then click on **Yes**



- If you have not installed DISKLESS a window will appear that informs you to log out and log back in before you start TOPSPIN, click **Seen**



19. The window wherein you entered the root password informs you that you have to press *Enter* to close this window
 - Press **Enter**
20. If you have started the installation manually by clicking the file 'install' in the Konqueror window, you have to close this Konqueror window.
21. To eject the DVD click right on the DVD icon on the desktop and choose **eject**
22. To start TOPSPIN log out and log back in

4.3.2 What do you want to do with TOPSPIN?

- it should control a spectrometer or you want to use a special spectrometer configuration
 - you have to execute *cf* and *expinstall*
- it should be used for data processing only?
 - you have to execute *expinstall*

For more details please refer to the next chapter

4.4 Perl and GCC

The packages Perl and GCC are automatically installed in the directory `<TOPSPIN_HOME>/gnu` and `<TOPSPIN_HOME>/perl` where `<TOPSPIN_HOME>` is the directory where *TOPSPIN* is installed (default `/opt/topspin`). Perl and GCC are needed for compiling *TOPSPIN* AU programs.

4.5 Installing TopSpin patches

TOPSPIN patches are parts of **TOPSPIN** which contain bug fixes. They are available on the Bruker WWW and FTP server. We recommend to install new patches whenever the *patchlevel* on your system is older than the one on the WWW / FTP server.

Perform the following steps:

1. Open a shell
2. Type **patchlevel**

This will show you the *patchlevel* on your system.

Download from WWW server

1. Log in on a computer which is connected to the internet
2. Open your web browser
3. Go to <http://www.bruker-biospin.de>
4. Click on **NMR** → **tech support** → **Software** → **Downloads** → **Linux PC** and select the package you want to download

Download from ftp server

Type **ftp ftp.bruker.de**

Name: **ftp**

Password: enter your own full Email address

```
ftp> cd /pub/nmr/topspin/patches/topspin<vers.>/LINUX
```

```
ftp> bin
```

```
ftp> ls
```

(if patches are available, you will see a file `topspin-<vers.>-patchX.tar.gz` where X is the patchlevel)

```
ftp> get topspin-<vers.>-patchX.tar.gz
```

```
ftp> get topspin-<vers.>-patchX.readme
```

```
ftp> bye
```

Install the patches according to the file `topspin-<vers.>-patchX.readme`

You can also get **TOPSPIN** patches from the American FTP server:

<ftp.brucker.com:/pub/nmr/mirror.brucker.de/topspin>

Check our WWW pages for *known bugs* :

<http://www.brucker-biospin.de/NMR/nmrsoftw/passwd/bugs/index.html>

4.6 The usage of the NMR-GUIDE

NMR-GUIDE 4.0 that comes with TOPSPIN 1.3 uses an own web server and the TOPSPIN internal java package. Therefore it is not necessary anymore to have these packages installed/configured separately!

4.6.1 Requirements for NMR Guide 4.0

The using of NMR-GUIDE 4.0 requires the installation of TOPSPIN 1.3. Also be sure that the following components are installed:

1. Web browser *Mozilla*
2. *PDF viewer*

NMR-GUIDE offers among other things a lot of PDF documents. To display these documents it is necessary to install the acrobat reader. For an installation description see chapter 9.2.1.2

3. NMR-GUIDE license (contact license@bruker.de for a three month demo license or order a official license from your local Bruker representative). Note that the NMR-GUIDE license is already part of a full or demo **TOPSPIN** 1.3 license.

Mozilla is an official part of the operating system RED HAT ENTERPRISE LINUX WS 3. The WWW browser gets documents from the locally installed WWW server. The client - server communication is established over a network socket connection, but the communication partners only use the localhost TCP/IP host identifier. These pages are shown in the right-hand side window of the browser.

The WWW browser also gets dynamically created pages from the localhost. These workstation specific pages are coming from a *gserver*. This is a platform independent Java web server written by Bruker BioSpin and acts as a daemon running on the locally installed JAVA virtual machine. The *gserver* is included in the NMR Guide&Encyclopedia package. The table of contents is generated dynamically and displayed in the left-hand side window of the browser.

TOPSPIN commands can be invoked from the WWW browser. This process is also

managed by the gserver. The gserver is an essential component of the NMR Guide&Encyclopedia and requires an NMR-GUIDE license.

4.6.2 Java Requirements for NMR-GUIDE 4.0

For using NMR-GUIDE 4.0 that comes with the TopSpin 1.3 version it is not necessary to install/configure a special Java package. NMR-GUIDE uses the TopSpin internal Java package that is automatically installed and configured during installation of TopSpin.

4.6.3 Working with the NMR-GUIDE

In order to work with the NMR-GUIDE, the following procedure must be followed.

4.6.3.1 Start up of the NMR-GUIDE

1. Start TOPSPIN and enter the command:

ghelp

or enter the command (this allows to start a search within NMR-GUIDE:


help

or select from the *Help* pull-down menu of TOPSPIN the entry:

NMR-GUIDE

2. Your standard web browser will automatically starts and shows the page

<http://localhost:6600/guide/index.html>



You can now use the functionality of NMR-GUIDE

4.6.3.2 Functionality test of NMR-GUIDE

If you want to test whether all functions and auxiliary programs that are used in connection with the NMR-GUIDE are correctly installed, then enter the URL

<http://localhost:6060/guide/test.html>

and follow the instructions.

4.6.3.3 Using the NMR-GUIDE search engine

One of the very important and useful tools of NMR-GUIDE is its search engine. It

offers the possibility to search within thousands of html pages that contain a huge amount of information about NMR spectroscopy.

For using this search engine you can start NMR-GUIDE and enter the keyword therein (or select an entry of the list of commonly used keywords). With NMR-GUIDE 4.0 you have also the possibility to search directly from within TOPSPIN, just enter the command:

help <keyword>

Chapter 5

The TOPSPIN license

5.1 Questions and answers about the TOPSPIN license

1. *Which programs need a license?*

TOPSPIN, TOPSPIN PLOT EDITOR, NMR-SIM and NMR-GUIDE need a license. However, if you order a license for TOPSPIN, you automatically get a license for all four programs.

If you do not have a license for TOPSPIN, it will start up in an emergency mode if your PC is connected to a spectrometer. However, NMR-SIM and NMR-GUIDE will not run without a valid license.

ICON-NMR is always used in connection with TOPSPIN and does not need an individual license.

AURELIA, AMIX, MAXENT, GLP and PARAVISION do need a license which must be ordered for each program separately.

2. *Which license types are available and what programs may be run with each?*

There are six different license types, details of each are given in the following table. For more information please refer to the Bruker license order page on the web server:

<https://www.bruker-biospin.de/NMR/nmr-softw/licenses/index.html>

or contact your local Bruker BioSpin representative.

License type	Contents
Full	Acquisition, Processing, Automation (ICON-NMR), Plotting, Experiment Simulation and NMR-GUIDE
	license period: 15 years
	Floating
Demo	Acquisition, Processing, Automation (ICON-NMR), Plotting, Simulation and NMR-GUIDE
	license period: 3 months
	Node locked, uncounted, free of charge
Processing-Only	Processing and Plotting
	license period: 15 years
	no Acquisition, no NMR-SIM, no NMR-GUIDE
Presenter	Viewing and Plotting
	license period: 15 years
	no Acquisition, no Processing, no NMR-SIM, no NMR-GUIDE
Teaching	NMR-SIM and restricted Processing and Plotting
	license period: 15 years
	TopSpin can only process datasets created with XWIN-NMR \geq 3.1 that are older than 4 weeks or datasets created by NMR-SIM. One TopSpin- plot session (without portfolio editor) can be started, Node locked, count 1.
Developer	Special license for extending TOPSPIN with own functionalities. Allows viewing of data.
	license period: unlimited
	no Acquisition/Processing/Plotting/NMR-SIM/NMR-GUIDE, free of charge

3. What kind of licenses do I need?

If you want to use TOPSPIN 1.x you need a FEATURE line 'TOPSPIN1' in your license.dat file. If you have 'XWINNMR3' you can only start XWIN-NMR 3.x.

If more than one license is available, as many copies of the program as are cov-

ered by the license can be started simultaneously.

4. *Which licensing program is used for TOPSPIN?*

The Flexlm license manager is used and is delivered on the TOPSPIN DVD. The Flexlm package occupies about 3 MByte of disk space. It contains the 'bruker_lmgr' license manager, the 'lmutil' program, Help files and example license.dat file.

5. *What are the requirements for the Flexlm to function:*

Your PC must have an Ethernet card

6. *What is a Floating license?*

Floating licenses are licenses that are available for more than one computer in a network. *Floating* licenses are issued based on the Host-ID of one particular computer. This computer becomes the 'license server'. All computers that can communicate with this 'license server' computer can use the licenses managed by this computer. For this to work, the same license file must be installed on all participating computers, including the 'license server'. The Flexlm license manager software must only run on the 'license server' computer.

A *Floating* license can have just one license for a program like TOPSPIN. In such a case, the respective program can be started once on one computer in the network. If more than one license is available, as many copies of the program as are covered by the license can be started simultaneously. For example, if you have a license file with 5 licenses for TOPSPIN, TOPSPIN can be started 5 times from any 5 computers in the network.

The number of licenses is stated in the FEATURE for the respective program.

7. *What is the differenz between 1-Server / 3-Server Licenses?*

Floating licenses can be ordered as 1-server or 3-server licenses. In a 1-server license you have only one license server, and in a 3-server license there are three license servers defined in the license file.

- The advantage of a 1-server license is that it is only necessary to have one computer active at a time to make the license available.
- The advantage of a 3-server license is that any of the three license servers can be inactive and the license management is then done by the remaining two servers.
- One requirement for the three license servers is that all of them have TCP/IP active.

- The operating system on the three computers can be the same or can be different (Linux, Windows).

8. What is a Node Locked license?

A node-locked license allows you to start the respective program only on the computer with the Host-ID for which that license was generated. Typically, this license form is used for Demo and Teaching licenses. No other computer can use the license in the way described above for floating licenses.

9. How do I know if my license is Floating or Node Locked?

You can recognize this from the syntax of the license FEATURES in the `license.dat` file:

```
/usr/local/flexlm/Bruker/licenses/license.dat
```

If the entry '`HOSTID=`' appears after the encrypted password, the license is Node locked and can only be used on the local PC

```
FEATURE TOPSPIN1 bruker_ls 0.0 6-aug-2018 3 9B1EA0113CD53A883974 HOSTID=006008d244fe vendor_info=...
```

→ *Node locked* license for three parallel sessions of TOPSPIN that can be started on the local PC

```
FEATURE TOPSPIN1 bruker_ls 0.0 6-aug-2018 3 9B1EA0113CD53A883974 vendor_info=...
```

→ *Floating* license for three parallel sessions of TOPSPIN that can be started on every PC in the local network

10. If I have a Floating license and I have a dual-boot PC (e.g Windows XP and Linux), can I use the license on both operating systems?

Yes, the license is based on the physical address of the Ethernet card and not on the operating system. Copy the license file `license.dat` to:

Windows:

```
c:\flexlm\Bruker\licenses\license.dat
```

Linux:

```
/usr/local/flexlm/Bruker/licenses/license.dat
```

11. How can I determine the hostid of my PC?

See chapter 5.2.2.

12. How do I install a license (full or demo) on my PC?

See chapter 5.2.

13. Does a Demo license require SERVER/DAEMON lines in license.dat?

No, a Demo license is *Node Locked uncounted* and therefore only requires FEATURE lines. If SERVER and/or DAEMON lines exist, they are ignored.

14. *I have a Floating license. Must the file license.dat be identical on all hosts, server and clients?*

Yes, you have to install the `license.dat` on the server and copy it to all clients.

15. *Must the Flexlm license manager run on all hosts in the network?*

No, it only needs to run on the license server. If the license manager also runs on a license client, this is simply ignored.

16. *Must the Flexlm license manager run on a host with a Demo license?*

No, it only needs to run on the server for counted licenses, a demo license is uncounted.

17. *What can I do if TOPSPIN does not start after I installed the license?*

See chapter 12.1

18. *How is it possible that TOPSPIN starts even though I do not have a license?*

If you do not have a license for TOPSPIN, it will start up in an emergency mode if your PC is connected to a spectrometer. Also one TOPSPIN PLOT EDITOR session (without portfolio editor) can be started from within TOPSPIN. However, NMR-SIM and NMR-GUIDE will not run without a valid license.

19. *What can I do if TOPSPIN cannot find the license server?*

Please try to restart the license server:

```
service bruker_lmgr start
```

5.2 Ordering and installing the TOPSPIN license

5.2.1 Ordering a license

You can order a license from Bruker BioSpin Germany. Demo licenses are free of charge, for other licenses you must specify your order number.

Licenses can be ordered from the German Bruker BioSpin web server:

<https://www.bruker-biospin.de/NMR/nmr-softw/licenses/index.html>

5.2.2 How do I determine the correct hostid

The hostid is a 12-digit hexadecimal number. It is the physical address (also called Mac address or Node address) of the Ethernet card. After installing TOPSPIN, including the Flexlm license manager, you can determine the hostid in two ways:

- Start TOPSPIN; a command prompt will pop up, which will list the hostid, if no licence.dat file is available.
- Open a Linux shell and enter:
`/sbin/ifconfig eth0`
 take the numbers of the *HWaddr*.

5.2.3 Installing a Demo license

A Demo license can be installed as follows:

1. Log in as Administrator
2. Install the Flexlm license manager as described in chapter 4.3
3. Open a Linux shell and enter:

```
kedit /usr/local/flexlm/Bruker/licenses/license.dat
```

enter the FEATURE lines, save and exit the file

5.2.4 Installing a full license on a license server

The Flexlm license manager must always be installed and run as a service on the license server. It is delivered on the TOPSPIN DVD and can be installed as described in chapter 4.3. In case the license file identifies the PC as license server, Flexlm is automatically installed as a service which is then automatically started during boot time

On the license server, the file license.dat must contain a SERVER line, a DAEMON line and FEATURE lines. An example would be:

```
SERVER tulip 0060080e830d 1700
DAEMON bruker_ls /usr/local/flexlm/Bruker/bruker_ls
FEATURE TOPSPIN1 bruker_ls 0.0 6-aug-2018 3 0B0E4011FF0CD32FBFF1 \
vendor_info=" for hostid(s) : 0060080e830d" ISSUER=00047573d35e
FEATURE TOPSPIN_ID bruker_ls 0.0 6-aug-2018 3 9B3EA021818829E2716A \
vendor_info=" for hostid(s) : 0060080e830d" ISSUER=00047573d35e
FEATURE TOPSPIN_2D bruker_ls 0.0 6-aug-2018 3 AB3E80014D54F0FB3887 \
```



```
vendor_info=" for hostid(s) : 0060080e830d" ISSUER=00047573d35e
FEATURE XWINPLOT bruker_ls 0.0 6-aug-2018 3 DB4E50F1618A2B2383D3 \
vendor_info=" for hostid(s) : 0060080e830d" ISSUER=00047573d35e
FEATURE TOPSPIN_ACQU bruker_ls 0.0 6-aug-2018 3 9BDE3051D4CB63255DCC \
vendor_info=" for hostid(s) : 0060080e830d" ISSUER=00047573d35e
FEATURE NMRSIM bruker_ls 0.0 6-aug-2018 3 2B4EC041E98523510FCB \
vendor_info=" for hostid(s) : 0060080e830d" ISSUER=00047573d35e
```

where *tulip* is the hostname and *0060080e830d* the hostid of the computer. Note that the second part of the FEATURE lines, the `vendor_info`, does not necessarily exist. In that case a FEATURE line would look like:

```
FEATURE TOPSPIN1 bruker_ls 0.0 6-aug-2018 3 0B0E4011FF0CD32FBFF1"
```

To inform FlexIm about the new license enter in a Linux shell:

```
/usr/local/flexlm/Bruker/lmreread -c /usr/local/flexlm/Bruker/licenses/license.dat
```

5.2.5 Installing a full license on a license client

On a license client, license can be installed in two different ways, as described for a Demo license chapter 5.2.3. The only difference is that you do not have to setup a new file `license.dat`. You can just copy it from the license server which can be another PC (see question 10 in chapter 5.1).

5.2.6 Modifying the license file

If you have modified a license file, you have to inform FlexIm about the modified license with the command:

```
/usr/local/flexlm/Bruker/lmreread -c /usr/local/flexlm/Bruker/licenses/license.dat
```

Chapter 6

Configuration of TOPSPIN

The idea of this chapter is to give you some guidelines to full spectrometer or workstation configuration. We do not discuss the commands in detail here. For a full description refer to the TOPSPIN manual.

6.1 Take over the configuration of XWIN-NMR/TOPSPIN into TOPSPIN 1.3

The decision how to install TOPSPIN is very important in case you have already XWIN-NMR or TOPSPIN installed on the same or another PC. If so, you might want to take over the configuration of a previous XWIN-NMR/TOPSPIN into TOPSPIN 1.3 because of:

- in case the XWIN-NMR installation was used for controlling a spectrometer there are very important configuration files like the spectrometer configuration, shim files, your own pulse programs, parameter sets, etc.
- if XWIN-NMR was used for processing only there can be also very important configuration files like AU programs, macros, plot layouts etc.



Note: Transferring files could cause permission problems especially if they are transferred from one PC to another.

6.1.1 You have a XWIN-NMR configuration that you like to use for TopSpin (on the same PC)

You can take over the configuration during the Installation of TOPSPIN 1.3 (see chapter 4.3) or after the installation with the AU program *nmr_save*.

During the installation you will be asked if you want to import a spectrometer configuration. If you want to take over the configuration after installation is finished, start TOPSPIN

- and type in the command:

nmr_save

- or click

Options → ***Spectrometer Tools*** → ***Save/Restore configuration***

A window 'Password request' will appear. Please enter the NMR administration password and click ***OK***.

A window 'NMR_Save' appears.

- Accept or modify the 'Location of Backup file' (default: <TOPSPIN-HOME>\nmr_backup).
- Enter the location of the installation to be saved (e.g. /opt/xwinnmr).
- Click on ***Save***.

This will create a tar file of all necessary information from the TOPSPIN installation. The tar-file is called *nmr_backup.tar* and will be stored in the directory you selected.

To take over this information into TOPSPIN, start *nmr_save* again and ...

- Make sure that 'Location of Backup file' contains the correct path to your backup file
- Make sure that the restore destination is correct (previous settings might be overwritten in the restore destination!)
- Click on ***Restore***.

.After this process is finished, execute ***cf***.

If you like to get more information about the command *nmr_save* refer to chapter 6.6

6.1.2 You have a XWIN-NMR configuration that you like to use for TOPSPIN (on another PC)

Before you start installation of TOPSPIN on PC-A you have to copy the two directories:

```
<xwinnmrhome>/conf and <xwinnmrhome>/exp
```

from the XWIN-NMR installation on the PC-B into the directory where TOPSPIN will be installed. The default installation directory of TOPSPIN is for Linux:

```
/opt/topspin
```

In case you will use the default installation path you have to create these directories and as a second step copy the two directories `conf` and `exp` into them. After these two directories are available on PC-A you can start the TOPSPIN installation.

6.2 Configuring TOPSPIN

There are two ways to start TOPSPIN:


1. Click the icon TOPSPIN on the desktop
2. Click **Start** → **All Programs** → **Bruker TOPSPIN** → **TOPSPIN<vers.>**

After the first startup of TOPSPIN you have to configure the software (depends on the hardware of your spectrometer). If you do this, you are mainly in one of the four following situations:

- You have installed the new software on the same PC as a previous XWIN-NMR/TOPSPIN installation, so the old spectrometer configuration is available on the same PC → go to chapter 6.2.1
- You want to configure a processing-only PC like your spectrometer PC, so you have the spectrometer configuration on another computer → go to chapter 6.2.2
- You have to configure the software on a new hard disk without a backup of the spectrometer configuration, so you have no configuration files anymore → go to chapter 6.2.4
- You have a new spectrometer hardware component, so you have no configuration file so far → go to chapter 6.2.5

6.2.1 If your PC controls a spectrometer


1. Be sure that the spectrometer is booted (see chapter 9.3.0.3 and 12.3.4.2)
2. Open a dataset
 - Type **new** and create a new dataset.
 - You can also select a dataset (e.g. *exam1d_1H'*) acquired with another TOPSPIN or XWIN-NMR version.
3. Enter **cf** (The command **cf** now configures the software according to the hardware of your spectrometer):
 - a) A window will appear:
Password request. Please enter the NMR ADMINISTRATION PASSWORD
 Type in the NMR ADMINISTRATION PASSWORD. Click **OK**.
 - b) A window Cf will appear:



Most likely you will select an existing configuration and press **Edit** to confirm this configuration.

- c) A window will appear:
 The selected configuration is displayed. Choose the type of your spectrometer e.g. **AVANCE** → **ENTER**. Check or type in the 1H frequency of the magnet e.g. **500.13** → **ENTER**. Use the debug mode only in a case

of trouble-shooting.



d) A window will appear:

Specify the channel to which external devices are connected. The RS232/485 channel for external devices-Table appears.

Check if the tty-numbers of the hardware components are set correctly.

e) A window will appear:

Additional configuration. Security check. Enable peak power check (POWCHK). Enable puls power check (PULCHK).

Choose the options if you need this security options.

f) Only if you have a sample changer! A window will appear:

Should the Sample Changer control the Lift? (yes/no)

The answer depends on how the air flow is connected, normally you will answer: **yes** → **ENTER**

In this cases the question appear:

Delay between SX and next command [sec]?

Type in the length of the delay e. g. **10** → **ENTER**

g) The *nuclei-table* appear.


Check if the frequency of the nuclei are set correctly. If in doubt click **Re-store** then click **Save**

h) The results of *cf* appears in a text window. Check this list carefully and compare the result with an list from an earlier configuration (if available) click **Print** and store the paper output with your other spectrometer documentation, then click **OK**. *These data are saved in uxnmr.info.*




If during 'cf' an error message occurs that invites you to do a firmware update see chapter 6.3

After finishing *cf* the window *configure* will appear showing a list of configuration steps. Accept the default selection and click **Start** (You will be asked for the NMR ADMINISTRATION PASSWORD a couple of times):




- Continue with ***expinstall***

- a) Type in the NMR ADMINISTRATION PASSWORD and **ENTER**
- b) Click **Next** if you either have no modified Bruker BioSpin-Parameterfiles,




AU-programs and Pulse-programs (If you have some, you have to archive them in a different directory before executing *expinstall*). Your individual Parameterfiles, AU-programs and Pulse-programs that have different names to the original Bruker BioSpin files will not be destroyed during the installation


- c) In the window that appears → Check *Installation for Spectrometer* → Click *Next*




- d) In the window that appears → Check *High Resolution Systems* → Click *Next*




- e) In the dialog box that appears → Select the spectrometer configuration name → Click *Next*



In the dialog box that appears → Check the items you want to install or, accept the default selection → Click *Next*




- f) In the dialog box that appears → Select the default printer and plotter and the desired paper format → Click *Next*




- g) In the window that appears → Enter the desired spectrometer frequency, acquisition mode and pre-scan-delay or, accept the default values

→ Click *Next*



- h) In the information box that appears → Check the configuration selection and, if it is correct → Click *Finish*



While *expinstall* is running you can go ahead with the other configuration commands.

- Continue with **edsolv**

Solvent	Description	Reference Number
Acetic	acetic acid-d4	[01]
Acetone	acetone-d6	[02]
CD3CN	acetonitrile-d3	[03]
C6D6	benzene-d6	[04]
CDCl3	chloroform-d	[05]
D2O	deuteriumoxide	[06]
pC6D4Br2	p-dibromobenzene-d4	[07]
oC6D4Cl2	o-dichlorobenzene-d4	[08]
DEE	diethylether-d10	[09]
DME	dimethylether-d6	[10]
DMF	dimethylformamide-d7	[11]
DMSO	dimethylsulfoxide-d6	[12]
Dioxane	dioxane-d8	[13]
EtOD	ethanol-d6	[14]
MeOD	methanol-d4	[15]
CD2Cl2	methylenechloride-d2	[16]
Pyr	pyridine-d5	[17]
THF	tetrahydrofurane-d8	[18]
Tol	toluene-d8	[19]
C6F6	hexafluorobenzene	[20]
TFA	trifluoroacetic acid-d	[21]
C6H5CF3	trifluorotoluene	[22]
C6H5F	fluorobenzene	[23]
pC6H4F2	p-difluorobenzene	[24]
H2O+D2O	90%H2O and 10%D2O	[25]
CH3CN+D2O	HPLC Solvent (Acetonitril/D2O)	[26]
CH3OH+D2O	HPLC Solvent (Methanol/D2O)	[27]
None	no Solvent	[28]

Print New Entry Delete Restore Save Close

Check the list of lock solvents and click **Save** if it is correct → Click **Close**

- Continue with **edhead**

Probe		
No current probe defined.		
Please select a probe then press the button 'Define as current probe'.		
5 mm 1H		[01]
6 mm 1H		[02]
5 mm Dual 13C/1H		[03]
5 mm Dual 19F/1H		[04]
5 mm QNP 1H/15N/13C/31P		[05]
5 mm QNP 1H/13C/31P/19F		[06]
5 mm QNP 1H/13C/31P/19F Z-grad		[07]
5 mm QNP 1H/29Si/13C/31P		[08]
5 mm Multinuclear		[09]
5 mm Multinuclear inverse		[10]
5 mm Multinuclear inverse Z-grad		[11]
5 mm TXI 31P Z-grad		[12]
5 mm TXI 13C Z-grad		[13]
10 mm 1H		[14]
10 mm Multinuclear inverse		[15]
10 mm 2H		[16]
10 mm 13C		[17]
10 mm 15N		[18]
10 mm 31P		[19]
10 mm QNP 1H/15N/13C/31P		[20]

Define as current probe Edit Probe Parameters Exit

Select the current probe, click the button **Define as current probe** → Click **Exit**

- Continue with *edprosol1*


Check the parameters for your favorite nuclei, change them and calculate them, if necessary → Click *Save* to store to disk → Click *Exit*

- Continue with *edlock*

Solvent	BSMFSI.	LockPo.	LoopGain	LoopTime	LoopFit	LockPh.	Distance	Ref	Width	Ref-Shift
Acetic	171	-40.0	-32	0.136	200	-1	2.0300	0.0000	0.5000	0.0000
Acetone	171	-40.0	-32	0.136	200	-1	2.0400	0.0000	0.5000	0.0000
Aceton	171	-40.0	-32	0.136	200	-1	2.0400	0.0000	0.5000	0.0000
CDCl3	171	-25.0	-32	0.136	200	-1	7.2400	0.0000	0.5000	0.0000
CD2Cl2	171	-30.0	-32	0.136	200	-1	5.3200	0.0000	0.5000	0.0000
CD3CN	171	-40.0	-32	0.136	200	-1	1.9300	0.0000	0.5000	0.0000
C6D6	171	-26.0	-32	0.136	200	-1	7.2800	0.0000	0.5000	0.0000
D2O	171	-20.0	-32	0.136	200	-1	4.7000	0.0000	0.5000	0.0000
H2O+D2O	171	-23.0	-32	0.136	200	-1	4.7000	0.0000	0.5000	0.0000
DEE	171	-30.0	-32	0.136	200	-1	1.0700	0.0000	0.5000	0.0000
DME	171	-35.0	-32	0.136	200	-1	3.3000	0.0000	0.5000	0.0000
DMF	171	-25.0	-32	0.136	200	-1	2.9100	0.0000	0.5000	0.0000
DMSO	171	-25.0	-32	0.136	200	-1	2.4900	0.0000	0.5000	0.0000
Dioxane	171	-30.0	-32	0.136	200	-1	3.5300	0.0000	0.5000	0.0000
EtOD	171	-30.0	-32	0.136	200	-1	1.1100	0.0000	0.5000	0.0000
MeOD	171	-35.0	-32	0.136	200	-1	3.3000	0.0000	0.5000	0.0000
THF	171	-25.0	-32	0.136	200	-1	1.7300	0.0000	0.5000	0.0000
Tol	171	-40.0	-32	0.136	200	-1	2.0900	0.0000	0.5000	0.0000
Pyr	171	-25.0	-32	0.136	200	-1	8.7100	0.0000	0.5000	0.0000
CH3CN+D2O	171	-30.0	-32	0.136	200	-1	4.7000	4.7000	0.5000	0.0000
MeOH+D2O	171	-30.0	-32	0.136	200	-1	4.7000	4.7000	0.5000	0.0000

Check the lock table values (field, lock phase, ...), change them and click **Save** if necessary, otherwise click **Abort**

- Continue with **edscon**



Check the spectrometer constants, change it and click **Save** if necessary, otherwise click **Cancel**


After **cf**, **expinstall**,... have finished, insert a sample and make sure that you can lock on the lock solvent. Type **ii** and read with **rpar** a standard proton parameter file, type **getprosol** and **rga** and collect a normal ^1H NMR spectrum.

6.2.2 If you want to configure a processing-PC like your spectrometer-PC:


If the PC does not control a spectrometer but is used in connection with a particular spectrometer, e.g. for setting up experiments, processing and or plotting, you should configure it as that spectrometer (see FAQ 210 on www.bruker-biospin.de): for example the spectrometer is a AV300 configured with the name 'av300'. Do the following.

- Copy the configuration directory '`<TOPSPINHOME>/conf/instr/av300`' from the spectrometer AV300 to the workstation.
- On the workstation paste this 'av300' directory to '`<TOPSPINHOME>/conf/instr`'.
(most likely this is the path: `/opt/topspin/conf/instr`)
- Start TOPSPIN as normal user.
- Perform **expinstall** as follows:
Type in the command line **expinstall** or
Click **Open** → **Setup** → **Experiment Installation [expinstall]**
- Type in the NMR ADMINISTRATION PASSWORD and **ENTER**


- f) In the appearing information box → Click **Next** if you either have no modified Bruker BioSpin-Parameterfiles, AU-programs and Pulse-programs (If you have some, you have to archive them in a different directory before executing **expinstall**). Your individual Parameterfiles, AU-programs and Pulse-programs that have different names to the original Bruker BioSpin files will not be destroyed during the installation




- g) In the appearing window → Check *Installation for Datastation (Customize)* → Click **Next**




- h) In the appearing window → Check *High Resolution Systems* → Click *Next*




- i) In the appearing dialog box → Select the spectrometer configuration name → Click *Next*




- j) In the appearing dialog box → Check the items you want to install or, accept the default selection → Click *Next*




- k) In the appearing dialog box → Select the default printer and plotter and the desired paper format → Click *Next*



- l) In the appearing window → Enter the desired spectrometer frequency, acquisition mode and pre-scan-delay or, accept the default values → Click *Next*



- m) In the appearing information box → Check the configuration selection and, if it is correct → Click **Finish**



The installation of the selected items, will start now. Wait until this process has finished.

For more details on **expinstall**, please refer to the description of this command in the Acquisition Reference manual.

6.2.3 If you want to use a default configuration on a datastasion

A default configuration can be used on a datastasion. It is automatically performed (no **cf** required) during the installation of TOPSPIN on a new computer, a new disk or in a new TOPSPIN installation directory. The default configuration name is *Bruker_default_av500* and corresponds to a Avance 500 MHz spectrometer.

For manual or interactive data processing, the automatic default configuration is sufficient. If, however, you want to use AU programs, you must execute **expinstall** once, selecting *Installation for Datastasion (Default)*.

6.2.4 If you do not know the configuration of your spectrometer

If you do not have any spectrometer configuration information, e. g. after a head-crash, then you have to recreate the configuration information.

In this situation there are two possibilities:

- a) your spectrometer has a standard configuration
- b) your spectrometer has a non-standard configuration

Case a)

If your spectrometer has a standard configuration, then you can easily do **cf** and if you give the correct information for the '*Type of Spectrometer*' you will get the correct list of your hardware in the '*RS232 table*'. You only have to type in the correct tty's of the interfaces.

Case b)

If your D*X spectrometer has a non-standard configuration, then the hardware configuration is read during the **cf** from the file 'hardware_list' file in the directory

```
<TOPSPINHOME>\conf\instr\<spectrometer-name>
```

TOPSPIN 1.3 offers the default file 'hardware.exam' which contains a list of all available hardware components which are important for the **cf** command. This file was generated in XWIN-NMR with the command **cf makelist**, now it is available in:

```
<TOPSPINHOME>\conf\instr
```

- Open this file with a text editor, e. g. kedit, and delete all lines except those that describe the hardware components of your spectrometer. You can get this information in the following way:

Check the list in the file 'hardware.exam' with the components of your spectrometer and/or with the Bruker BioSpin delivery note(s). Please carefully check the PN-C/D/E and ECL numbers of the components

- **Save and exit** from 'hardware.exam'
- Copy the modified file 'hardware.exam' into the directory

```
<TOPSPINHOME>\conf\instr\<spect>
```

- Rename the file

```
<TOPSPINHOME>\conf\instr\<spect>/hardware.exam to
```

```
<TOPSPINHOME>\conf\instr\<spect>/hardware_list
```

6.2.5 If you added a new hardware component to your spectrometer



If you want to add a new hardware component it is recommended to order somebody of your local Bruker BioSpin office to do this for you.

If you have to add a new hardware component to the configuration file of the software, do the following:

TOPSPIN 1.3 offers the default file 'hardware.exam' which contains a list of all available hardware components which are important for the **cf** command. This file was generated in XWIN-NMR with the command **cf makelist**, now it is available in:

```
<TOPSPINHOME>\conf\instr
```

- Open this file with a text editor, e. g. kedit, and search and copy these line(s) that describe the new hardware component(s). Please carefully check the PN-C/D/E and ECL number(s) of the component(s).

- Copy the file

```
<TOPSPINHOME>/conf/instr/spect/hardware_list to
```

```
<TOPSPINHOME>/conf/instr/spect/hardware_list_sav
```

- Open the file

```
<TOPSPINHOME>/conf/instr/spect/hardware_list
```

and paste the new entries to the correct positions.

6.3 Firmware update of the spectrometer hardware

This chapter describes how to update the firmware of certain spectrometer components. This is only necessary if the **cf** command has crashed and popped up a message window which informs that you have to do a firmware update of a component.



It is strongly recommended to do the firmware update only in this case, because the hardware component could get unusable if the update would be done incorrectly.

If you answer the question about the setting of the default path during the installation of the NMR Suite with **yes**, you can start the following service tools in a shell by typing their names:

```
acb
bsms
grad
hpcu
hppr
rx22
```

If you answer **no** you have to type in the name with the entire path:

```
<TOPSPIN_HOME>/prog/bin/scripts/<service-tool-name>
```

The following example shows the check of the *bsms service tool* with the entire pathname

1. Open a shell
2. Type `<TOPSPIN_HOME>/prog/bin/scripts/bsms`
 - An automatic interactive procedure starts in this shell
3. Type **y** → **ENTER**
 - A dialog comes up in the shell
4. To check and eventually download a new firmware (from *TOPSPIN*) type **2** → **ENTER**
5. The check starts and produce some messages e. g.:
 - if there is a problem with a component you get an information to contact your service office
 - if you need a firmware update you will get a message together with the length of time for the download. It is recommended to do the download as soon as possible.



During the download it is absolutely necessary that spectrometer and workstation are not disturbed! Make sure that no one pulls a cable out of the spectrometer or will crash the Linux system while using it in parallel! (Yes we know this seems very improbable, but if the download would

crash, your board might be unusable afterwards)

6. press \mathcal{Y} → ENTER to start the download

A 'normal' error during downloading the new firmware (you will see a comment in the shell) requires to do the download once again!

If you have done a download of newer firmware, you have to do \mathcal{CF} again!

If you do not need any downloads, you can follow the instructions in the command prompt to exit this tool.

6.4 Configuring the BSMS-Keyboard

For a full description how to configure the BSMS-Keyboard see your

- hard copy of *BSMS-User Manual* or
- the respective entry on the *BASH CD* (Order Number Z36541)

A configuration step that should be done periodically is the spin rate calibration. For this you have to enter the **Menu** mode of your BSMS keyboard:

press simultaneously **2nd** and \mathcal{Y}^3 , then rotate the control knob as long as you see **1.Sample** in the display, press **2nd** to go into this directory. The display will show **1.1 Spin calib.** press **2nd** to start this process and wait for ~ 3 minutes.

In this Menu mode, you can also define the helium warning level. If the liquid helium reaches the critical point, the BSMS-Keyboard starts beeping. A special password, the so called security code, is required to be able to change the setting. If in doubt about the setting, contact your service representative.

For an additional possibility to get a warning about a low helium level see chapter 16.1.

6.5 A log of the entire system configuration

The *showconf* command provides you with various information about your hardware and software, both of the PC and of the spectrometer. Here is how showconf works: You type *showconf* in a shell - and wait the about 2 minutes until it has finished. After that showconf has produced a file called `confout` that contains all you need. `Confout` is located in your home directory. The size of this report is

usually less than 100 kbytes.

DO KEEP A COPY OF THIS FILE IN A SAFE PLACE!

6.6 Backup of VIFs (=Very Important Files) with `nmr_save`

It is not only necessary to create some backups for the operating system like Boot Disks and to store the acquired datasets.



It is strictly recommended to create a frequently backup of the 'VIFs' - the 'Very Important Files' of the configuration of TOPSPIN.

There are many reasons like e. g. head crash or theft of the hard disk that makes it important that a weekly or monthly backup is made of important data, so that you would be able to restore the spectrometer configuration in a very short time.

TOPSPIN offers an AU program¹ (`nmr_save`) for saving those important files that are not delivered by BRUKER. This AU program can save and extract user specific data from the TOPSPINHOME directory and puts them into a tar-file. Files that start with an "." are ignored.

The AU program can be started as follows:


- Start TOPSPIN and type in the command:

`nmr_save`

- or click

Options* → *Spectrometer Tools* → *Save/Restore configuration


A window 'Password request' will appear.



Please enter the NMR administration password and click **OK**.

1. BRUKER also offers a script (`xwinnmr.save`) which can be used for saving the 'VIF's (see chapter 16.2).

A window 'NMR_Save' will appear.



Accept or modify the 'Location of Backup file' (the default is: <TOPSPIN-HOME>\nmr_backup). Enter the location of the installation which should be saved, select the respective spectrometer configuration and click on **Save**.

This will create a tar file of all necessary information from the selected installation/spectrometer configuration. The tar-file is called `nmr_backup.tar` and is stored in the directory:

```
<TOPSPINHOME>\nmr_backup
```

Files and directories will be saved into the special tar-file `nmr_backup.tar` which is, by default, in the directory `TOPSPINHOME/nmr_backup`.

The following directories are saved completely (hidden files are ignored!):

```
<TOPSPINHOME>/exp/stan/nmr/lists/group/
<TOPSPINHOME>/exp/stan/nmr/lists/scm/
<TOPSPINHOME>/exp/stan/nmr/lists/bsms/
<TOPSPINHOME>/exp/stan/nmr/lists/mac/
<TOPSPINHOME>/exp/stan/nmr/lists/fl/
<TOPSPINHOME>/exp/stan/nmr/lists/ds/
```

```

<TOPSPINHOME>/exp/stan/nmr/py/
<TOPSPINHOME>/prog/server/
<TOPSPINHOME>/conf/instr/probehead/
<TOPSPINHOME>/conf/instr/<spect>/prosol/
<TOPSPINHOME>/conf/instr/<spect>/cortab/
<TOPSPINHOME>/conf/instr/<spect>/inmrusers/
<TOPSPINHOME>/conf/instr/autoshim/refmaps/
<TOPSPINHOME>/data/final/nmr/protocolfiles/
<TOPSPINHOME>/QTP/

```

<spect> represents the current spectrometer name or the name specified explicitly when **nmr_save** is started.

From the following directories every file is saved that does not belong to the genuine Bruker distribution:

```

<TOPSPINHOME>/exp/stan/nmr/lists/pp/
<TOPSPINHOME>/exp/stan/nmr/lists/cpd/
<TOPSPINHOME>/exp/stan/nmr/lists/gp/
<TOPSPINHOME>/exp/stan/nmr/lists/wave/
<TOPSPINHOME>/exp/stan/nmr/au/scl/
<TOPSPINHOME>/exp/stan/nmr/au/src/
<TOPSPINHOME>/exp/stan/nmr/par/
<TOPSPINHOME>/plot/layouts/
<TOPSPINHOME>/prog/tcl/xwish3_scripts/

```

Additionally, the following files are also saved:

```

<TOPSPINHOME>/conf/instr/<spect>/uxnmr.par
<TOPSPINHOME>/conf/instr/<spect>/uxnmr.info
<TOPSPINHOME>/conf/instr/<spect>/specpar
<TOPSPINHOME>/conf/instr/<spect>/2Hlock
<TOPSPINHOME>/conf/instr/<spect>/nuclei
<TOPSPINHOME>/conf/instr/<spect>/best-nmr/CurrPushSolvent
<TOPSPINHOME>/conf/instr/<spect>/best-nmr/hardware.par

```


```
<TOPSPINHOME>/conf/instr/<spect>/best-nmr/hardware.src  
<TOPSPINHOME>/conf/instr/<spect>/best-nmr/racks.src  
<TOPSPINHOME>/conf/instr/<spect>/best-nmr/scanner.src  
<TOPSPINHOME>/conf/instr/<spect>/best-nmr/solvents.add  
<TOPSPINHOME>/conf/instr/<spect>/best-nmr/solvents.add  
<TOPSPINHOME>/conf/instr/<spect>/best-nmr/solvents.use  
<TOPSPINHOME>/conf/instr/<spect>/rs232_device/best  
<TOPSPINHOME>/conf/instr/<spect>/rs232_device/bestscan  
<TOPSPINHOME>/conf/instr/probehead  
<TOPSPINHOME>/exp/stan/nmr/lists/solvents  
<TOPSPINHOME>/exp/stan/nmr/lists/probeheads  
<TOPSPINHOME>/conf/instr/<spect>/hardware_list  
<TOPSPINHOME>/prog/logfiles/heliumlog  
<TOPSPINHOME>/prog/logfiles/heliumlog.err  
<LM_LICENSE_FILE>1
```

For safety reasons the tar-file will not be replaced by executing ***nmr_save*** once again. So if you like to use this AU program several times, you have to move the `nmr_backup.tar` file from the directory `<TOPSPINHOME>/nmr_backup` on a backup medium (floppy, CD ...) and then start the script once again or you have to select the 'Overwriting existing backup file' option.

Example for transferring configuration information from a XWIN-NMR installation in `/opt/xwinnmr` and the spectrometer configuration 'spect' into a TOPSPIN installation in `/opt/topspin`:

1. This is the environment variable for the license file, in a default installation this is:
`/usr/local/flexlm/Bruker/licenses/license.dat`

Backup :



Restore:



If 'nmr_backup.tar' does not exist in the backup directory, the backup file 'xwin_backup.tar' will be taken (if existing). After this process is finished, execute **cf**.

Chapter 7

Plotting and Printing with TopSpin

7.1 Printing / plotting data

7.1.1 Plotting and printing


Under RED HAT ENTERPRISE LINUX WS 3, all common printer drivers can be used for plotting and printing. When a printer is installed under RED HAT ENTERPRISE LINUX WS 3, it can be used for plotting from TOPSPIN and TOPSPIN PLOT EDITOR.

7.1.2 Print/Plot from the Menu

The current data window can be printed as follows:

1. From the TOPSPIN menu:
 - Click the print button in the upper toolbar
 - or Click **File** → **Print**
 - or Enter *print* or *Crtl-p*

All these actions are equivalent; they open the Print dialog box.



2. In the Print dialog box:

- a) Select ***Print active window [prnt]***
- b) Click ***OK***

Before printing starts, the operating system print dialog box will appear to select e.g. the printer and printer properties.

The Print dialog box contains two further options:

- ***Print with layout - start Plot Editor [plot]***

If you select this option and click OK, the Plot Editor will be started. This option is equivalent to entering plot on the TOPSPIN command line.

- ***Print with layout - plot directly [autoplot]***

Selecting this option activates the Plot Editor layout list box. Select the desired layout and click ***OK*** to print. Standard layouts are delivered with TOPSPIN. They use the Windows default printer. User defined layouts use the printer defined in the TOPSPIN PLOT EDITOR. On a 1D dataset, only 1D layouts are listed, on a 2D dataset only 2D layouts are listed etc.

For each of the last two options one of the following parameters is required.

Use plot limits:

- ***from screen / CY***

The plot limits and maximum intensity are used as they are on the screen

(processing parameter F1P, F2P and CY, respectively)

- **from Plot Editor Reset Actions**

The plot limits and maximum intensity are set according to the TOPSPIN PLOT EDITOR Reset Actions (right-click inside the TOPSPIN PLOT EDITOR data field and choose Automation to set the Reset Actions).

- **as saved in Plot Editor**

The plot limits and maximum intensity are set in the specified layout

Fill dataset list:

- **from your default portfolio**

The portfolio contains the current TOPSPIN dataset plus the data from the default Plot Editor portfolio

- **from portfolio saved in dataset**

The portfolio contains the current TOPSPIN dataset plus the data from the portfolio stored in this dataset

7.1.3 Plot data from the Processing guide

Printing/plotting data can be done from the Processing guide by clicking the **Plot/Print** button. If **Automatic mode** is checked, the active data window will be printed as it appears in the screen. If **Automatic mode** is unchecked, you will get the dialog box.

7.1.4 Plot data with the TOPSPIN PLOT EDITOR

Although the **print** command allows you to use or start the TOPSPIN PLOT EDITOR, you can also start it from the command line (enter **plot**). This will open the TOPSPIN PLOT EDITOR from which you can create layouts and plot data. The complete TOPSPIN PLOT EDITOR functionality is described in its online manual that can be opened from the TOPSPIN PLOT EDITOR Help menu.

7.1.5 Store (Export) a Data Window as Graphics File

The clipboard and metafile formats are resizable vector formats. In addition to this, TOPSPIN allows you to save the contents of a data window in a graphics file of selectable type, e.g. .png, .tif, .wmf etc. To do that:

1. Click **Edit** → **Export**

2. Navigate to the storage folder
3. Enter the destination filename and extension
4. Click ***Export***

The resolution of such a screen dump equals the resolution of your screen. When you import a graphics file into an other program, you may loose information when resizing the graphics.

7.1.6 TOPSPIN PLOT EDITOR

1. TOPSPIN PLOT EDITOR uses the printer driver system of the operating system for any printer output. This means that any printer for which a system driver is correctly installed can be selected as print device. Note that the printer name in a layout file written by the Windows XP version is probably not known by the Linux version of TOPSPIN PLOT EDITOR because of the different printer handling. Apart from that, layout files can be freely interchanged between the two platforms.
2. TOPSPIN PLOT EDITOR for Linux offers an internal Clipboard. You can copy objects from one Plot Editor window into another by using copy and paste.
3. Another way of exporting graphics is the menu option '*File/Print*' which allows you to print the complete layout into a file. The format can be chosen from a list (e. g. PostScript, Epsi, PDF, JPEG, PCX, PNG or TIFF). Most current graphical-oriented windows applications can import this format. This provides a convenient way of importing TOPSPIN PLOT EDITOR graphics into programs like Word or Powerpoint.
4. With the TOPSPIN PLOT EDITOR you can directly send your current plot as Email. Just click **File** → **Send** and TOPSPIN PLOT EDITOR will pop-up your default mail client and add the current plot as an EMF attachment.

Chapter 8

User Management

8.1 Adding a new user account


If you want to add a new user¹ and this new user already exists on other hosts in the network, make sure that you use the same *User ID* and *User Name* on all hosts in the network. It is also useful to use the same *Group ID* and *Group Name*, especially, if you want to mount directories via NFS between IRIX and Linux systems.

Perform the following steps:

1. Click <red hat> → **System Settings** → **User and Groups**
2. If you are not logged in as root, a window *run as root - KDE su* will appear
 - enter root password


1. The default configurations of a new user are copied from `/etc/skel/` into the home directory of the new user. If you want to change the default informations, e. g. by adding an icon that starts **TOPSPIN** copy the respective file into `/etc/skel/`

3. A window *KDE User Manager* will appear:



- Click **Add User**

4. A window *Create New User* will appear



- Enter the new **User Name**
- Enter **Full Name** of the new user
- Enter and confirm the **password**
- Make sure that **Create home directory** is enabled and that a legal path is entered in the field **Home directory**. The default `/home/` will create the home directory:

/home/<USER_NAME>

- The field **Create a private group of this user:**
 - a) enabled (default):

This configuration in conjunction with the default umask (002) will lead to the situation, that only the user is able to write in those files he created with the TopSpin. That means no other user could work with his spectroscopic data.
 - b) disabled:


New user will be assigned to the group **users**

This configuration in conjunction with the default umask (002) will lead to the situation, that the user and all members of the group *users* are able to write in those files he created with the TopSpin. That means all other users of his/her group can work with his spectroscopic data.
- Enable the field *Specify user ID manually* if necessary (see introduction of this chapter) and type in the desired ID of the new user
- click **OK**


8.1.1 Change group membership of an existing user

1. Click <red hat> → **System Settings** → **Users and Groups**
2. If you are not logged in as root, a window *run as root - KDE su* will appear
 - enter root password
3. A window *KDE User Manager* will appear:
 - Click on an existing user and click on **Properties**

4. A window *User Properties* will appear:



- Click on **Groups**



You will see in the appearing window the field Primary Group that shows the original group of the user.

- In the field *Select the groups that the user will be a member of* you can now add the user to other groups simply by selecting the respective field of the list of group(s)
- Click **OK**

Part II

RED-HAT Enterprise Linux WS 3 Installation-Tips and Configuration

Chapter 9

RED HAT

Enterprise Linux WS 3

9.1 Hints for installation of Red Hat

A Linux NMR workstation purchased from Bruker comes with a pre-installed Linux distribution, currently RED HAT ENTERPRISE LINUX WS 3.

If you want to install Red Hat by yourself, Bruker recommends to use the installation guide of the Red Hat Documentation CD. This offers, in different languages, a very easy and detailed description of all necessary steps. Those steps which are recommended by Bruker are listed below, for all other configuration steps you can choose your personnel settings.

9.2 Scratch installation of Desktop computer

1. Start the installation of RED HAT ENTERPRISE LINUX WS 3 as described on the Red Hat Documentation CD
2. In the window *Mouse Configuration*:
select your mouse (the default *generic mouse* might not support all features) (e.g

select *Microsoft Intellimouse* if you have a wheel mouse, otherwise the middle button will not work correctly)

3. In the window *Disk Partitioning Setup* choose **Manually partition with Disk Druid**

4. In the window *Partitioning*:

Brucker BioSpin recommends to create the following partitions on a 80GB hard-disk. For all these partitions you should select the button: **Force to be a primary partition** and as **File System Type** ,*ext3*¹:

- a) a partition */boot* with 125MB size
- b) a partition */* with 16GB size
- c) define the swap partition */swap* with 2GB size by selecting *swap* as **File System Type**
- d) a partition */opt* should use the remaining space (enable **Fill to maximum allowable size**)

5. In the window *Boot Loader Configuration* accept the default *Use GRUB as the boot loader*²

6. In the window *Network Configuration* disable *Active on Boot* for *eth1* if your PC does not control a spectrometer (e. g. notebook or workstation).

If the PC is controlling a spectrometer you should select *Active on Boot* for *eth1*. Additionally you can click on *Edit* of *eth1* and configure it for spectrometer control. For this you can type in the IP address and subnet mask:

IP-address (eth1): 149.236.99.1

Subnet mask (eth1): 255.255.255.0

The hostname should be created by DHCP (if available). Changing the hostname to another name can be done in the installed system, see chapter 10.1

7. In the window *Firewall*, select the firewall configuration you prefer.

If you select other than *No Firewall* and the installation is for spectrometer con-



1. *ext3* file system is a journaling filesystem. The advantage of a journaling filesystem is, that an accidental power off would not end in most cases to a file system check and missing data.
2. GRUB is more convenient than LILO especially for dual boot computer. if you want to use Norton GHOST 2003 for creating images of the harddisk you have to use LILO. PCs delivered from Bruker Germany have Norton Ghost 2003 attached. For more details see chapter 15.8

trol, you have to enable your spectrometer communication for full access. Click on **eth1** in the section:

If you would like to allow traffic from a device, select it below

8. In the window *Timezone Selection*:

it is necessary to configure the *Location*, but do not enable *System clock uses UTC*

9. In the window *Package Defaults*:

we recommend to select: *Customize the set of packages to be installed*

and in the following window *Package Group selection*

scroll down to the end of the list and select **Everything**¹

10. During the installation process, you will be asked to insert all four installation CDs one after the other . After finishing with the fourth CD you will be asked to insert the first CD once again.

11. The window **Congratulations** appears and informs you about a successful installation. You have to remove the CD and press **Exit**, this will reboot the PC.

12. When the system is coming up again, a window *Welcome* appears

13. You have to accept the license agreement, to adjust date and time and you will be asked for creating a non-administrative user account.

14. In the window *Red Hat Network* you can select to register at the Red Hat Network right now (**Yes, ...**) or you can choose **No...** in case you have currently no network connection etc.

15. A window *Additional CDs* appear

On HP computers it is necessary to install the drivers/update CD that can be downloaded from the HP web side

<http://welcome.hp.com/country/us/en/support.html>

Enter there the type of your PC (e. g. xw4100 for the currently delivered HP xw 4100) and follow the instructions of the following pages.

If you have the correct CD for your PC available or any other CD that should be installed now (e. g. a Red Hat Extras CD), you have to insert the respective CD then click on the **Install...** button right to the text *Additional CDs*

16. You will be guided trough the installation process

1. This installation requires about 4GB space on the harddisk

17. After finishing the installation process of the HP driver/update CD you will be asked to reboot now or later. In case you want to install further CDs it is recommended to click on reboot later, because otherwise this process is stopped.
18. After finishing this installation you can start TopSpin installation (see chapter 2.2 for a brief description or chapter 4.3 for a detailed description)

9.2.1 Necessary additional packages

Bruker recommends to use the full installation of Red Hat Enterprise Linux WS 3 (see chapter 9.2 step 9). In addition to this complete installation you need for some special cases (see table below) the two packages *bootparamd* and Acrobat Reader. For an installation description see the following chapters.


Package	available on	Necessary for
bootparamd	<i>TopSpin 1.3 DVD</i>	acquisition
acrobat reader	www.gurulabs.com/downloads	documentation viewer

9.2.1.1 Additional packages for spectrometer control - bootparamd

The software packages *bootparamd* is necessary if the PC controls a spectrometer.


1. Check if *bootparamd* is installed so far
 - `rpm -q bootparamd`
 - is the output like this
bootparamd-<version-number>
 - Yes? You can skip this chapter
 - No? Go to step 2.
2. Insert DVD *TopSpin 1.3*
3. Click on the DVD icon on the desktop
4. A window `file:/mnt/cdrom - Konqueror` will appear
 - click **Linux** → **rpm**
5. Click on `bootparamd-<version_number>.i386.rpm`

6. A window *Query* will appear



- enter root password

7. A window *Completed System Preparation* will appear



- click **Continue**

9.2.1.2 Additional packages for documentation view - Acrobat Reader

The software packages *acroread* is necessary if you wants to view manuals in the *Help* menu of TOPSPIN.

1. Check if *acroread* is installed so far

- Open a shell and type in:

```
rpm -q acroread
```
- is the output like this

```
acroread-<version-number>
```
- Yes? You can skip this chapter
- No? Go to step 2.

2. Open Mozilla and go to <http://www.gurulabs.com/downloads.html>

3. Search for the entry `acroread-<version-number>.i386.rpm`
4. Download `acroread-<version-number>.i386.rpm`
5. Download `acroread-plugin-<version-number>.i386.rpm`
6. Install the rpm packages. you can do this either...
 - on a shell, type in the command


```
rpm -Uvh <path_of_the_rpm>/acroread-<version-number>.i386.rpm
```
 - click on the icon *Home* on the desktop. Go to the directory the rpm packages were stored and click on the file. Installation will start immediately or you will be asked for the root password.

9.3 Installation of additional Red Hat packages on already installed systems

If you have done the full installation of RED HAT ENTERPRISE LINUX WS 3 as it is described in chapter 9.2, the packages described in this chapter are already installed and you can skip this chapter.

If you want to prepare an already installed RED HAT ENTERPRISE LINUX WS 3 system for TOPSPIN, you have to make sure that the following packages are available:

Package	on CD	Necessary for
nfs-utils	RED HAT ENTERPRISE LINUX WS 3 <i>Operating System CD 2 of 4</i>	acquisition
redhat-config-nfs	RED HAT ENTERPRISE LINUX WS 3 <i>Operating System CD 2 of 4</i>	acquisition
uucp	RED HAT ENTERPRISE LINUX WS 3 <i>Operating System CD 2 of 4</i>	acquisition
glibc-kernheaders	RED HAT ENTERPRISE LINUX WS 3 <i>Operating System CD 3 of 4</i>	compiling AU programs
glibc-devel	RED HAT ENTERPRISE LINUX WS 3 <i>Operating System CD 3 of 4</i>	compiling AU programs
mutt	RED HAT ENTERPRISE LINUX WS 3 <i>Operating System CD 2 of 4</i>	email

For an installation description see the following chapters. For more information

about useful packages, please refer to chapter 15.3.1

9.3.0.3 uucp, nfs-utils and redhat-config-nfs

The software packages *nfs-utils*, *redhat-config-nfs* and *uucp* are necessary if the PC controls a spectrometer.

1. Check if *nfs-utils*, *redhat-config-nfs* and *uucp* is installed so far

- Open a shell and type in:

```
rpm -q nfs-utils
```

```
rpm -q redhat-config-nfs
```

```
rpm -q uucp
```

- is the output like this

```
nfs-utils-<version-number>
```

```
redhat-config-nfs-<version-number>
```


```
uucp-<version-number>
```

- Yes? You can skip this chapter
- No? Go to step 2.

2. Insert RED HAT ENTERPRISE LINUX WS 3 *Operating System CD 2 of 4*

3. Click <red hat> → **System Tools** → **More System Tools** → **File Manager - Super User Mode**

4. A window *run as root - KDE su* will appear



- enter root password

5. A window *file:/root - Konqueror* will appear

- go to **file:/mnt/cdrom/RedHat/RPMS**

6. Search for the entry *uucp*-<version-number>.i386.rpm

7. Double click **uucp<version-number>.i386.rpm**
8. A window *Completed System Preparation* will appear
 - click **Continue**
9. In the window *file:/root - Konqueror* search for the entry *redhat-config-nfs-<version-number>.i386.rpm*
10. Double click **redhat-config-nfs-<version-number>.i386.rpm**
11. A window *Completed System Preparation* will appear
 - click **Continue**
12. In the window *file:/root - Konqueror* search for the entry *nfs-utils-<version-number>.i386.rpm*
13. Double click **nfs-utils-<version-number>.i386.rpm**
14. A window *Completed System Preparation* will appear
 - click **Continue**
15. Close the *Konqueror* window *file:/mnt/cdrom/RedHat/RPMS*
16. To eject the CD click right on the CD icon on the desktop and choose **eject**

9.3.1 Additional packages for the compilation of AU programs

The packages *glibc-kernheaders* and *glibc-devel* are necessary for the compilation of AU programs.

9.3.1.1 glibc-kernheaders

1. Check if *glibc-kernheaders* is installed so far
 - **rpm -q glibc-kernheaders**
 - is the output like this
glibc-kernheaders-<version-number>
 - Yes? You can skip this chapter
 - No? Go to step 2.
2. Insert RED HAT ENTERPRISE LINUX WS 3 *Operating System CD 2 of 3*
3. Click **<red hat>** → **System Tools** → **More System Tools** → **File Manager - Super User Mode**
4. A window *run as root - KDE su* will appear

- enter root password
5. A window *file:/root - Konqueror* will appear
 - go to **file:/mnt/cdrom/RedHat/RPMS**
 6. Search for the entry *glibc-kernheaders-<version-number>.i386.rpm*
 7. Double click **glibc-kernheaders<version-number>.i386.rpm**
 8. A window *Completed System Preparation* will appear
 - click **Continue**
 9. Close the *Konqueror* window *file:/mnt/cdrom/RedHat/RPMS*
 10. To eject the CD click right on the CD icon on the desktop and choose **eject**

9.3.1.2 glibc-devel

1. Check if *glibc-devel* is installed so far
 - **rpm -q glibc-devel**
 - is the output like this
glibc-devel-<version-number>
 - Yes? You can skip this chapter
 - No? Go to step 2.
2. Insert RED HAT ENTERPRISE LINUX WS 3 *Operating System CD 2 of 3*
3. Click **<red hat>** → **System Tools** → **More System Tools** → **File Manager - Super User Mode**
4. A window *run as root - KDE su* will appear
 - enter root password
5. A window *file:/root - Konqueror* will appear
 - go to **file:/mnt/cdrom/RedHat/RPMS**
6. Search for the entry *glibc-devel-<version-number>.i386.rpm*
7. Double click **glibc-devel-<version-number>.i386.rpm**
8. A window *Completed System Preparation* will appear
 - click **Continue**
9. Close the *Konqueror* window *file:/mnt/cdrom/RedHat/RPMS*
10. To eject the CD click right on the CD icon on the desktop and choose **eject**

9.3.2 Additional package for sending data by email

The package *mutt* is necessary for the *smail* command in *TOPSPIN*, the email notification in *ICON-NMR* and for the automatic notification for helium refill. If you want to work with these tools, you have to install the *mutt* package.

9.3.2.1 Mutt

1. Check if *mutt* is installed so far
 - `rpm -q mutt`
 - is the output like this
`mutt-<version-number>`
 - Yes? You can skip this chapter
 - No? Go to step 2.
2. Insert RED HAT ENTERPRISE LINUX WS 3 *Operating System CD 2 of 3*
3. Click **<red hat>** → **System Tools** → **More System Tools** → **File Manager - Super User Mode**
4. A window *run as root - KDE su* will appear
 - enter root password
5. A window *file:/root - Konqueror* will appear
 - go to **file:/mnt/cdrom/RedHat/RPMS**
6. Search for the entry `mutt-<version-number>.i386.rpm`
7. Double click **mutt-<version-number>.i386.rpm**
8. A window *Completed System Preparation* will appear
 - click **Continue**
9. Close the *Konqueror* window *file:/mnt/cdrom/RedHat/RPMS*
10. To eject the CD click right on the CD icon on the desktop and choose **eject**

9.4 Check if all necessary components are installed and active

To check if all necessary components are installed and active, please type

```
/sbin/chkconfig --list
```


in a shell and compare the output with this example

```
nfs          0:off 1:off 2:on 3:on 4:on 5:on 6:off
nfslock      0:off 1:off 2:off 3:on 4:on 5:on 6:off
bootparamd   0:off 1:off 2:on 3:on 4:on 5:on 6:off
portmap      0:off 1:off 2:off 3:on 4:on 5:on 6:off
```

If an entry is off, you can activate it in general:

```
/sbin/chkconfig <service> on
```

If an entry is off, you can activate it for a specific runlevel:

```
/sbin/chkconfig --level <runlevel> <service> on
```

for example to start bootparamd at runlevel 5:

```
/sbin/chkconfig --level 5 bootparamd on
```

With **chkconfig** command you define (permanently) which service will be started/not started on which runlevel after the next reboot. It does NOT activate/deactivate the respective service, instantly.

For a temporary activation/deactivation of a process you have to use the command:

```
service <service> start/stop
```

There are mostly more possibilities for the different services like reload, restart To get a list which possibilities are available just type in

```
service <service>
```

9.5 Dual Boot computer

It is possible to configure a PC in a dual boot mode, so that you can decide during boot time if you like to start Windows or Linux. This is not a necessary feature for spectrometer PC's and Bruker does not recommend to configure the PC in dual boot mode.

If you like to configure a dual boot PC on your own, Bruker has some tips in the FAQ list on the Bruker web page:

<http://www.bruker-biospin.de/NMR/nmrsoftw/passwd/faq/faqalone/161.html>

Chapter 10

Network configuration

10.1 Configuring the laboratory network

The easiest way to add a Linux PC to a laboratory network is using 'DHCP' (**D**ynamic **H**ost **C**onfiguration **P**rotocol). It allows you to go online just by defining the hostname. The IP address, name- and mailserver will be created dynamically. Please contact your network administrator to find out whether DHCP is the right choice for you, or whether a fixed TCP/IP address is to be assigned to your PC. Because the licensing of TopSpin needs a correct network configuration there are three ways to do this:

1. Those PC's which are configured by Bruker BioSpin Germany are set to DHCP. If this configuration can be used the PC should be connected to the network before the boot process starts. After the PC is coming up the network works correctly and TopSpin can be started.
2. If the PC should be connected to the network by using a fixed IP address you have to ask your system administrator for it. After the correct configuration the network works correctly and TopSpin can be started.
3. If the PC should not get a network connection it is necessary to work with a neutral IP address e. g. 192.168.1.1. This address is not routed to any network computer by world-wide-agreement. If you configure this number you will not have a connection to any other computer (internet, intranet, network printers ...)

but TopSpin can be started

10.2 Some useful commands

- There are two commands that allows to easily check if network is set up correctly:

```
hostname
```

```
hostname -f
```

hostname has to show the short name of the PC, e. g. *nmrpc* and **hostname -f** has to show the full qualified name e. g. *nmrpc.chemistry.university.anywhere*. If the answers are not as described, you find out that something is wrong with your network configuration. Because the Linux system reacts very sensitive on wrong network configuration, you might get several problems, e. g. TopSpin does not start, spectrometer will not boot, the TopSpin license can not checked out ...

- A convenient command for checking the functionality of the ethernet cards is

```
/sbin/mii-tool
```

The output of this command shows a list of all ethernet cards and gives the information if the connection is working properly or not.

- If a ethernet card is not able to be activated in the graphical tool, you can use the following commands (in this example for ,eth0'):

```
service network restart
```

```
ifdown eth0
```

```
ifup eth0
```

10.2.1 Configuring the laboratory network by DHCP

1. Click <red hat> → **System Settings** → **Network** ¹


a window *Query* will appear and asks for the root password

- enter root password

1. The shell command for starting the network configuration is:


redhat-config-network or simply **neat**

2. A window *Network Configuration* will appear



- On the tab pane **Devices** click on the ethernet card you want to use for the laboratory network (typically this 'eth0') and click **Edit**

3. A window *Ethernet Device* will appear



- select the field:

Activate device when computer starts

- select the field:

Automatically obtain IP address settings with:


and click the little arrow to the right of this field and select:

dhcp

!

- enter the hostname of the PC in the field *Hostname (optional)*
- Contact your network administrator to get the information if you can check the field *Automatically obtain DNS information from provider*
- click **OK**


e) In the window *Network Configuration* click on the tab pane *Hosts*¹



- Either this entry is correct and was visible in previous Linux versions, you will not see here an entry *127.0.0.1 / Name: localhost.localdomain / Nick-name: localhost* The operating system controls this entry and does not show it here in this interface. Nevertheless you can check the file */etc/hosts* to get the information if this entry is there or not.
- Be sure that there is **NO** entry with the name of this PC. If there is one highlight it and click on **Delete**

1. This information is stored in the file */etc/hosts*

- The two entries that can be seen in the picture are necessary for the spectrometer connection. If this PC does not control a spectrometer they are not necessary.
- f) Click on the tab pane *DNS*



- enter the hostname without the domainname into the field *Hostname* ¹
- click **File** → **Save**
- click **File** → **Quit**


10.2.2 Configuring the laboratory network by using a fixed TCP/IP address

1. Click <red hat> → **System Settings** → **Network** ²
a window *Query* will appear and asks for the root password
 - enter root password

1. In this field you have to enter the hostname only, without domainnames, i. e. only *nmrhc* and not *nmrhc.chemistry.university.anywhere*


2. The shell command for starting the network configuration is:
redhat-config-network or simply **neat**

2. A window *Network Configuration* will appear




- On the tab pane **Devices** click on the ethernet card you want to use for the laboratory network (typically this 'eth0') and click **Edit**

3. A window *Ethernet Device* will appear



- select the field:
Activate device when computer starts
- select the field:
Statically set IP addresses
and type in the
IP address,
Subnet mask and
Default Gateway Address
that should be used from this PC (if you do not know them, contact your network administrator)
- click **OK**

g) In the window *Network Configuration* click on the tab pane *Hosts*¹




- Either this entry is correct and was visible in previous Linux versions, you will not see here an entry *127.0.0.1 / Name: localhost.localdomain / Nick-name: localhost*. The operating system controls this entry and does not show it here in this interface. Nevertheless you can check the file */etc/hosts* to get the information if this entry is there or not.

1. This information is stored in the file */etc/hosts*

- click **Add** and create an entry for the IP address of this PC. The Name has to be the full qualified domain name like this:
nmrpc.chemistry.university.anywhere
and the Nickname:
nmrpc
- The two entries that can be seen in the picture are necessary for the spectrometer connection. If this PC does not control a spectrometer they are not necessary.

h) Click on the tab pane *DNS*



- enter the hostname without the domainname into the field *Hostname*¹
- enter the domain name in the field *DNS Search Path*
- click **File** → **Save**
- click **File** → **Quit**

Remember, that the network administrator must enter your new host on the name-server.

1. In this field you have to enter the hostname only, without domainnames, i. e. only *nmrpc* and not *nmrpc.chemistry.university.anywhere*

10.2.3 Spectrometer PC that should not be connected to the laboratory network


In case the spectrometer PC should not be connected to the laboratory network, it is nevertheless necessary to configure a kind of 'virtual network'.

1. Click <red hat> → **System Settings** → **Network** ¹

a window *Query* will appear and asks for the root password

- enter root password


2. A window *Network Configuration* will appear



- On the tab pane **Devices** click on the ethernet card you want to use for the laboratory network (typically this 'eth0') and click **Edit**


1. The shell command for starting the network configuration is:
redhat-config-network or simply **neat**

3. A window *Ethernet Device* will appear



- select the field:
Activate device when computer starts
- select the field:
Statically set IP addresses
and enter in the field:
Address:
192.168.1.1
and in the field:
Subnet Mask:
255.255.255.0
- click **OK**


- i) In the window *Network Configuration* click on the tab pane *Hosts*¹



- Either this entry is correct and was visible in previous Linux versions, you will not see here an entry *127.0.0.1 / Name: localhost.localdomain / Nickname: localhost*. The operating system controls this entry and does not show it here in this interface. Nevertheless you can check the file */etc/hosts* to get the information if this entry is there or not.
- click **Add** and create an entry for the IP address of this PC. The Name has to be the full qualified domain name like this:
nmrpc.chemistry.university.anywhere
and the Nickname:
nmrpc
- The two entries that can be seen in the picture are necessary for the spectrometer connection. If this PC does not control a spectrometer they are not necessary.

1. This information is stored in the file */etc/hosts*

j) Click on the tab pane *DNS*



- enter the hostname without the domainname into the field *Hostname*¹
- enter the domain name in the field *DNS Search Path*
- leave the DNS fields empty
- click **File** → **Save**
- click **File** → **Quit**

10.3 Configuring the spectrometer network



The network configuration for spectrometer control will be done automatically in case the ethernet card is not yet configured. This configuration will be done by the *Diskless* package that is part of a TopSpin installation. The description of the following chapter can be used for checking the configuration.

10.3.1 Configuration of the spectrometer network


1. Click <red hat> → **System Settings** → **Network**²

1. In this field you have to enter the hostname only, without domainnames, i. e. only *nmrhc* and not *nmrhc.chemistry.university.anywhere*

a window *Query* will appear and asks for the root password

- enter root password


2. A window *Network Configuration* will appear



- click on the ethernet card you want to use for the spectrometer network (normally this 'eth1') and click **Edit**

2. The shell command for starting the network configuration is:
redhat-config-network or simply **neat**

3. A window *Ethernet Device* will appear



- be sure that the field *Activate Interface at boot time* is checked¹
- be sure that the field *Allow any user to (de)activate interface* is not checked
- select the field:

Statically set IP address

and enter in the field:

Address:

149.236.99.1

and in the field:

Subnet Mask:

255.255.255.0

- click **OK**

k) In the window *Network Configuration* click on the tab pane *Hosts*²

1. On notebooks it might be useful to disable the 'Activate Interface at boot time' function. This will avoid waiting for time outs during boot process, if there is no network connection.

- be sure that there are the following entries:

IP	Name	Nickname
149.236.99.1	ASP_ST2	<whatever>
149.236.99.99	spect	<whatever>

If one of them is missing click **Add** and create them

- click **File** → **Save**
- click **File** → **Quit**

10.3.2 Changing the hostname

For changing the hostname please use the network configuration tool as described in chapter 10.1.

Changing the hostname may have an effect on programs. For the Bruker BioSpin software you have to check these two places:

1. In case the PC is a license server the name of the PC is stored in the license.dat file that is stored in the flexlm installation directory. In a default installation this is:

```
/usr/local/flexlm/Bruker/licenses/license.dat
```

If this PC acts as a license server and you changed the license file, you have to be sure that the respective file is changed in the same way on all license clients.

2. The user specific configuration is stored in the directory:

```
<USER_HOME>/ .topspin-<PC_NAME>
```

To keep the individual TopSpin configuration in case the computer name was changed, rename these directories.

10.4 Email client

For using the *TOPSPIN* command *TOPSPIN smail* under Linux, it is necessary to install the program *MUTT*. *MUTT* is distributed by Red Hat and most Linux sys-

2. This information is stored in the file */etc/hosts*

tems. For installation please refer to chapter 9.3.2

10.5 Network diagnostic commands

Here is a list of the most important network diagnostic commands as they can be entered from the shell:

- **arp** - displays or modifies the IP to physical address translation tables
- **hostname** - displays the short name of the current host
- **hostname -f** - displays the full qualified name of the current host
- **host IP** - displays the name that is connected with this IP
- **ifconfig** - displays information about the current TCP/IP configuration
- **netstat** - displays protocol statistics and current TCP/IP connections
- **ping** - checks if a destination host is receiving TCP/IP packets
- **route** - maintains and displays routing tables
- **traceroute** - displays the route, packets take to a destination host
- **mii-tool** - list all ethernet cards and show if they work or not

For more details about commands you can enter either:

```
man <command>
```

or

```
<command> --help
```

10.6 Accessing remote files and directories


10.6.1 Sharing directories from a Linux host using Samba

If you are logged in on a Windows 2000/XP host, you can access data which reside on a Linux host, using the sharing mechanism. This, however, only works if the Samba package is installed on the Linux host (see chapter 15.3).

After installation of Samba is done, you can configure it with SWAT (Samba Web Administration). This tool can be started from your Web Browser, e.g.:

```
<red hat> → Internet → Mozilla Web Browser
```

by entering the address: *http://localhost:901/*



You will be asked for username and password. For configuration of SWAT it is necessary to login as root. Other users are neither able to see the whole information nor can change it!

There you have to define:

- the **password** and the **user** for configuration in the **Password** section
- the **workgroup** in the **GLOBALS** section
- the **directories** you want to share in the **SHARES** section

Reboot the PC

10.6.2 Mounting Windows 2000/XP partitions on a Linux host using Samba

This is more convenient than mounting via NFS because you do not need to have a NFS package installed on the Windows PC. The only thing that you need is the Samba package for the Linux PC, which is available on the Red Hat CDs.

To define a mounted directory you have to gather the following information:

1. Check a valid user and the name of the Windows PC and get the information which directories are available (=shared), enter the command:

```
smbclient -L <WIN_PVC_NAME> -U <WIN_PC_USER>
```

for example for the Windows PC *nmrpc.chemistry.university.anywhere* and the user *Administrator* the command looks like this:

```
smbclient -L nmrpc.chemistry.university.anywhere -U Administrator
```

You will be asked to enter the users password. The answer of the command is a list of all shared directories on the Windows PC, e. g.:

Sharename	Type	Coment
-----	-----	-----
E\$	Disk	Default share
tmp	Disk	

now you have four information available

- a) <WIN_PC_NAME>
- b) <WIN_PC_USER_NAME>
- c) <WIN_PC_USER_NAME_PASSWORD>
- d) <WIN_PC_SHARE>

2. Create a group called *nmruser* and add those user who should get access into the shared directory. To get information how to create a group please refer to chapter 8.1
3. As a second step you have to define the directory wherein the contents of the shared Windows directory should appear. It is recommended to create therefor a new directory:

```
mkdir <WIN_PC_DATA_DIR>
```

4. To allow an automatic mount of the directory of the Windows PC during boot of the Linux PC you have to create a text file (= <FILEX>) with the information of the user name and corresponding password. The syntax of the file has to be:

```
username = <WIN_PC_USER_NAME>
password = <WIN_PC_USER_NAME_PASSWORD>
```

Because the file contains a plain password, you have to make sure that it is only readable by root. You can do this as root on a LINUX shell with the two commands:

```
chown root.root <FILEX>
```

```
chmod 600 <FILEX>
```

Now you can put together the information and add a new line in the file

```
/etc/fstab:
//<WIN_PC_NAME>/<WIN_PC_SHARE> /<WIN_PC_DATA_DIR> smbfs
rw,gid=nmruser,dmask=775,credentials=<FILEEX>,user 1 2
```

This is one line between `smbfs` and `rw` is a space. The option `gid=nmruser,dmask=775` gives full control to the user and the group, others get read and execute.

10.6.3 Mounting Windows 2000/XP partitions from a Linux host using NFS

This requires that the NFS server is installed on the Windows PC. If this PC is controlling a spectrometer, the Hummingbird NFS Server is already installed. The partition `C:\Bruker\Diskless\clients\spect` is mounted via NFS by the spectrometer CCU (which is a UNIX host). In the same way, you can mount any Windows 2000/XP partition or directory on any UNIX or Linux host. Proceed as follows:

1. Export the Windows 2000/XP partition as described in *Installation Guide for Windows XP*.
2. On the Linux host:
 - a) Open a shell
 - b) Type `su` to become superuser
 - c) Create a mount directory, e.g. `mkdir /ntdata`
 - d) mount the partition, e.g.:

```
mount nthost:/c/data /ntdata
```

Note that if the pathname on the Windows 2000/XP host contains upper case letters, these must also be specified with the mount command.

10.6.4 FTP and Telnet

The security level of `ftp` and `telnet` are not very high. Because of this, these services are not installed on a typical RED HAT ENTERPRISE LINUX WS 3 installation. The modern and secure alternative is SSH.

If you want to install `ftp` and `telnet` services see chapter 15.3

Chapter 11

Installing a Printer

11.1 Introduction

In the following you can find the description of the print system ,CUPS'. In RED HAT ENTERPRISE LINUX WS 3 CUPS is the only print system, the older ,LPRng' is not available anymore.

Note: While using CUPS, TOPSPIN PLOT EDITOR 1.3 can print the Bruker logo which is part of the standard Bruker layouts. BUT previous XWIN-PLOT versions are not able to create a plot with the Bruker logo using CUPS. The workaround here is to remove one of the last lines in the file:

```
<TOPSPIN_HOME>/plot/epsi/bruker.eps
```

Find the following lines at the end of the file:

```
...
%%Trailer
cleartomark
countdictstack exch sub { end } repeat
restore
%%EOF
```

and remove the line ,%%Trailer' so that it looks like this:

```
...
cleartomark
```

```
countdictstack exch sub { end } repeat
restore
%%EOF
```

Another workaround is, to use the additional option `,-oraw'` in the print command of the TopSpin plot editor layout. This is only possible in case the required printer is able to understand PostScript code.

11.2 Print system ,CUPS‘ 1

11.2.1 Installation of a Local printer


In order to install a Local printer on a Linux PC, you have to take the following steps:

1. Click <red hat> → **System Settings** → **Printing**

A window *Query* will appear and asks for the root password

- enter root password

2. The window *Printer configuration* will appear



- Click **New**


1. To define the default printer of the system root can also execute the command:

`lpadmin -d <printer_name>`

To define a user specific default printer the user has to execute the command:


`lptions -d <printer_name>`

3. A window *Add a new Print Queue* will appear



- Click **Forward**


4. The window changes to *Queue name*



Enter name and short description of the printer (in this example *printer-nmr3* and *HP printer on NMR-3*)

- Click **Forward**


5. The window changes to *Queue type*



Select *Locally-connected* and the port where the printer is connected (in this example */dev/lp0*)

- Click **Forward**


6. The window changes to *Printer model*



Select the manufacturer and the model of the printer in this example *HP Laser Jet 1220*)


- Click **Forward**

7. The window changes to *Finish*, and create the new print queue



- Click **Apply**


8. The new window *Question* appears



In case the printer is already connected it is recommended to print a test page


- Click **Yes**

If it is not connected so far, click on **No** and back in the window *Printer configuration* on **Apply**




9. To configure the printer click on its name and on **Edit**

10. The new window *Edit a print queue* appears



11. To change the driver click on the tab pane **Printer driver**



If there are several printer drivers, 'Postscript' is always a good choice, and the old 'ljet4' works also in most cases very well)

12. Click **OK** and back in the window *Printer configuration* on **Apply**

13. The printer is ready to use

11.2.2 Installation of a Unix network printer


In order to install a Unix network printer on a Linux PC, you have to take the following steps:

1. Click **<red hat>** → **System Settings** → **Printing**

A window *Query* will appear and asks for the root password


- enter root password

2. The window *Printer configuration* will appear




- Click **New**

3. A window *Add a new Print Queue* will appear



- Click **Forward**


4. The window changes to *Queue name*



Enter name and short description of the printer (in this example *printer-nmr3* and *HP printer on NMR-3*)

- Click **Forward**

5. The window changes to *Queue type*




Select *Networked UNIX (LPD)* and enter the hostname/IP address¹ and the queue of the printer (in this example *viola.software.brucker.de* and *laserjet1200*)

- Click **Forward**

1. You can enter either a name or a IP address of the printer. Everything will work what can be used for a successfully **ping** to the respective PC


6. The window changes to *Printer model*



Select the manufacturer and the model of the printer in this example (*HP Laser Jet 1220*)


- Click **Forward**

7. The window changes to *Finish, and create the new print queue*



- Click **Apply**


8. The new window *Question* appears



In case the printer is already connected it is recommended to print a test page


- Click **Yes**

If it is not connected so far, click on **No** and back in the window *Printer configuration* on **Apply**




9. To configure the printer click on its name and on **Edit**

10. The new window *Edit a print queue* appears



11. To change the driver click on the tab pane **Printer driver**



If there are several printer drivers, 'Postscript' is always a good choice, and the old 'ljet4' works also in most cases very well)

12. Click **OK** and back in the window *Printer configuration* on **Apply**

13. The printer is ready to use

11.2.3 Installation of a Windows network printer


In order to install a Windows network printer on a Linux PC, you have to take the following steps:

1. Click <red hat> → **System Settings** → **Printing**

A window *Query* will appear and asks for the root password


- enter root password

2. The window *Printer configuration* will appear




- Click **New**

3. A window *Add a new Print Queue* will appear



- Click **Forward**


4. The window changes to *Queue name*



Enter name and short description of the printer (in this example *printer-nmr3* and *HP printer on NMR-3*)

- Click **Forward**


5. The window changes to *Queue type*



Select *Networked Windows (SMB)* and the required shared printer from the list. You will see all printers that are reachable without username and password. If the required printer is not in the list, ...

- Click **Specify**

The window Authentication appears, enter all entries manually




.. if it is on the list, click on the entry and the window Authentication will appear with those information that are available. Fill in the missing information about the user name and password.

- Click **Forward**

6. The window changes to *Printer model*

Select the manufacturer and the model of the printer




- Click **Forward**

7. The window changes to *Finish, and create the new print queue*

- Click **Apply**


8. The new window *Question* appears



In case the printer is already connected it is recommended to print a test page


- Click **Yes**

If it is not connected so far, click on **No** and back in the window *Printer configuration* on **Apply**




9. To configure the printer click on its name and on **Edit**

10. The new window *Edit a print queue* appears



11. To change the driver click on the tab pane **Printer driver**



If there are several printer drivers, 'Postscript' is always a good choice, and the old 'ljet4' works also in most cases very well)

12. Click **OK** and back in the window *Printer configuration* on **Apply**

13. The printer is ready to use

11.2.4 Installation of a Jet Direct printer


In order to install a Jet Direct printer on a Linux PC, you have to take the following steps:

1. Click **<red hat>** → **System Settings** → **Printing**

A window *Query* will appear and asks for the root password


- enter root password

2. The window *Printer configuration* will appear




- Click **New**

3. A window *Add a new Print Queue* will appear



- Click **Forward**

4. The window changes to *Queue name*




The screenshot shows a Windows dialog box titled "Add a new print queue" with a blue header bar. Below the header, the title "Queue name" is displayed. The main area contains the following text: "Please enter a name for this queue. Choose a short name that begins with a letter and contains no spaces." Below this is a text input field with "printer-nmr3" entered. Underneath is the "About" section: "If you like, you can enter a description of the printer to help you identify it more easily." Below this is another text input field with "HP printer on NMR-3" entered. At the bottom of the dialog are four buttons: "Help" (with a question mark icon), "Cancel" (with an X icon), "Back" (with a left arrow icon), and "Forward" (with a right arrow icon).

Enter name and short description of the printer (in this example *printer-nmr3* and *HP printer on NMR-3*)

- Click **Forward**

5. The window changes to *Queue type*

Select *Networked JetDirect* and enter name and port of the printer (in this example *printersw* and *9100*)




The screenshot shows the same "Add a new print queue" dialog box, but now the title is "Queue type". The main area contains the following text: "Select a queue type:" followed by a dropdown menu showing "Networked_JetDirect". Below this is a table with two columns: "Printer:" and "Port:". The "Printer:" field contains "printersw" and the "Port:" field contains "9100". At the bottom of the dialog are four buttons: "Help" (with a question mark icon), "Cancel" (with an X icon), "Back" (with a left arrow icon), and "Forward" (with a right arrow icon).

- Click **Forward**

6. The window changes to *Printer model*


Select the manufacturer and the model of the printer in this example *HP Laser*

Jet 1220)




- Click **Forward**

7. The window changes to *Finish*, and create the new print queue



- Click **Apply**


8. The new window *Question* appears




In case the printer is already connected it is recommended to print a test page


- Click **Yes**

If it is not connected so far, click on **No** and back in the window *Printer configuration on **Apply***



9. To configure the printer click on its name and on **Edit**
10. The new window *Edit a print queue* appears



11. To change the driver click on the tab pane **Printer driver**

If there are several printer drivers, 'Postscript' is always a good choice, and the old 'ljet4' works also in most cases very well)

- 12. Click **OK** and back in the window *Printer configuration* on **Apply****
- 13. The printer is ready to use**

Part III

Tips, Tricks and Troubleshooting

Chapter 12

Troubleshooting TopSpin

12.1 TopSpin does not start

If *TOPSPIN* does not start, you should first look for an error message in the window where the *TOPSPIN* startup messages appear. Then check if you have one of the following problems:

- A window *Choose Server or File* appears:

The file `license.dat` was not found. You might have one of the following problems:

`license.dat` has been deleted or renamed

→ Setup the file `license.dat` as described in chapter 5.2.

- *Invalid license key (inconsistent license key)*

The file `license.dat` might contain one of the following errors:

- The FEATURE line contains the wrong license key, the wrong date and/or the wrong number of licenses
- The `hostid` is appended at the end of the FEATURE line

- *Invalid host*

- The file `license.dat` contains additional characters at the end of the FEATURE line (after ““““)

12.2 TopSpin troubleshooting

TOPSPIN should not get into any trouble in general! Nevertheless, if your current *TOPSPIN* session has a problem, you can find in this chapter some instructions how to get rid of this situation.

In order to stop TopSpin in case of any problems, the first thing to know is, that TopSpin is working with a client / server structure. The server is the program itself and the client is the graphical interface. It is possible to exit each of them separately. It is also possible to start more than one client (graphical interface). This feature is used for the remote control but can also be very useful for troubleshooting (see below)!

- The client (=graphical user interface of *TOPSPIN*) can be stopped without touching a running acquisition.
- Stopping the program itself will also stop a running acquisition.

12.2.1 Troubleshooting while no acquisition is running

The most convenient way to get rid of a problem inside TOPSPIN is, to stop the complete program. For this you can do the following steps. These steps are increasingly drastic.

1. Enter **`exit`** in the TOPSPIN command line or press the **X** in the upper right corner of the program window.

Is the TOPSPIN window closed?

- Yes? Start TOPSPIN once again, in case there will be still some problems go ahead with step 2
- No? go ahead with step 2

2. Stop the current TOPSPIN window and start a new one

a) Enter in the TOPSPIN command line:

- **`logoff`**

b) Open another Linux shell and start a second TOPSPIN:

- **`topspin -client`**

Has the first TOPSPIN window been closed and does the second one work fine?

- Yes? The problem is solved you can go ahead with TOPSPIN
- No? go ahead with step 3

3. Stop all TOPSPIN windows and end the program itself

a) Try to stop all TOPSPIN windows by enter the command:

- ***logoff***

in the TOPSPIN command line

b) To close all the TOPSPIN windows that could not be closed in step a), open a Linux shell and enter

- **`ps -ef | grep <TOPSPIN_HOME>/jre/bin/java`**

You have to replace the <TOPSPIN_HOME> with the path of your *TOPSPIN* installation. In a default installation this is /opt/topspin, in this case the command would be:

- **`ps -ef | grep /opt/topspin/jre/bin/java`**

The answer is something like this:

```
<USER> <xxx> <yyy> 0 <date/time> pts/<z> 00:00:00 /opt/topspin/jre/bin/java ...
```

c) Kill all those sessions by entering the command in a Linux shell:

- **`kill <xxx>`**

d) In case the ,normal' kill does not work use the kill -9 command:

- **`kill -9 <xxx>`**

e) Enter the command in a Linux shell:

- **`shrm`**

12.2.2 Troubleshooting while an acquisition is running

The difference to the situation that no acquisition is running is, that you can try to stop only the graphical user interface without stopping the acquisition. Steps 1-3 do not stop the acquisition, step 3-5 will stop also the acquisition.

1. Stop the current TOPSPIN window and start a new one

a) Enter in the TOPSPIN command line:

- ***logoff***

b) Open another Linux shell and start a second TOPSPIN:

- **topspin -client**

Has the first TOPSPIN window been closed and does the second one work fine?

- Yes? The problem is solved you can go ahead with TOPSPIN
- No? go ahead with step 2

2. Stop all TOPSPIN windows and end the program itself

a) Try to stop all TOPSPIN windows by enter the command:

- **logoff**

in the command line of the TOPSPIN windows

b) Open a Linux shell and enter

- **ps -ef | grep <TOPSPIN_HOME>/jre/bin/java**

You have to replace the <TOPSPIN_HOME> with the path of your *TOPSPIN* installation. In a default installation this is /opt/topspin, in this case the command would be:

- **ps -ef | grep /opt/topspin/jre/bin/java**

The answer is something like this:

```
<USER> <xxx> <yyy> 0 <date/time> pts/<z> 00:00:00 /opt/topspin/jre/bin/java .....
```

c) Kill all those sessions by entering the command in a Linux shell:

- **kill <xxx>**

d) In case the ,normal' kill does not work use the kill -9 command:

- **kill -9 <xxx>**

Has the first TOPSPIN window been closed and does the second one work fine?

- Yes? The problem is solved you can go ahead with TOPSPIN
- No? go ahead with step 4

3. Enter the command in a Linux shell:

- **shmrn**

4. To get a list of all TOPSPIN processes open a Unix shell and enter

- **uxproc**

5. Kill any hanging *TOPSPIN* processes with the command:

- **kill -9 <xxx>**

and the respective ID of the process you got in step 4.

12.2.3 Send debug information to Bruker BioSpin

If you are in a situation, that TOPSPIN problems appear periodically you can send the information about this situation to Bruker BioSpin. To get the information go into the Linux shell where TOPSPIN was started and press the keys **Ctrl** and **Backslash** simultaneously¹. This will create a debug output inside the shell. Please copy this output into an email and send it to this address:

nmr-software-support@bruker.de

For a complete debugging it is useful to attach the history file of the TOPSPIN session. For this attach the following two files:

```
<TOPSPIN_HOME>/prog/curdir/<USER_NAME>/history
```

```
<USER_HOME>/topspin-<PC_NAME>/prop/protocol.txt
```

12.3 Communication problems between PC and CCU

If the acquisition (**zg**, **gs**, **wobb**) and/or the configuration (**cf**) fail, you might have a communication problem between the Linux PC and the spectrometer CCU. Typically you will get error messages from **TOPSPIN** like:

```
Connect to spect: unable to connect to remote host
Connection to spect. aqport0, failed
Acquisition hardware not present or defect
```

Check the communication between the PC and the CCU in the following way:

- a) Open a shell
- b) Type **telnet spect**

Depending on whether **telnet spect** works (you get a login prompt) or not, continue with paragraph 12.3.1 or 12.3.2.

12.3.1 Communication problems but telnet spect works

If **telnet spect** works, the CCU has booted. This means that the spectrometer network software is installed and running. The communication problems obviously have other causes. Two of the most common causes are mentioned below.

1. If this key combination does not work, you might have a German keyboard configuration, in this case you have to switch it into English.

In each of the following steps you must be logged in on the PC as root and type the indicated commands in a shell.

1. Check if the process *startd* is running on the CCU:

- a) Type **telnet spect** and log in on the CCU as root
- b) Type **ps -efw | grep startd**

Does this command show you 2 processes */etc/startd*?

YES? Then *startd* is running, type **exit**, go to step 2

NO? Then you found the problem, go to step 1c.

- c) Type **/etc/startd**

Trying to start *startd* in this way might show you the reason why *startd* was not started automatically. You might have a hardware problem and starting the *startd* manually indicates which board has a failure. Run the corresponding hardware test (see chapter 12.8.1).

2. Check if the spectrometer has an alias name which is not specified in */etc/hosts*:

- a) Type **cd /<TOPSPIN_HOME>/conf/instr**

<TOPSPIN_HOME> is the directory where *TOPSPIN* is installed (default is */opt/topspin*)

- b) Type **more curinst**

What is the output of this command?

- *spect*? Then the spectrometer does not have an alias name
- different from *spect*, e.g. *av600*? Then do the following:

type **grep spect /etc/hosts**

What is the output of this command?

149.236.99.99 spect av600? Then the alias name is specified

149.236.99.99 spect? Then you found the problem, go to c)

- c) Edit the file */etc/hosts*:

append the alias name, e.g. *av600* after *spect*

If these two steps did not solve your problem, copy the contents of the two history files into an email to your local Bruker office. Where do you can find the files, can be seen by entering the command *hist* in TopSpin. These files contain the neces-

sary information:

```
<TOPSPIN_HOME>/prog/curdir/<USER_NAME>/history  
<USER_HOME>/ .topspin-<PC_NAME>/prop/protocol.txt
```

12.3.2 Communication problems and telnet spect does not work

If **telnet spect** does not work, the CCU has probably not booted up completely. This can have several causes. In each of the following steps you must be logged in on the PC as root and type the indicated commands in a shell. After each step, type **telnet spect** and see if you get a login prompt.

12.3.2.1 Software related problems

1. Open a Linux shell and become superuser; type **su -**
2. Check if there are firewall rules that prohibit to contact the spectrometer

a) go to the directory */etc/sysconfig*

- **cd /etc/sysconfig**

b) what is the configuration of *iptables*¹

- **iptables -L**

Did you get a information that contains three times (*policy ACCEPT*)?

- Yes ? No problem, go to step 3
- No, the answer is at great length? go to step c

c) to check if this is the problem stop the firewall temporarily, type in :

- **service iptables stop**

Check the spectrometer connection once again. Does it work now?

- No ? No problem, go to step 3
- Yes ? You might have found the problem. The system administrator has to make sure that the ethernet card for spectrometer control gets full access or he has to disable firewall package in general. This can be done by the command:

```
chkconfig iptables off
```

Note: You have to reboot the PC.

1. If you use the firewall package *iptables*, you have to make sure that it allows the ethernet card for spectrometer control full access (see chapter 15.2.1)

3. Reboot the PC and CCU using *reviveccu*¹:

type **reviveccu**

→ follow the instructions:

- a) You will be asked to switch off the acquisition rack
- b) *reviveccu* will automatically reboot the PC
- c) Switch on the acquisition rack after the PC has rebooted

4. Check the nameserver configuration:

- a) Type **grep hosts /etc/nsswitch.conf**

What is the output of this command?

Nothing? No problem, go to step 5.

- *Cannot open nsswitch.conf*? No problem, no nameserver; go to step 5.
- *hosts: files dns nisplus*? No problem, this entry is correct; go to step 5
- *hosts dns....* or *hosts nisplus....*? You found the problem: go to step b.

- b) Edit the file */etc/nsswitch.conf* and check if the line *hosts* looks like this: *hosts: files dns nisplus*

where *files* refers to the local */etc/hosts* file and *dns* to the DNS (bind) and *nisplus* refer to nameservers.

This specifies that network commands like **telnet** should first check the local file */etc/hosts* then, if the host was not found, contact the Domain Name System (*dns*) and finally the Network Information System (*nisplus*). The entries *dns* and *nisplus* are only useful if you have the corresponding network systems. Make sure the entry *files* comes first, the existence and the order of entries *nisplus* and *dns* depends on your network.

5. Check if NFS is running:

- a) Type **ps -efw | grep nfsd**

Does the output of this command show you 8 processes */usr/etc/nfsd* ?

YES? Then NFS is running, continue with step 7

NO? Then you found the problem, continue with step c

- b) Check the file

1. The command **reviveccu** executes a shell script which will first remove the two files: */usr/diskless/clients/spect/root/etc/mtab* and *.../utmp* and then reboot spectrometer-CCU and Linux-PC

- `/etc/sysconfig/network`

Does the file contain the entry

`NETWORKING=yes`

YES? Then NFS can be started, continue with step 7

NO? Then you found the problem, change the entry to yes and restart NFS manually. Open a shell and type in:

`service nfs restart`

- c) Type **`/sbin/chkconfig nfs on`**

- d) Type **`reviveccu`**

- e) Run step 5a again and if NFS is still not running continue with step 6

6. Check if NFS is installed:

- a) Type **`rpm -q nfs-utils`**

Does the output of this command show you the package `nfs-utils`?

YES? Then NFS is installed, go to step 7

NO? Then you found the problem, install the package `nfs-utils`. For a installation description please refer to chapter 9.2.1.1

7. Check if the process `bfscd` is running:

- a) Type **`ps -efw | grep bfscd`**

Does the output of this command show you the process

`/usr/diskless/bfscd.linux`

YES? Then `bfscd` is running, continue with step 8

NO? Then you found the problem, continue with step 7b

- b) Check the configuration of `bfscd`. Correct it if it is wrong.

- type in **`/sbin/chkconfig --list bfscd`** and check if the output show `on` in runlevel 5, e. g.¹:

`bfscd` `0:off 1:off 2:on 3:on 4:on 5:on 6:off`

- c) Start the `bfscd` manually.

- type in **`service bfscd start`**

and check if the output of step 7a is correct now. If it is not, go ahead.

1. If a process like this is not started, start it on a shell by typing **`service <process> start`**

- d) Check the output of the command **hostname**

Is the output of this command:

a plain **hostname**? No problem, continue with step 8.

a **hostname.domainname**? You found the problem. Change the hostname as described in chapter 10.3.2. Then continue with step 7h.

- e) Check the output of the command **hostname -f**

Is the output of this command:

a **hostname.domainname**? No problem, continue with step 8.

a plain **hostname**? You found the problem. Change the hostname as described in chapter 10.2. Then continue with step 7h.

- f) Check the file `/etc/hosts`¹ - at least the following entries have to exist for spectrometer control:

```
127.0.0.1      localhost.localdomain    localhost
149.236.99.1  ASP_ST2
149.236.99.99 spect
```

If the ethernet card for the laboratory network is not configured as DHCP, also the name of the own PC has to be listed

```
xyz.xyz.xyz.xyz name_of_own_pc.domain  name_of_own_pc
```

Does this entry exists?

Yes? No problem, continue with step 8.

No? You found the problem. Add the name of the PC as described in chapter 10.2.2. Then continue with step 7h.

- g) Check if the ethernet card for the laboratory network is configured as DHCP (see chapter 10.2.1). If so check if the DHCP server is still active:

- open a Linux shell
- type in **su -** to become superuser
- type in **ifconfig**

you will get a description of the configurations of the ethernet cards. Check the IP address of the PC, this is the number behind *inet addr*: , typ-

1. It is also possible to use the graphical tool for checking / configuring the entry, see chapter 10.2.2

ically of *eth0*

- take the IP address and type in **host <IP-address>**

is the answer like this:

```
<IP-address>.in.addr.arpa. domain name pointer <PC-Name.domainname>
```

Yes? No problem, continue with step 8.

No? You found the problem. Contact your network administrator to check the nameserver and DHCP server.

h) Type **reviveccu**

8. Check if the `/usr/diskless` partitions are exported:

a) Type **/usr/sbin/exportfs -v**

Does the output of this command show you both partitions with correct options :

```
/usr/diskless/dl_usr
spect(ro,async,wdelay,root_squash)
```

and

```
/usr/diskless/clients/spect
spect(rw,async,wdelay,no_root_squash)
```

YES? Then the partitions are exported, continue with step 9

NO? Then you found the problem, continue with step 8b

b) Check the file `/etc/exports` and correct it if it is wrong (see chapter 12.3.3.1)

c) Type **reviveccu**

9. Check if the spectrometer network is correctly configured:

Please refer to chapter 10.3.1

10. Reconfigure the spectrometer network:

- You can use the automatic procedure which will be started by typing:

```
/<TOPSPIN_HOME>/prog/bin/install.net/install.net
```

where `<TOPSPIN_HOME>` is the directory where *TOPSPIN* is installed (default `/opt/topspin`). This command reconfigures the entire spectrometer network.

- You can check it manually, see chapter 10.3

11. Reconfigure the spectrometer network by re-installing the DISKLESS package:
 Insert the **TOPSPIN** DVD into the DVD drive, select the Diskless package and deselect all other packages. For more details see installation description in chapter 10.3.1

12.3.2.2 Hardware related problems

1. Check if the ethernet cable between PC and CCU is properly connected.
 → Remove and re-connect the ethernet cable and the transceiver or minihub.
2. Check if the ethernet interface is broken.
 - a) Physically connect the CCU to the first ethernet interface (*eth0*)¹
 - b) Become superuser; type **su -**
 - c) Change the configurations between eth0 and eth1 as described in chapter 10.2.2 and 10.3.1.

! Note: Write down all configuration parameters like IP addresses, nameserver etc.
3. Replace the ethernet cable, transceiver or minihub between the PC and CCU to find out if either of these components is broken.
4. If a new CCU was installed:
 - a) login as root
 - b) Type **arp -d spect**
 - c) Type **reviveccu**

12.3.3 Files which are needed for the spectrometer network

12.3.3.1 File entries made for the spectrometer network

- /etc/hosts


```
149.236.99.1    ASP_ST2
149.236.99.99  spect
```
- /etc/bootparams


```
spect root=ASP_ST2:/usr/diskless/clients/spect/root
swap  =ASP_ST2:/usr/diskless/clients/spect/swapfile
```

1. assuming that the CCU is currently connected to the second ethernet interface.

```
dump= ASP_ST2:/usr/diskless/clients/spect/dump
```

Note that this is one line with a blank between the three entries.

- /etc/exports
/usr/diskless/dl_usr spect(ro)
/usr/diskless/clients/spect spect(rw,no_root_squash)
- /etc/sysconfig/network
NETWORKING=yes
HOSTNAME=hostttest
GATEWAYDEV=
GATEWAY=
- /etc/sysconfig/network-scripts/ifcfg-lo
DEVICE=lo
IPADDR=127.0.0.1
NETMASK=255.0.0.0
NETWORK=127.0.0.0
If you're having problems with gated making 127.0.0.0/8
a martian, you can change this to something else
(255.255.255.255, for example)
BROADCAST=127.255.255.255
ONBOOT=yes
NAME=loopback
BOOTPROTO=none
- /etc/sysconfig/network-scripts/ifcfg-eth1
BROADCAST=149.236.99.255
NETWORK=149.236.99.0
NETMASK=255.255.255.0
IPADDR=149.236.99.1
DEVICE=eth1
ONBOOT=yes
BOOTPROTO=none
USERCTL=no
- /usr/diskless/clients/spect/root/etc/fstab

```

ASP_ST2:/usr/diskless/clients/spect/root / nfs rw,bg,hard,retry=3,timeo=20 0 0
ASP_ST2:/usr/diskless/dl_usr /usr nfs ro,bg,hard,retry=3,timeo=20 0 0
ASP_ST2:/usr/diskless/clients/spect/var /var nfs rw,bg,hard,retry=3,timeo=20 0 0

```

- /etc/xinetd.d/bootparamd¹

```

service bootparamd
{
type= RPC
socket_type= dgram
user= root
server= /usr/sbin/rpc.bootparamd
wait= yes
protocol= udp
rpc_version= 1
disable= no
instances= 1
log_type= FILE /usr/diskless/bootparamd.log
}

```

12.3.3.2 Files created for the spectrometer network

- /usr/diskless/clients/spect/swapfile

This file can also be created manually with the command:

```
dd if=/dev/zero of=/usr/diskless/clients/spect/swap-file bs=8000000 count=1
```
- /etc/rc2.d/S95bfzd

This file is a symbolic link to the file /etc/init.d/bfzd. From this file the *bfzd* daemon /usr/diskless/bfzd is started.
- /usr/diskless/clients/spect/root/dev

This directory contains a large number of special files such as *tty00*, *tty01*, *tty02*, etc. These files are created during the installation of the DISKLESS package from the *TopSpin* DVD. They can also be created manually:

 - a) Become superuser: type **su -**

1. If the file is modified you can reload it with the shell command: **service xinetd reload**

```
b) cd /usr/diskless/clients/spect/root/dev
c) ./MKDEV -m RS3330
```

12.3.4 View the output of the CCU console

In order to view the output of the CCU you must connect serial port 2¹ of the PC with `ty00` of the CCU, you need a cable with a round rs232 plug on the CCU side (order number HZ04161/A). The program `cu` can be used to make the output of the CCU visible in a shell on the PC.

Check if `cu` is installed:

```
rpm -q uucp
```

if it is installed and active the answer should be:

```
uucp-<versions-no>
```

If it is not installed yet, please refer to the installation description in chapter 9.2.1.1

In case of any problems with this package check the permissions of the program.

1. Type `ls -l /dev/ttyS1`

The permissions, owner, and group of this file must be: `crw-rw-rw-`, `root` and `uucp` respectively. If they are different, change them as follows:

```
chmod a+rw /dev/ttyS1
chown root.uucp /dev/ttyS1
```

1. Type `ls -l /var`

The group of the directory `lock` must be `uucp` and has to have write permissions. If they are different, change them as follows:

- login as root
- type in:

```
chgrp uucp /var/lock
```
- type in:

```
chmod g+rw /var/lock
```

Type `/sbin/reboot`

1. If you connect the cable with serial port 1, you have to use `tyS0` instead of `tyS1`

12.3.4.1 View/change the prom parameters

1. Open a shell and become root

2. Type `cu -l ttyS1`

→ You will get the prompt *connected* and you can see the CCU output messages:

a) Push the reset button of the CCU

b) Press **Control-C** to interrupt the boot process

The CCU will be in Monitor Mode now, the prompt will change to `>>`, type:

```
>> printenv
```

You will get a list of environment variables, the most important are:

```
netaddr=149.236.99.99
```

```
netmask=0xffffffff00
```

```
bootfile=bfs()/usr/diskless/clients/spect/root/unix.r4600_std1
```

```
bootmode=c
```

These values must be shown as specified above². You can set the environment variables as follows:

```
>> setenv netaddr 149.236.99.99
```

```
>> setenv netmask 0xffffffff00
```

To boot the CCU again, type:

```
>> u
```

c) To exit the *cu* program type:

```
>> ~.Enter (press the keys: tilde, dot, Enter)
```

12.3.4.2 Viewing the boot process

When all environment parameters are correct, you can view the boot process and see where it hangs. This might give you an indication as to what the problem is.

1. For CCU 4-7 the bootfile entry is:

```
bootfile=bfs()/usr/diskless/clients/spect/root/unix.r3004_std
```

2. If the *netaddr* is different, the spectrometer network files must be adjusted. However, when you install the DISKLESS package the entries are set back to *149.236.99.99* and your CCU will not boot anymore. Therefore we strongly recommend not to change the *netaddr*.

1. Open a shell and become root
2. Type **cu -l ttys1**
→ you will get the prompt *connected* and after typing **Enter** you can see the CCU output
3. Push the reset button of the CCU

The CCU will try to boot now. The output of the boot process is listed below. The marks ‘---?---’ indicate positions where the boot process possibly hangs or goes into an endless loop:

```
Autoboot: Waiting to load
bfs(/usr/diskless/clients/spect/root/unix.r4600_std
(CTRL-C to abort)
loading
---1---

Obtaining /usr/diskless/clients/spect/root/unix.r4600_std from server c85
851056+115728+824448 entry: 0x80021000
CPU: MIPS R4700 Processor Chip Revision: 1.0
FPU: MIPS R4700 Floating Point Unit [CP1] Revision: 1.0
RISC/os Release 4_52 mips Version UMIPS
Total real memory = 16777216
Available memory = 14356480
---2---

Root on nfs file : Swap on nfs file :
---3---

hostname: spect
domainname: dummy
---4---

Root fstype nfs
Available memory = 12705792
---5---

Root on nfs file : Swap on nfs file :
---6---

The system is coming up. Please wait.
ASP_ST2:/usr/diskless/clients/spect/var mounted on /var
ASP_ST2:/usr/diskless/dl_usr mounted on /usr
```

Internet daemons: portmap inetd.

NFS daemons: nfsd biod lockd statd.

The system is ready.

Booting might hang, stop (and print an error message) or go into an endless loop at one of the positions 1 through 5. This might be caused by one of the problems mentioned below. Check the indicated files and correct them if necessary or perform the indicated steps.

- a) If booting stops at position ---1--- you might have one of the following problems:
- the *bfsd* is not running; go to 12.3.2.1, step 7
 - the DISKLESS packages is not installed; see chapter 12.3.2.1, step 11
 - the network is not set up correctly; check the configuration (see chapter 10.3.1)

If the problems above do not exist, please enter the following command in a shell:

```
/sbin/arp -a
```

The output should contain a line like:

```
spect (149.236.99.99) at a.b.c.d.e.f[ether] on eth1
```

where *a.b.c.d.e.f* is a hexadecimal hardware ethernet address of the CCU

- b) If booting stops at position ---2--- please try the following:
- push the CCU reset button
 - type Ctrl-C to interrupt the boot process and enter:


```
>> setenv bootmode m
>> setenv bootmode c
```
 - push the CCU reset button
- c) If booting stops at position ---3--- one of the following problems exist:
- the CCU (hostname *spect*) is searched for on a nameserver first:
The search order is determined by *hosts*:

```
/etc/nsswitch.conf: hosts: files dns nisplus
```

Make sure the first location to be searched is the local */etc/hosts* file as specified by the entry *files* as specified above.

- the network is not set up correctly; check `/etc/hosts`, `/etc/boot-params`
 - the swapfile does not exist; check the file
`/usr/diskless/clients/spect/swapfile`
If the file does not exist please refer to chapter 12.3.3.2
 - start of `bootparam` fails; check if `bootparamd` is active right now¹:

```
ps -ef | grep bootparamd
```


Does the answer show the two processes `rpc.bootparamd` and `grep bootparamd`?:
Yes? Continue with the troubleshooting list
No? Do the following:

start `bootparamd` right now:

```
service bootparamd start
```


and check what is the general configuration for `bootparamd`:

```
/sbin/chkconfig --list bootparamd
```


If the entry for 5 is off, activate it in general:

```
/sbin/chkconfig --level 5 bootparamd on
```
 - edit the file `/etc/hosts.allow` and add the line
`ALL: spect`
- d) If booting stops at position ---4--- one of the following problems exist:
- NFS is not running; see chapter 12.3.2.1, step 5
 - diskless partitions are not exported; see chapter 12.3.2.1, step 8
 - the network is not set up correctly; check the file `/etc/hosts`
- e) If booting stops at position ---5--- or ---6--- with the following error message please do like described:
- ```
mount root ASP_ST2:/usr/diskless/clients/spect/root
failed, rpc status 15
```

---

1. Note that in case the laboratory network is not active, starting of the `bootparamd` might fail. For checking the network you can use the command `/sbin/mii-tool`. The output shows if the ethernet cards are connected and working properly

```
PANIC: vfs_mountroot: cannot mount root'
you had to check the file (see chapter 12.3.3.1)
/etc/xinetd.d/bootparamd
```

f) If booting is very slow or hangs at an arbitrary point you might have a hardware problem (see chapter 12.3.2.2).

4. To exit the *cu* program type:

```
>> ~.Enter (press the keys: tilde, dot, Enter)
```

### 12.3.4.3 Running CCU diagnostic tests

1. Open a shell and become root

2. Type **cu -l ttys1**

→ you will get the prompt *connected* and you can see the CCU output:

a) Push the reset button of the CCU

b) Press **Control-C** to interrupt the boot process

The CCU will be in monitor mode now, the prompt will change to >>

```
>> setenv bootmode m
```

c) Push the CCU reset button

The CCU diagnostic tests will start automatically. If one or more diagnostic tests fail, the CCU could be broken. If they are all passed successfully, your CCU seems to be okay. The bootmode is automatically set back to its normal value **c**.

d) To exit the *cu* program press:

```
>> ~.Enter (press the keys: tilde, dot, Enter)
```

---

## 12.4 General steps that you can try to locate or solve problems

---

### 12.4.1 Check the WWW for FAQ's, known bugs or known problems

The Bruker WWW pages contain a large amount of information about known errors and problems:

1. Start your WWW browser

2. Go to the Bruker home page: [www.bruker-biospin.de](http://www.bruker-biospin.de)
3. Click **Analytical** → **NMR** → **Software**
4. Click on **FAQ** or **Bugs**

Note that you need an account to enter these pages. If you do not have one, click on **WWW account** and fill out the form to get a free account.

### 12.4.2 RED HAT help

1. Manual Pages: help on Linux commands
  - Click **<red hat>** → **Documentation**

a conqueror window will appear that offers a description of Linux commands

- Type **man <command>**, e.g. **man find**

### 12.4.3 Checking the history of TopSpin

The history function keeps track of all the commands you have entered in *TOPSPIN* and of all error messages. This can be very useful if you discuss a problem with your Bruker service or support person.

- Enter **hist** in *TOPSPIN*

you will get a window with the information about the two history files and the contents of these files. You can use this window to copy-and-paste error messages into your email to Bruker BioSpin.

The history file:

```
<TOPSPIN_HOME>/prog/curdir/<USER_NAME>/history
```

is automatically re-initialized (overwritten) when *TOPSPIN* is started and the previous one is store as history.old

### 12.4.4 Rebooting the PC

Rebooting the PC workstation is a quick and efficient way to kill and restart any hanging processes. Perform the following steps:

1. Check who is logged in; type **who**
  - warn all users who are logged in that you are going to shutdown the system
2. Shutdown the computer; type **/sbin/shutdown**<sup>1</sup>

3. When the computer is down, switch it off; push the power button
4. Switch the computer on, it will boot automatically

### 12.4.5 Rebooting the PC and the CCU with ‘reviveccu’

If you have any kind of communication problem between the PC and the spectrometer CCU, it is always a good idea to reboot both.

1. Check who is logged in; type **who**  
→ warn all users who are logged in that you are going to shutdown the system
2. Become superuser: type **su -**
3. Type **reviveccu**  
→ follow the instructions:
  - a) You will be asked to switch off the acquisition rack
  - b) `reviveccu` will automatically reboot the PC
  - c) Switch on the acquisition rack after the PC has rebooted

---

## 12.5 Various problems in Linux

---

### 12.5.1 The entire screen is frozen

If you have lost mouse control on the entire screen, there are several things you can do. The following steps are increasingly drastic; perform the next higher step only if the previous one didn't help.

1. Press the Escape key: **Esc**
2. Press the three keys: **Ctrl-Alt-F1** simultaneously
  - this will open a text console, you may login as root and enter **init 3** and then **init 5**, if this does not work ...
  - ... press the three keys: **Ctrl-Alt-Delete** simultaneously  
→ this will reboot the system
3. If the system is connected to the network, log in as root via ssh from a different host and type **/sbin/reboot**

- 
1. To force a file system check after the reboot you can use the command **shutdown -F**

→ this will reboot the PC

4. Push the power button, this is not dangerous in case the Linux system uses an ext3 file system

Important: before you reboot or shutdown the system, first check who is logged in if this is still possible: type **who**. Warn all users who are logged in that you are going to shutdown the system.

### 12.5.2 Programs like **shrm**, **reviveccu** and **touser** are not found

These programs are delivered with *TOPSPIN* and are executed by the superuser from a shell. However, if you become superuser with the command **su**, the programs may not be found. You can solve this problem in two different ways:

1. Become superuser with the command '**su -**' or
2. Log out and log in as root

Now the programs **shrm**, **reviveccu**, **touser** etc. can be typed in from any directory and will always be found.

### 12.5.3 Xserver does not start after reboot / not enough disk space available

In case the Xserver does not start after a reboot of the Linux PC, the cause could be that the harddisk is 100% full. In every user's home directory there is a file called

`.xsession-errors`

that could contain a lot of error messages and become very big (xGB!). It is possible to delete this file with the shell command:

```
rm -rf .xsession-errors
```

Note: After deleting this file you have to log out to free the reserved disk space.

It is helpful to check these files and of course the core file which are named as `core.<pid>` under Linux. To find them, you can use the command:

```
find / -name core.* -print
```

---

## 12.6 Trouble under Linux

---

A selection of some important problems that could occur under Linux are described in this chapter.

### 12.6.1 You cannot log in as normal user

There are several possible reasons why you cannot login as normal user. First find out if other normal users can log in. If they can, the problem is user specific:

1. Your password is wrong
2. One of your personal setup files is corrupt. Check the files:
  - /home/<user>/.desktop-<hostname>/log
  - /var/log/messages

---

## 12.7 Useful configuration hints

---

In this chapter you can find some useful configurations of the XF86Config file.

### 1. Definition of the font path

In the file /etc/X11/XF86Config you can add an information about the location of the fontpath, which make it easier for programs to find the necessary fonts:

```
...
Section "Files"
 RgbPath"/usr/X11R6/lib/X11/rgb"
 FontPath "unix:/7100"
EndSection
```

...

2. Definition of the two possibilities USB and PS/2 to connect the mouse. If this configuration is used, you should not use anymore the hardware configuration tool *kudzu* for configuration of the mouse.

```
Section "ServerFlags"
```

...

```
AllowMouseOpenFail
```

```
EndSection
...
Section "InputDevice"
 Identifier"Mouse1"
 Driver"mouse"
 Option"Device""/dev/input/mice"
 Option"Protocol""IMPS/2"
 Option"Emulate3Buttons""off"
 Option"ZAxisMapping""4 5"
EndSection
Section "InputDevice"
 Identifier"DevInputMice"
 Driver"mouse"
 Option "Device" "/dev/psaux"
 Option "Protocol""IMPS/2"
 Option"Emulate3Buttons""off"
 Option"ZAxisMapping""4 5"
EndSection
...
Section "ServerLayout"
 Identifier "AGP"
 Screen "Screen AGP"
 InputDevice "Mouse1" "CorePointer"
 InputDevice "DevInputMice" "SendCoreEvents"
 InputDevice "Keyboard1" "CoreKeyboard"
EndSection
```

## 12.8 Spectrometer hardware tests

---

The *TOPSPIN* DVD contains various programs to test the spectrometer hardware components.

### 12.8.1 Spectrometer internal boards

Test programs for spectrometer internal boards like RCU, TCU, FCU etc. are started on the CCU.

1. Type `telnet spect` and log in as root
2. Type `/u/systest/rcu/rcutest` (to test the rcu)

After starting a test, type `h` to get help or type `auto` to run a complete test. The tests for the other boards are started in the same way, just replace the names accordingly.



---

# Chapter 13

## Important Linux features

---

### 13.1 The Boot Floppy

---

#### 13.1.1 Creating a boot floppy

During Installation process you were asked if you wanted to create a boot floppy. If you want to create a boot floppy after the installation, you have to get the information about the 'operating system release' first. The command `mkbootdisk` uses this information and will create a boot floppy.

1. Open a shell and enter `uname -r`
2. The output of this command looks like `w.x.y-z`, e. g. `2.4.2-2`
3. Enter `mkbootdisk <output_of_uname-r>` , in the example above `mkbootdisk 2.4.2-2`



**Store your boot floppy in a safe place!**

---

## 13.2 Important Commands

---

### 13.2.1 rpm

calls the Red Hat package manager which manage the installation of software packages.

- **rpm -q ABC** → list version of installed package *ABC*
- **rpm --help** → list options and their description for the rpm command
- **rpm -qa** → list all installed programs
- **rpm -Uvh abc** → install or upgrade the rpm package *abc*. This action is silent and shows #-marks as display of progress
- **rpm -Fvh abc** → upgrade the rpm package *abc*. Does not install any new packages. This action is silent and shows #-marks as display of progress
- **rpm -Va** → list all files of the default installation which are modified since the installation
- **rpm -qi ABC** → list informations of version and installation of package *ABC*
- **rpm -ql ABC** → list informations of installed files of package *ABC*

### 13.2.2 cd

**cd:** "change directory"

- **cd** → change to your home directory
- **cd dira** → change to directory *dira*
- **cd /u/data/guest/nmr** → change to directory */u/data/guest/nmr*
- **cd datasetsa** → change to directory */u/data/guest/nmr/datasetsa*
- **cd ..** → change to directory */u/data/guest/nmr*
- **cd ../../usera/nmr** → change to directory */u/data/usera/nmr*

### 13.2.3 ls

**ls**: list the contents of a directory

- **ls** → list files and subdirectories in current directory
- **ls -l** → list complete file information
- **ls -al** → also list files starting with '.'
- **ls -Rl** → also list the contents of all subdirectories
- **ls -lt** or **ls -lrt** → sort by time of last modification

### 13.2.4 chmod

**chmod**: change the permissions of a file or directory

(To understand the permissions and type of access please refer to chapter 14.1)

1. using the character representation: **chmod [ugo][+-][rwx] filename**

- **chmod a+r fila**

make fila readable for all users

```
-r--r--r-- user group 102 Apr 10 12.20 fila
```

- **chmod a+w fila**

make fila writable for all users

```
-rw-rw-rw- user group 102 Apr 10 12.20 fila
```

- **chmod ug+x fila**

make fila executable for the owner and the group

```
-rwxrwxrwx- user group 102 Apr 10 12.20 fila
```

- **chmod o-w fila**

make fila not writable for others

```
-rwxrwxr-- user group 102 Apr 10 12.20 fila
```

- **chmod -R +w dira**

make the directory tree dira, including all files and sub directories, writable for all users

## 2. using the numerical representation: **chmod xxx filename**

- **chmod 444 fila**

make fila readable for all users

```
-r--r--r-- user group 102 Apr 10 12.20 fila
```

- **chmod 222 fila**

make fila writable for all users

```
--w--w--w- user group 102 Apr 10 12.20 fila
```

- **chmod 640 fila**

make fila read/write for the owner, read for the group

```
-rw-r----- user group 102 Apr 10 12.20 fila
```

- **chmod 751 fila**

make fila read/write/execute for the owner, read/execute for the group and execute for others

```
-rwxr-x--x user group 102 Apr 10 12.20 fila
```

- **chmod -R 222 dira**

make the directory tree dira, including all files and sub directories, writable for all users

### 13.2.5 chown / chgrp

**chown /chgrp**: change the owner/group of a file or directory. The output of:

```
ls -l filea or
```

**ls -ld dira** shows the result.

- change the owner of the file *fila* to *usera*

```
chown usera fila
```

```
-r--r--r-- usera system 102 Apr 10 12.20 fila
```

- change the group of the file *fila* to *grpa*

```
chgrp grpa fila
```

```
-r--r--r-- usera grpa 102 Apr 10 12.20 fila
```

- change the owner of directory *dira* to *usera*

```
chown usera dira
```

- change the owner of directory tree *dira* to *usera*, including all files and subdirectories

```
chown -R usera dira
```

- change the group of directory tree *dira* to *grpa*

```
chgrp -R grpa dira
```

- set the owner and group of the directory tree *dira* to *usera* and *grpa*

```
chown usera.grpa dira
```

```
-rwxr-xr-x usera grpa 102 Apr 10 12.20 dira
```

### 13.2.6 **grep**

**grep**: search for a string in a text file and show the lines containing this string

- search for string *stringa* in the file *fila*

```
grep stringa fila
```

- search for *stringa* in all files in the current directory  
**grep stringa \***
- search for *stringa* in all files AND directories in the current directory  
**grep -r stringa \***
- search for *stringa* in all files which starts with *a, b, c* or *d*  
**grep stringa [a-d] \***
- print the lines in *fila* not containing *stringa*  
**grep -v stringa fila**
- ignore upper/lower case in *stringa*  
**grep -i stringa fila**
- search for *stringa* or *stringb* in the file *filea*  
**grep -i "stringa|stringb" filea**
- search for *stringa* in the output of a command, e. g. **ps -ef**  
**ps -efw | grep stringa**

### 13.2.7 find

**find**: search for files and directories in a directory tree

- search for entries (files and directories) called *fila* in the directory tree *dira*  
**find dira -name fila -print**
- search for *fila* on all local partitions, not on NFS mounted partitions

```
find / -local -name fila -print
```

- search for entries that start with *a*

```
find dira -name "a*" -print
```

- search for entries not called *fila*

```
find dira ! -name fila -print
```

- search for entries newer than the file *fila*

```
find dira -newer fila -print
```

- search for files bigger than 1000 blocks

```
find dira -size +1000 -print
```

- run a Linux command (e.g. *rm*, *lp*) on all entries *fila*

```
find dira -name fila -exec lp {} \;
```

for example to delete 2D processed data: *2rr*, *2ri*, *2ir* and *2ii* in */opt/topspin*:

```
find /opt/topspin -name "2[ir][ri]" -exec rm {} \;
```

### 13.2.8 ps

ps: show the Linux process table

- show all processes

```
ps -efw
```

```
UID PID PPID C STIME TTY TIME CMD
root 0 0 0 14:01:12 ? 0:01 sched
root 1 0 0 14:01:12 ? 0:00 /etc/init
root 164 1 0 14:01:48 ? 0:00 /usr/etc/nfsd 4
```

```

root 270 1 0 14:01:56 ? 0:00 /sbin/cron
lp 254 1 0 14:01:54 ? 0:00 /usr/lib/lpsched
root 308 1 0 14:01:58 ? 0:00 /opt/topspin/prog/topspin/bfsd
guest 312 309 0 14:01:59 pts/2 0:00 /opt/topspin/prog/top-
spin/cpr/cpr

```

- show all *TopSPIN* processes

```
ps -efw | grep prog
```

```

UID PID PPID C STIME TTY TIME CMD
guest 1143 686 0 17:07:14pts/1 0:00 /opt/topspin/prog/cpr/cpr
guest 1151 1143 0 17:07:24 pts/1 0:00 /opt/topspin/prog/mod/go
guest 1149 1143 0 17:07:15 pts/1 0:01 /opt/topspin/prog/cpr/xcpu
guest 1170 1168 0 17:09:42 pts/0 0:00 grep prog

```

### 13.2.9 kill

kill: terminate a Linux process

(for a graphic tool use <red hat> → **System** → **KDE System Guard**)

- kill processes with PID 1149, 1151 and 1143

```
kill 1149 1151 1143
```

- kill directly processes with PID 1149, 1151 and 1143

```
kill -9 1149 1151 1143
```

Use **kill -9 PID** only if **kill PID** does not kill the process!

### 13.2.10 compress, gzip

**compress, gzip, bzip**: compression to reduce the size of a file

File compression can be used to:



- save disk space on rarely used files
- reduce network traffic (e.g. **r**cp, **f**tp, E-mail)
  
- compress a file, *fila* is replaced by *fila.Z*  
**compress fila**
  
- uncompress a file, *fila.Z* is replaced by *fila*  
**uncompress fila**
  
- compress a file, *fila* is replaced by *fila.gz*  
**gzip fila**
  
- uncompress a file, *fila.gz* is replaced by *fila*  
**gunzip fila**
  
- compress a file, *fila* is replaced by *fila.gz*  
**bzip2 fila**
  
- uncompress a file, *fila.gz* is replaced by *fila*  
**bunzip2 fila**

Some remarks:

- Uncompression does not require the suffix *.Z* or *.gz*
- compression factor depends on the file type (20-70%)
- **gzip** gives a better reduction in size than **compress**
- many files on the Bruker FTP server have the extension *.gz*, which means they have been compressed with **gzip**

### 13.2.11 tar

**tar**: This command allows to convert a directory (tree) into a file. This is a convenient tool for attaching many files to an email or transferring them via ftp

- copy the directory-tree *dira* into the file *dira.tar*  
**tar cvf dira.tar dira**
- extract *fila* from *dira.tar*  
**tar xvf dira.tar fila**
- copy the directory-tree *dira* into the file *dira.tar* and compresses (gzip) the file  
**tar cvzf dira.tar.gz dira**
- copy the directory-tree *dira* into the file *dira.tar* and compresses (bzip2) the file  
**tar cvjf dira.tar.bz2 dira**
- extract *fila* from compressed file *dira.tar.gz*  
**tar xvzf dira.tar.gz fila**
- list the occurrences of *fila* in *dira.tar*  
**tar tvf dira.tar fila**
- append *filb* to the end of the archive  
**tar rvf dira.tar filb**
- append if *fila* does not exist in *dira.tar* or if it exists but has been modified  
**tar uvf dira.tar fila**
- you can copy with relative pathnames e.g.

```
cd /u/data/guest/nmr
tar cvf dataseta
```

- or with absolute pathnames e.g.

```
tar cvf /u/data/guest/nmr/dataseta
```

---

## 13.3 Important Files

---

### 13.3.1 /etc/shadow

the file /etc/shadow:

- contains the complete password information
- has one line for each user (account) on the system, each line contains 9 fields separated by colons
- example

```
root:1Oz.sVA5j$yMP3jpxgNTIjLT06CZ7QK1:11432:0:99999:7:::
bin:*:11432:0:99999:7:::
daemon:*:11432:0:99999:7:::
adm:*:11432:0:99999:7:::
lp:*:11432:0:99999:7:::
sync:*:11432:0:99999:7:::
shutdown:*:11432:0:99999:7:::
halt:*:11432:0:99999:7:::
mail:*:11432:0:99999:7:::
news:*:11432:0:99999:7:::
uucp:*:11432:0:99999:7:::
operator:*:11432:0:99999:7:::
games:*:11432:0:99999:7:::
gopher:*:11432:0:99999:7:::
ftp:*:11432:0:99999:7:::
nobody:*:11432:0:99999:7:::
```

```
nscd:!!:11432:0:99999:7:::
mailnull:!!:11432:0:99999:7:::
ident:!!:11432:0:99999:7:::
rpc:!!:11432:0:99999:7:::
rpcuser:!!:11432:0:99999:7:::
xfs:!!:11432:0:99999:7:::
gdm:!!:11432:0:99999:7:::
apache:!!:11535:0:99999:7:::
testuser:oaHKvUemEH3TE:11561:0::7:::0
```

meaning of the fields:

1. *Login name*
2. *Encrypted password*
3. *Days since Jan 1, 1970 that password was last changed*
4. *Days before password may be changed*
5. *Days after password must be changed*
6. *Days before password is to expire that user is warned*
7. *Days after password expires that account is disabled*
8. *Days since Jan 1, 1970 that account is disabled*
9. *A reserved field*

### 13.3.2 /etc/passwd

the file `/etc/passwd`:

- contains the complete login account information
- has one line for each user (account) on the system, each line contains 7 fields separated by colons

example

```
root:x:0:0:root:/root:/bin/bash
bin:x:1:1:bin:/bin:
```

```
daemon:x:2:2:daemon:/sbin:
adm:x:3:4:adm:/var/adm:
lp:x:4:7:lp:/var/spool/lpd:
sync:x:5:0:sync:/sbin:/bin/sync
shutdown:x:6:0:shutdown:/sbin:/sbin/shutdown
halt:x:7:0:halt:/sbin:/sbin/halt
mail:x:8:12:mail:/var/spool/mail:
news:x:9:13:news:/var/spool/news:
uucp:x:10:14:uucp:/var/spool/uucp:
operator:x:11:0:operator:/root:
games:x:12:100:games:/usr/games:
gopher:x:13:30:gopher:/usr/lib/gopher-data:
ftp:x:14:50:FTP User:/home/ftp:
nobody:x:99:99:Nobody:/:
nscd:x:28:28:NSCD Daemon:/:bin/false
mailnull:x:47:47:/:var/spool/mqueue:/dev/null
ident:x:98:98:pident user:/:bin/false
rpc:x:32:32:Portmapper RPC user:/:bin/false
rpcuser:x:29:29:RPC Service User:/var/lib/nfs:/bin/false
xfs:x:43:43:X Font Server:/etc/X11/fs:/bin/false
gdm:x:42:42:/:home/gdm:/bin/bash
apache:x:48:48:Apache:/var/www:/bin/false
testuser:x:1001:1003:testuser:/home/testuser:/bin/bash
```

meaning of the fields:

1: *Login name*

2: *password*

3: *User Id*; must be a unique number for each user

4: *Group Id*; must be a unique number for each group

5: *Full Name*; the users full name, may contain blanks

6: *home directory*; usually `/home/"user"`, but can be different

7: *shell*; must contain full pathname e.g. `/bin/sh` or `/bin/bash`

In the field *password* there are several entries possible:

- the field is empty  
Security leak: the user has no password, everybody is able to login on this account
- the field contains an encrypted password  
Security leak: the user has a password, everybody is able to read and probably able to decrypt this password
- the field contains: `*`  
There is no password which translates to the encrypted password: `*`  
Nobody is able to login on this account
- the field contains: `x` (default configuration)  
The account has a shadow password. The encrypted password is written in the file `/etc/shadow` that is only readable for root (see chapter 13.3.1). This is the default setting on Red Hat Enterprise Linux WS 3

### 13.3.3 `/etc/group`

The file `/etc/group`

- contains the complete information for groups of login accounts
- has one line for each group on the system

Each line contains 4 fields separated by colons:

example

```
root::0:root
```

```
users::100:guest, jim
```

meaning of the fields:

1. *the Group Name*
2. *encrypted group password; usually empty*
3. *group Id; must be unique for a group and the same on each host*
4. *the list of users who belong to this group, separated by commas*

Why are users divided into groups?

You can share your files with users within your group and protect them from all other users.

The file `/etc/group` can only be changed by root.

### 13.3.4 `.bashrc`

The file `.bashrc`:

- sets your personal Linux environment
- is executed every time you log in
- resides in your home directory
- is used in connection with the shells `sh` and `bash`

You can edit the file `.bashrc` and add your personal preferences, e.g.:

- `umask 002`
- `PATH=$PATH:/home/guest`
- `NMR=/opt/topspin/data/guest/nmr`
- `export PATH NMR`
- `lr() { find $@ -name lr -ls }`
- `alias ts=topspin`

Now open a new shell:

- *umask* will be set to 002
- guests home directory is part of the Linux search path
- you can type `cd $NMR` instead of  
`cd /opt/topspin/data/guest/nmr`
- you can type `lr <dirA> <dirB>` for searching all files called *lr* in the in the directory trees *dirA* and *dirB*
- you can type `ts` instead of `topspin`

Note the difference between *NMR*, an environment variable and *xw*, a Linux function.

For *csh* and *tcsh* the file `.cshrc` is used. Note that the setting of environment variables is handled differently.

### 13.3.5 umask

umask: set the initial permissions for file creation

the initial permissions of a file are set according to:

$$permission = creation\ permission - umask$$

- the creation permissions are:  
666 for files created by an editor (e.g. vi, emacs)  
777 for an executable program, created by a compiler
- each user can set his own umask, e.g in his `.profile`<sup>1</sup>

examples

1. create a file *fila*, e.g. with vi, for different values of umask

---

1. umask can be set for all users in the file `/etc/profile` and `/etc/csh.login` (*sh* and *bash*) or `/etc/csh.cshrc` (*csh* and *tcsh*)



## 2. show the output of `ls -l fila`

- `umask 000` → permissions 666  
`-rw-rw-rw- 1 owner group 10 Jun 7 12:20 fila`
- `umask 002` → permissions 664  
`-rw-rw-r-- 1 owner group 10 Jun 7 12:20 fila`
- `umask 024` → permissions 642  
`-rw-r---w- 1 owner group 10 Jun 7 12:20 fila`
- `umask 666` → permissions 000  
`----- 1 owner group 10 Jun 7 12:20 fila`

## 13.4 Linux shells

---

A Linux shell:

- is a layer between the operating system and the user.
- offers the user a way to communicate with the OS.
- is selected upon installation of a new user

Type `echo $SHELL` to find out which shell you use!

different Linux shells and the files executed during login:

- `.profile`  
`sh` and `bash` (also executes `.bashrc` if it exists)
- `.cshrc`  
`csh` and `tcsh`

Important differences between shells <sup>1</sup>:

- How do you set environment variables

*sh* and *bash*:

1. **export variable=value**

*csh* and *tcsh*:

1. **setenv variable value**

- *tcsh* and *bash* offer history substitution  
you can use the arrow keys to repeat and modify previously entered commands.

### 13.4.1 How to use the bash

- Open a shell and type in commands
- How to repeat a previously entered command
  1. use the *arrow keys* to get a command back
  2. hit the *Enter key* to execute the command again
- How to repeat and change a previously entered command
  1. use the *arrow keys* to get a command back
  2. use the *arrow keys* to move around in the line  
use **Control-a** to move to the beginning of the line  
use **Control-e** to move to the end of the line
  3. change the command according to your needs
  4. hit the *Enter key* to execute the command
- How to search for previously entered commands:  
type **Control-r** and enter a string or keyword to search for the last command containing that string by typing **Control-r** again, you can search for earlier

- 
1. type **man sh** or **man csh** for more information on these shells

commands containing the same string

- Use <TAB> key for filename/command completion:

Not sure how File or command is written? Just enter the first characters and press <TAB> key, e. g. :

```
vi /etc/X11/XF<TAB> -> vi /etc/X11/XF86Config
```

- How to set the prompt in the bash

The prompt can be changed by defining the environment variable PS1.

1. type **cd** to change to your home directory
2. edit (e.g. with vi) the file `.profile` and set PS1; you can include several items into the prompt, e.g.:
  - a) the hostname of the computer:

```
PS1=' \h'
```
  - b) the logname of the user:

```
PS1=' \u'
```
  - c) the current directory:

```
PS1=' \w'
```
  - d) the last part of the current directory:

```
PS1=' \W'
```
  - e) certain characters, e.g.:

```
PS1=hello or PS1=% or PS1='>' 1
```
  - f) a combinations of the items above, e.g.:

```
PS1=' \h: \w>'
```

in case of the last example the prompt will look like:

```
<hostname> : /u/data/guest >
```

---

1. characters which have a special meaning in Linux, e.g. <, >, \, and ; must be quoted

3. save the file and leave the editor
4. log out and log in again to see the effect of the change

The second prompt can be set by defining the variable `PS2`.

If you want to read the manual pages for the `bash`, type: `man bash`

---

## 13.5 Useful tips for Linux

### 13.5.1 Configure the screen resolution with Xconfigurator

*NOTE:*

It is possible to define the resolutions of the screen depending on the colour depth by the software package *Xconfigurator*. This tool removes all old entries in the respective configuration file

```
/etc/X11/XF86Config-4
```

and stores just those one, you are defining with *Xconfigurator*. This might cause some trouble, so it is strongly recommended to create a backup of the configuration file, e. g.:

```
/etc/X11/XF86Config-4.orig
```

before starting *Xconfigurator* (see below).

Bruker recommends not to use the software package *Xconfigurator*, manually adding the screen resolution entries in the configuration file is more safe (see e. g. description in chapter 12.7).

Starting *Xconfigurator*:

- Press **'Ctrl'-'Alt'-'F1'** to switch to text mode
- Login as root
- Go to the directory `/etc/X11`  
`cd /etc/x11`
- create a backup copy of all files starting with `XF86Config`<sup>1</sup>

```
cp XF86Config XF86Config.orig
cp XF86Config-4 XF86Config-4.orig
```

- Type **Xconfigurator**
- Confirm a couple of questions and choose your monitor
- A window *Probing finished* occurs
- Choose the entry: **Let Me Choose**
- In the window *Select Video Modes*, you have to select all the different monitor settings that you want.
- Confirm the next questions and finish the Xconfigurator
- Back in the text mode type **init 3**
- Type **init 5**

With **Ctrl-Alt+** you can switch between the different monitor settings

---

## 13.6 How to avoid much typing in Linux

---

### 13.6.1 General

1. use environment variables, e.g.:

```
NMR=/opt/topspin/data/usera/nmr
export NMR
```

these 2 lines can be entered in your file `.bashrc`

Open a new shell and type `cd $NMR` to go to the corresponding directory.

2. use Linux functions, e.g.:

```
data()
{
 cd /opt/topspin/data/usera/nmr
}
```

this line can be defined in your file `.bashrc`

- 
1. If you will have problems after the configuration, remove the new files `XF86Config` and `XF86Config-4` and rename the backup copy of the original files to `XF86Config` resp. `XF86Config-4`

Open a new shell and type **data** to execute the function.

3. use copy and paste
  - mark text by holding the left mouse button down
  - paste this text somewhere by clicking the middle mouse button
  
4. use the Linux shells bash or tsh
  - use the arrow keys up/down to repeat and/or modify previously entered commands
  - use the Tab key to complete commands and paths
  
5. use a Linux shell script
  - They are described in the next chapter 13.6.2.

### 13.6.2 Linux shell scripts

A shell script executes (a series of ) Linux commands:

- is an executable text file
- can contain a series of Linux commands
- can contain control loops and variables

example

1. *scripta* displays the date, current directory and contents:

edit a file *scripta* and enter:

```
date
```

```
pwd
```

```
ls
```

make scripta executable: **chmod a+x scripta**

execute the script by typing **./scripta** <sup>1</sup>

When executing scripta a possible output is:

```
Wed Jan 3 03:51:07 PST 1996
/home/guest
scripta
fila
```

Note: the first line of the script should be:

```
#! /bin/sh
```

2. *scriptb* displays the contents of a pulse program:

- edit the file *scriptb* and enter:

```
dir=/opt/topspin/exp/stan/nmr/lists/pp
cd $dir
cat $1
```

\$dir = the value of a parameter which is set in the script

\$1 = an argument given to the script when it is executed

- make the file executable: **chmod +x scriptb**
- execute the script by typing its name and argument: **./scriptb zg30**

3. Example of a control and loop structure in a shell script: <sup>1</sup>

```
for i in `ls *.c` (for all files with the extension .c)
do
 if grep strn $i (if the file contains the string strn)
 then
 rm $i (remove the file)
 fi
done
```

---

1. "." means: the command resides in the current directory. You can skip the "." if the current directory is part of the Linux search path. Type `echo $PATH` to see the Linux search path.

1. For more information on control and loop structures type **man sh** or **man test**.

### 13.6.3 cron - crontab

Any Linux command or program can be executed automatically on a regular basis.

check if crond is running: **ps -efw | grep crond**

if crond is not running: **service crond start**

create a crontab file using the following syntax:

*min hour daymo month daywk command*

*min*: minute (0-59)

*hour*: hour (0-23)

*daymo*: day of the month (1-31)

*month*: month of the year (1-12)

*daywk*: day of the week (0-6 with 0=Sunday)

*command*: the command or program to be executed

Suppose you want to:

- backup your data at 4.30 am on Tuesday to Friday
- delete all core files at 2 am on Sunday

Note that the second task needs root permissions.

- a) type **cd** to go to your home directory
- b) type **su** to become root
- c) edit the file `fila` and enter the following lines:

```
30 04 * * 2-5 tar cv /u/data/<user>/nrm
00 02 * * 0 find / -name core -exec rm {} \;
```
- d) type **crontab fila**

You can skip step b) if all tasks in your crontab file can be done as normal user. Type **man crontab** for more information.



---

# Chapter 14

## File and directory handling

---

### 14.1 Permissions and Type of Access

---

#### 14.1.1 Permissions of a file

The permissions of a file:

- determine what each user can do with the file
- are represented :
  - by a set of characters  $r$ ,  $w$  and  $x$  or
  - by a set of numbers between  $0$  and  $7$

$1 = x$  = file is *executable*

$2 = w$  = file is *writable*

$4 = r$  = file is *readable*

The numbers  $3$ ,  $5$ ,  $6$  and  $7$  are combinations of  $1$ ,  $2$  and  $4$ :

e.g.  $7 = 4 + 2 + 1 = rwx$

### 14.1.2 Permissions of a directory

The permissions of a directory are represented as follows:

*l* = *x* = *search*

*2* = *w* = *write*

*4* = *r* = *read*

### 14.1.3 Types of Access

Each file has 3 sets of permissions for 3 categories of users:

- the actual *owner* (shortened 'u' like *user*)
- users who have the *same group* ID as the owner <sup>1</sup> (shortened 'g' like *group*)
- users who have a *different group* ID than the owner (shortened 'o' like *other*)

### 14.1.4 Example

- Type **ls -l filename** to see the 3 sets of file permission

e.g. **ls -l fila**

```
-rwxrw-r-- 1 jim sys 987 Sep 15 15:08 fila
```

The permissions of the file *fila* are *764*:

*7* = *rw**x* for the *user jim*

*6* = *rw*- for all users of the *group sys*

*4* = *r*- - for *other* users

- Use the command **ls -ld** to see the directory permissions.

e.g. **ls -dl dira**

```
drwxr-xr-- 1 jim sys 987 Sep 15 15:08 dira
```

---

1. the group attribute of a file can be changed with **chgrp**, the permissions of a file always refer to the group ID which is shown with **ls -l filename**

The permissions of the directory *dira* are 754:

- everybody can read (= see) the directory *dira*, e.g. use **ls** on *dira*
- only the owner *jim* and members of the group *sys* can read and search e.g. use **cd**, **pwd**, **ls**, **cat**, **more**, **du**, **file** on or in *dira*
- only the owner *jim* can also write or remove files in *dira*, e.g. use **rm**, **cp**, **mv**, **touch** on files in *dira*
- In the following example the directories *dira*, *dirb* and *dirc* are assumed to be *readable* and *searchable*.

- creating files

If the user *guest* wants to create the file *fila*, e.g. :

```
touch /dira/dirb/dirc/fila
```

the directory *dirc* must be *writable* for *guest*, the parent directories *dirb* and *dira* do not have to be *writable* for *guest*

- removing files

If the user *guest* wants to remove the file *fila*, e.g. :

```
rm /dira/dirb/dirc/fila
```

the directory *dirc* must be *writable* for *guest*, the parent directories *dirb* and *dira* do not have to be *writable* for *guest*. The file *fila* can have any permissions, even 000, and any owner or group, you can always remove it if its directory is *writable* !



---

# Chapter 15

## Useful tools for Red Hat Enterprise Linux WS 3

---

This chapter contains a description of some tools we found useful. We realize that many other tools are available for various purposes and we are very interested in your experiences. Just send your comments to:

`nmr-software-support@bruker.de`

### **15.1 How to use a memory stick under Linux**

---

It is possible to use a USB-memory stick under RED HAT ENTERPRISE LINUX WS 3. In this chapter you can find a description how to activate this functionality.

All USB memory sticks are handled by the scsi disk driver, therefore you must detect which device is attached to the stick. The following steps should give you access to the USB memory:

#### **15.1.1 PC with SCSI drives**

1. Connect the USB stick.
2. Enter the command :

**dmesg**

You will get an answer like the following:

```

hub.c: new USB device 00:1f.4-2, assigned address 17
usb.c: USB device 17 (vend/prod 0xc76/0x5) is not claimed by any
active driver.
Initializing USB Mass Storage driver...
usb.c: registered new driver usb-storage
scsi2 : SCSI emulation for USB Mass Storage devices
 Vendor: Model: Rev:
 Type: Direct-Access ANSI SCSI revision: 02
Attached scsi removable disk sdc at scsi2, channel 0, id 0, lun 0
SCSI device sdc: 256000 512-byte hdwr sectors (131 MB)
sdc: Write Protect is off
 sdc: sdc1
WARNING: USB Mass Storage data integrity not assured
USB Mass Storage device found at 17
USB Mass Storage support registered.

```

The bold line give you a hint for the next step. This system has already two SCSI disks installed sda and sdb and your stick get the next: sdc.

**3. Enter the command:**

**fdisk -l**

The answer is like the following:

```

Disk /dev/sda: 255 heads, 63 sectors, 2213 cylinders
Units =ylinders of 16065 * 512 bytes
Device Boot Start End Blocks Id System
/dev/sda1 * 1 6 48163+ 83 Linux
/dev/sda2 7 2083 16683502+ 83 Linux
/dev/sda3 2084 2213 1044225 82 Linux swap

Disk /dev/sdb: 255 heads, 63 sectors, 4427 cylinders
Units =ylinders of 16065 * 512 bytes
Device Boot Start End Blocks Id System
/dev/sdb1 * 1 523 4200966 83 Linux
/dev/sdb2 524 1046 4200997+ 83 Linux
/dev/sdb3 1047 1569 4200997+ 83 Linux
/dev/sdb4 1570 4427 22956885 83 Linux

```

```
Disk /dev/sdc: 50 heads, 32 sectors, 160 cylinders
Units =cylinders of 1600 * 512 bytes
Device Boot Start End Blocks Id System
/dev/sdc1 * 1 160 127952 6 FAT16
Partition 1 has different physical/logical endings:
phys=(26, 49, 32) logical=(159, 47, 32)
```

The **/dev/sdc1** line is the needed info for the next step.

4. Create a corresponding entry in `/etc/fstab`:


```
/dev/sdc1 /mnt/memstick auto noauto,user 0 0
```

5. Don't forget to create the mount point used in `fstab` (in the example from step 4 this is `/mnt/memstick`) and allow everybody to use that directory. This can be done with the following two commands:

```
mkdir /mnt/memstick
chmod a+rwX /mnt/memstick
```

6. For your convenience create a desktop icon.

- a) Click right on the desktop and choose **Create New** → **Hard Disc**



- b) a window *Properties for HD.desktop* will appear

- c) in section *General*:

- choose a name for the new icon, e. g. `memstick`
- click on the icon symbol, in the upcoming window you can choose one that you enjoy

- d) Select the Device tab and enter the correct path, in this example `/dev/sdc1`

### 15.1.2 PC with IDE drives

1. Connect the USB stick.
2. Enter the command :

**dmesg**

You will see something like this:

```
USB Mass Storage device found at 2
USB Mass Storage support registered.
Attached scsi removable disk sda at scsi1, channel 0, id 0, lun 0
SCSI device sda: 1015808 512-byte hdwr sectors (520 MB)
sda: Write Protect is off
sda: sda1
```

The bold line give you a hint for the next step. The stick on this system get: sda.

3. Enter the command:

**fdisk -l**

The answer is something like this:

```
Disk /dev/sda: 16 heads, 62 sectors, 1024 cylinders
Units = cylinders of 992 * 512 bytes

 Device Boot Start End Blocks Id System
/dev/sda1 * 1 1024 507632 6 FAT16
Partition 1 has different physical/logical beginnings (non-Linux?):
 phys=(0, 1, 1) logical=(0, 0, 33)
Partition 1 has different physical/logical endings:
 phys=(958, 15, 32) logical=(1023, 7, 46)
Partition 1 does not end on cylinder boundary:
 phys=(958, 15, 32) should be (958, 15, 62)

Disk /dev/hda: 255 heads, 63 sectors, 9964 cylinders
Units = cylinders of 16065 * 512 bytes

 Device Boot Start End Blocks Id System
/dev/hda1 * 1 6 48163+ 83 Linux
/dev/hda2 7 2046 16386300 83 Linux
/dev/hda3 2047 2111 522112+ 82 Linux swap
/dev/hda4 2112 9964 63079222+ 83 Linux
```

The **/dev/sda1** line is the needed info for the next step.



4. Create a corresponding entry in `/etc/fstab`:


```
/dev/sda1 /mnt/memstick auto noauto,user 0 0
```

5. Don't forget to create the mount point used in `fstab` (in the example from step 4 this is `/mnt/memstick`) and allow everybody to use that directory. This can be done with the following two commands:

```
mkdir /mnt/memstick
chmod a+rxw /mnt/memstick
```

6. For your convenience create a desktop icon.

- a) Click right on the desktop and choose **Create New** → **Hard Disc**



- b) a window *Properties for HD.desktop* will appear

- c) in section *General*:

- choose a name for the new icon, e. g. `memstick`
- click on the icon symbol, in the upcoming window you can choose one that you enjoy

- d) Select the Device tab and enter the correct path, in this example `/dev/sda1`

## 15.2 Firewall configuration

RED HAT ENTERPRISE LINUX WS 3 has a firewall package inside. You can activate and configure it during installation (see chapter 9.1) or in the running system.



Note that for acquisition the ethernet card for spectrometer control has not to be restricted by the firewall! For this it is strongly recommended to free the ethernet

card after activating the firewall (see below)!

The rules for the firewall, checked and executed during boot time, are stored in the file<sup>1</sup>: `/etc/sysconfig/iptables-config`

The rules in this file are:

- defined during installation (see chapter 9.1)
- added manually (see troubleshooting, chapter 12.3.2.1)
- added manually by using a graphical user interface (see the next chapter)

### 15.2.1 Firewall configuration


RED HAT ENTERPRISE LINUX WS 3 offers a graphical user interface that allows you to define some firewall configurations. The following example describes a firewall configuration that allows *SSH* traffic on *eth0* and frees *eth1* completely.

1. Open a Linux shell
2. Become superuser
  - `su -`
  - enter superuser password
3. Type in
  - `redhat-config-securitylevel`or click on
  - `SystemSettings` → `Security level`


---

1. In RED HAT ENTERPRISE LINUX WS 3 there is also another possibility for the firewall settings: `/etc/sysconfig/iptables`


a window *Security Level Configuration* will appear



4. Select under *Security level* the entry **Enable firewall**
5. Select under *Trusted Services* the entry **SSH**
6. Select under *Trusted devices* the entry **eth1**
7. The window now looks like this, click on **OK**



8. The window *Warning* will appear



Click on **Yes**

Now the configuration of the Firewall is finished and your system is protected right now.

Note that this tool described here allows not a very specific configuration, but it is useful for a quick and secure setting. In case you like to free another service like smb you can use the start command

**redhat-config-securitylevel-tui**

(note that the securitylevel-tui is one word and contains no space)

A window appears wherein you can define the same rules as above described but you get an additional field called 'Other ports'.

Another useful tool is the freeware guarddog that can be downloaded from:

<http://www.simonzone.com/software/guarddog/>

All these three ways that were described in this chapter uses the firewall tool of RED HAT ENTERPRISE LINUX WS 3 *iptables*. With *iptables* it is possible to define a very complex firewall installation.

### 15.3 Tools for internal communication

---

If the Linux NMR workstation should be used as server, the following packages should be installed:

| Package              | on CD                                                           | Necessary for                     |
|----------------------|-----------------------------------------------------------------|-----------------------------------|
| <i>xinetd</i>        | RED HAT ENTERPRISE LINUX WS 3 <i>Operating System CD 2 of 4</i> | "superserver" for several servers |
| <i>telnet-server</i> | RED HAT ENTERPRISE LINUX WS 3 <i>Operating System CD 3 of 4</i> | for telnet logins                 |
| <i>vsftpd</i>        | <i>on TopSpin 1.3 Linux CD in /linux/rpm</i>                    | for ftp sessions                  |
| <i>rsh-server</i>    | RED HAT ENTERPRISE LINUX WS 3 <i>Operating System CD 3 of 4</i> | for rlogin and rsh                |
| <i>samba</i>         | RED HAT ENTERPRISE LINUX WS 3 <i>Operating System CD 2 of 4</i> | to export directories to Windows  |
| <i>samba-swat</i>    | <i>on TopSpin 1.3 Linux DVD in /linux/rpm</i>                   | to configure samba                |

### 15.3.1 Installation description of an additional package: `xinetd`<sup>1</sup>

Other packages are installed similarly

1. Check if `xinetd` is installed so far
  - `rpm -q xinetd`
  - is the output like this  
`xinetd-<version-number>`
  - Yes? You can skip this chapter
  - No? Go to step 2.
2. Insert RED HAT ENTERPRISE LINUX WS 3 *Operating System CD 2 of 4*
3. Click **<red hat>** → **System Tools** → **More System Tools** → **File Manager - Super User Mode**
4. A window *run as root - KDE su* will appear
  - enter root password
5. A window *file:/root - Konqueror* will appear
  - go to **file:/mnt/cdrom/RedHat/RPMS**
6. Search for the entry `xinetd-<version-number>.i386.rpm`
7. Double click `xinetd<version-number>.i386.rpm`
8. A window *Completed System Preparation* will appear
  - click **Continue**
9. Close the *Konqueror* window *file:/mnt/cdrom/RedHat/RPMS*
10. To eject the CD click right on the CD icon on the desktop and choose **eject**

### 15.3.2 Check if all necessary components are installed

To check if all necessary components are installed, please type

```
/sbin/chkconfig --list
```

on a shell and if these lines are exists<sup>2</sup>

---

1. A detailed description of a software installation is described in chapter 9.2.1.1

2. If a process like this is not started, start it on a shell by typing `/sbin/chkconfig <process> on`

```

xinetd 0:off 1:off 2:off 3:on 4:on 5:on 6:off
smb 0:off 1:off 2:off 3:on 4:on 5:on 6:off
xinetd
based services:
rexec: on
rlogin: on
rsh: on
telnet: on
vsftpd: on
swat: on

```

## 15.4 Setting up permissions for rsh, rlogin and rcp

---

The commands **rcp**, **rlogin** and **rsh**<sup>1</sup> require special access permissions. There are two ways to set up these permission.

1. The system administrator can setup up these permission for all users:

- a) Open a shell
- b) Type **su** to become superuser
- c) **cd /etc**
- d) Edit the file `hosts.allow`

Make two entries for every remote host from which access is allowed, e.g.:

```

sunflower
sunflower.netx.lab.com

```

2. As a normal user, you can allow remote access from a remote host for yourself:

- a) Open a shell
- b) Go to your home directory; type **cd**
- c) Edit the file `.rhosts`

Make two entries for every remote host from which access is allowed, e.g.:

```

sunflower
sunflower.netx.lab.com

```

- d) **chmod 644 .rhosts**

---

1. For installation please see chapter 15.3

In both examples *sunflower* is the hostname and *netx.lab.com* is the domainname of the remote host.

To allow **r**cp, **r**sh and **r**login network wide in all directions, you must setup the file `/etc/hosts.allow` or `.rhosts` on each host.

---

## 15.5 Creating a screenshot of the desktop

---

There are several programs that offer the possibility for creating a screenshot of the desktop or of individual windows.

- Click on <red hat> → **Graphics** → **Screen Capture Program**

or the program *Gimp*:

- Open a Linux shell and type in **gimp**  
Inside *gimp* you has to go to **File** → **Acquire** → **Screen Shot ...**

---

## 15.6 Creating an icon on the desktop


---

For a more convenient using of Red Hat, it is recommended to create icons of frequently used commands on the desktop


### 15.6.1 Create an icon for opening a pdf file with Acrobat Reader

1. Click right mouse button on the desktop
  - a popup window appears

- choose **Create New** → **Link To Application**



- a window *Properties for Program.desktop* will appear



- in section *General*:
  - choose a name for the new icon, e. g. Installation Guide
  - click on the icon symbol, in the upcoming window you can choose one that you enjoy
- in section *execute*:
  - type in the command for opening the manual in a shell:

acroread /opt/topspin/prog/docu/english/xwinproc/pdf/install\_guide\_linux.pdf

(if case you have installed TOPSPIN in the default directory /opt/topspin)



## 15.7 Red Hat Network: Hotfixes to keep the system secure

The license agreement of RED HAT ENTERPRISE LINUX WS 3 demanded to register the PC to the Red Hat Network service. A description how to register and how to use the Red Hat Network can be found in the following chapters.


You have to do the activation of the license key first and the registration at Red Hat Network as a second step!

### 15.7.1 Activate your Linux license at Red Hat Network

1. Start mozilla
2. Enter the URL

*http://redhat.com/now*

you will see the following web side



3. If you are not registered so far or you want to use an additional account, click on

- **New Log In**

#### 4. The window changes to

**Activate Your Product**

**Create A New Login**

Take a moment and complete your personal profile. Once registered, you will enjoy special site features, valuable free offers, and faster checkout in our store.

It is Red Hat's policy to ask how you want us to use your contact information and then respect your wishes as to its use. For more information, see our [Privacy Statement](#).

For a detailed list of what is covered under Red Hat Installation Assistance, see our [Service Level Agreement](#).

Fields marked with an asterisk (\*) are required.

**Login Information**

\*Username:  Your username must be alphanumeric.

\*Password:  Your password must be at least 5 alphanumeric characters.

\*Confirm Password:

Always log in automatically  
Note: You should not use automatic log in if you are using a public computer. [Go to it](#)

**Corporate Account**

I am registering using my company information.

Enter all required information (fields marked with a small red asterisk). At the bottom of the window click on

- **Create Login**

#### 5. The window changes to

**Activate Your Product**

All customers with valid, registered Red Hat Enterprise Linux subscriptions are covered under the Intellectual Property Warranty. To register your Enterprise Linux subscription, please follow the steps below. Learn more about [Red Hat's Intellectual Property Warranty](#).

**Red Hat Enterprise Linux Users**

**STEP 1: Confirm Your log in**

Your product will be registered to the following log in:

**User Name:**

If you would like to register this product to a different log in, log out and re-enter your log in information. [log out](#)

**STEP 2: Enter Your Product Information to complete the process**

Activate your product to become entitled to Red Hat services.

**Product ID Number:**  example: XXXXXXXXXXXXXXXXXXXX (slashes optional)


I am proceeding with product activation, I acknowledge that I have read the [Applicable Agreement](#) for my Red Hat product and consent to its terms.

PLEASE NOTE: You will be able to access your Red Hat services in approximately 15 minutes. This delay is required to properly process your subscription.

6. In the field Product ID Number enter your license key for RED HAT ENTER-

## PRISE LINUX WS 3

## 7. The window changes to



8. you have successfully activated your license code, you can now go to next chapter and register the PC to the Red Hat Network


## 15.7.2 Register the PC to the Red Hat Network

1. Log in as root
2. Open a Linux shell and enter
  - Click **up2date**
3. If you are not registered to Red Hat Network so far, the window *Red Hat Network Registration* will appear



- Click **Register with RHN**

4. The window *Red Hat Network Configuration* will appear




- Click on the tab pane **Retrieval / Installation**

5. The window changes to *Retrieval / Installation*

- Disable the entry


*Do not upgrade packages when local configuration file has been modified*



- Click on the tab pane **Package Exceptions**


6. The window changes to *Package Exceptions*

- Click on the entry *kernel* and click on *Remove*




- Click **OK**

## 7. The window *Question* will appear




- Click **Yes**

## 8. The window *Welcome to Red Hat Update Agent* will appear




- Click **Forward**

## 9. The window will change to *Review the Red Hat Privacy Statement*



- Click **Forward**


## 10. The window will change to *Login*



Enable *Use Existing Account* and enter the username and password that you have specified in chapter 15.7.1


- Click **Forward**

## 11. The window will change to *Register a System Profile - Hardware*




- Click **Forward**

## 12. The window will change to *Register a System Profile - Packages*




- Click **Forward**

## 13. The window will change to *Send Profile Information to Red Hat Network*



- Click **Forward**

If you have already activated your license code as it is described in chapter 15.7.1, the process will go on with updating your system. Just follow the instructions given. If you have not done the activation so far, the following information window would appear.



### 15.7.3 Activate the Red Hat Alert Notification Tool

RED HAT ENTERPRISE LINUX WS 3 offers a convenient tool that checks automatically the available security fixes that can be downloaded. If there are any of those packages available or if the connection to Red Hat Network not working, he shows a blinking, red exclamation mark at the lower right side of the desktop.



If the system is up-to-date and the connection to the Red Hat Network is fine, it shows a friendly white hook in a green circle




To activate this tool proceed as follows:



1. Log in as root
2. Click on the Icon at the lower right side of the desktop




3. The window *Red Hat Network Alert Notification Tool* will appear




- Click **Forward**

4. The window *Terms of Service* will appear



- Click **Forward**


5. The window *Proxy Configuration* will appear



If you use a proxy server enter the required information, if not just leave the fields empty

- Click **Forward**

6. The window *Configuration Complete* will appear



- Click **Apply**

## 15.8 The Emergency Recovery CD / Image CD

Each PC which is delivered from Bruker Germany has a license for the backup tool *Symantec Norton Ghost 2003*. With this software and the CD-Writer of the PC it is possible to create a one-to-one backup of the harddisk of the system. Because *Norton Ghost 2003* does not support GRUB boot loader, the PCs are still delivered with LILO.

It is recommended to periodically create a backup of the harddisk. In case of any problems (harddisk crash, virus, hacker attack ...) you are able to rebuild the situa-

tion of your last backup within a few minutes.

In order to create the Emergency Recovery DVD under Linux it is necessary to produce a bootable floppy disk, called Rescue Disk.

### 15.8.1 Creating a boot floppy

If you create an Emergency Recovery Disk, it is recommended to also create a boot floppy. For this follow the instructions given in chapter 13.1.1

### 15.8.2 Creating the Rescue Disk

1. Login as root (only root has the icon *Ghost2003* on the desktop)
2. Preposit one empty floppy disk
3. Click on the icon *Ghost2003* on the desktop

This will activate the script that creates the Rescue Disk. You will be asked for inserting the floppy disk.

### 15.8.3 Creating the Emergency Recovery CD

1. Insert the Rescue Disk into the PC
2. Reboot the PC
3. Nowadays it is not the default that a PC is booting first from floppy. In this case you might have to perform two steps to activate this feature.
  - change the boot order to boot first from floppy (during the boot-time of the PC you can choose to switch into this configuration)
  - change the BIOS settings to allow to boot from floppy (during the boot-time you can choose to switch into BIOS setup)

On the HPxw4100 currently delivered by Bruker Germany you can switch into the BIOS setup by pressing F10 (F9 would offer the boot sequence)

4. A window *About Norton Ghost* appears
  - click **OK**
5. The window *Norton Ghost 2003* appears
6. Click on **Local** → **Disk** → **To Image**
7. The window *Select local source ...* appears
  - Choose your harddisk (typically there is only one)

- Click **OK**
8. Insert a CD into the CD-Writer (number the CD's you use, starting with "1")
  9. The window *File name to copy image to* appears
    - Click on the arrow down under *Look in* and choose your CD-Writer
    - Click **Save**
  10. The window *Compress Image* appears
    - Choose **High**
  11. The window *Copy a bootable floppy to the CD/DVD disc?* appears
    - Click **Yes**
  12. You will be asked *Is the floppy disk ready in drive a:?*
    - Check whether the boot floppy has been inserted and click on **Yes**
  13. The window *Proceed with Drive Backup to CD/DVD?*
    - Click **Yes**
  14. Norton Ghost is going to write the first CD. Wait for completion, insert the next (numbered) CD when you will get asked for it and click OK, this is to be repeated until all data is stored on CD
  15. If the backup process is finished, a window *Image Creation Completed Successfully* appears
    - Click **Continue**
    - Click **Quit**
  16. You will be asked *Are you sure you want to quit?*
    - Click **Yes**
  17. Remove the CD and the Floppy, store them in a safe place (not too close to the PC!) and reboot the PC (by switching it off and on for example)

#### 15.8.4 Using the Emergency Recovery CD

This should be done by experienced people only! All data currently on the disk will be deleted!

1. Insert the first Image-CD into the CD-drive and reboot the PC
2. Ghost will automatically be started. If the PC doesn't boot from CD, you have to change the boot-order in the BIOS of the PC (see previous chapter)

3. Select **Local** → **Disk** → **From Image**
  - choose DVD
  - click **Open** → **OK**
4. The window *Destination Drive Details* appears. If you want to keep the original partition sizes, change the numbers under *New Size* to the same values as under *Old Size*.
5. If you want to change the partition sizes in order to use the whole disk (no unpartitioned space left), then you can change the values under *New Size* according to your needs.
6. To continue, follow the instructions.

---


## 15.9 Using the CD-Writer

---


In previous Linux versions there were the CD burning program ‚XCDroast‘. This is not part of a Red Hat Enterprise Linux WS 3. In this chapter you can find a description how to burn a CD under Red Hat Enterprise Linux WS 3 in a very easy and convenient way. A major difference to the behaviour of XCDroast is that you do not need the root password to burn a CD.

1. Open a LINUX shell and type in the command:  
**nautilus**
2. A browser window appears that shows the content of your home directory. Click on  
**Go** → **CD Creator**
3. The window changes to the Location

*burn:///*




4. Click on **File** → **New Window**
5. An additional nautilus window appears, starting again in your home directory
6. Go inside the second window to the information you like to burn on CD
7. Transfer the information what do you want to burn on CD from the second window into the first one by drag-and-drop
8. When you have transferred all files you want to burn into the first nautilus window, click on **Write to CD**
9. A window *Write files to a CD recorder* will appear, click on **Write files to CD**



10. The window *Writing CD* appears



11. If successfully written to CD, the window *Writing CD* changes to



12. Click **Close**
13. The CD has been burned, close all windows

## 15.10 Add a new harddisk



*NOTE: Using the command `fdisk` is recommended only for seasoned persons, because formatting the old harddisk will destroy all data on this disk!*

To buy a new harddisk on the Linux PC is quite easy because you can buy a usual harddisk in a PC store. After installing the new disk you have to do the following:

1. Start the PC
2. Open a shell and become root
3. If the new harddisk is a scsi-disk, type in the command

```
cat /proc/scsi/scsi
```

If it is a ide-disk use the command

```
cat /proc/ide/hdX/model
```

(the *X* in *hdX* have to be replaced with the correct number of the disk, which is *a*, *b*, *c* or *d*)

4. Does the command show the harddisk which was built in?
  - Yes? Go to step 5
  - No? There was something wrong during build in, or the harddisk is damaged.

5. Type in the command **fdisk -l**

the output looks like this :

| <i>Device</i>    | <i>Boot</i> | <i>Start</i> | <i>End</i> | <i>Blocks</i> | <i>Id</i> | <i>System</i> |
|------------------|-------------|--------------|------------|---------------|-----------|---------------|
| <i>/dev/sda1</i> | *           | 1            | 3          | 24066         | 4         | FAT16 <32M    |
| ...              |             |              |            |               |           |               |

A second harddisk would show up as device '*sdbx*', where *x* is a number. All Devices *sda<sub>x</sub>* are part of harddisk number 1, *sdb<sub>x</sub>* of number two and so on.

6. Does the output show the new disk?
  - Yes and the partitioning is OK? Go to step 9
  - Yes but the partitioning was not done so far? Go to step 7
  - No? Go to step 7
7. To partition the new disk, type in the command **fdisk /dev/sdb**  
In case you need help just type *m* and choose the possibility you want
8. After partitioning is done you have to create the filesystem on the new partitions. This can be done with the command  
**mkfs /dev/sdb<sub>x</sub>**  
where *x* is the partition you want to configure
9. Now you can add the new harddisk to the file */etc/fstab*, open this file with a text editor and add a line like this

```
/dev/sdb1 /X ext2 defaults 1 2
```

where */X* is an example for the mount point of the new harddisk. Be sure that the mount point you have chosen exists. If you want to use one that does not exist so far, you have to create this directory with the command

**mkdir *x***

in case you want to use */X*.

10. After reboot of the PC you can mount the new disk with the graphical tool  
**<red hat> → System → Disk Management**



---

# Chapter 16

## Appendix

---

### 16.1 Some notes for the helium handling

---

In the *TOPSPIN* home directory (default is `/opt/topspin`) in the directory `prog/logfiles` are two files named `heliumlog` and `heliumlog.err`. The `heliumlog` saves the current helium level in a new line every night at 3.00 am. The `heliumlog.err` documents possible errors during detection of the helium level.

#### 16.1.1 Automatic notification for helium refill

*TOPSPIN* is able to send an email notification to a list of defined users if the helium level is equal or less than a predefined limit.

##### 16.1.1.1 Set minimum limit of helium level

You can define the value of the helium level which should trigger the automatic notification (see next chapter)

1. Click **<red hat>** → **Editors** → **Text Editor**
2. A window *New Document - Text Editor* will appear
3. type in the percentage (without '%') of the helium level, that should activate the

notification, e. g. 25

4. Click **File** → **Save as**
5. Save the file as `heliumlimit` in the installation directory of `TOPSPIN` (default is `/opt/topspin`) in the subdirectory `conf/instr/<SPECT>`

### 16.1.1.2 Create shell script 'heliumwarn' for email notification

In this chapter you find an example for the shell script that will be execute the email notification in case the helium level is less than the entry (in %) in the file `heliumlimit`

1. Click **<red hat>** → **Editors** → **Text Editor**
2. A window *New Document - Text Editor* will appear
3. Copy the following example in the file

```
#!/bin/sh
send an email when the helium level is too low
#####
#-----
user defined variables:
#-----
ADDRESSES="person1@domain person2@domain person3@domain"
INSTRUMENT="spectrometer identification"
MAIL="/usr/bin/mutt"
#####
no more changes are needed below this point. You may tailor
the variable MSG for specific needs though. Be careful!
#####
#-----
catch command line arguments:
#-----
MYSELF=$0
LEVEL=$1
LIMIT=$2
#-----
get the name of the spectrometer host:
#-----
```

```

HOST=`hostname -f`
user=`whoami`
#-----
prepare the email message:
#-----
MSG="`n\
Caution:`n\
The magnet of the spectrometer $INSTRUMENT `n\
has an helium level of $LEVEL % `n\
which is less than the limit of $LIMIT %: `n\
Please order liquid helium and refill this magnet.`n\
`n\
This email is automatically send from $MYSELF (user $USER) `n\
on $HOST to:`n\
- $ADDRESSES`n\
"
#-----
mail the message:
#-----
echo -e $MSG/$MAIL -s "$INSTRUMENT : Helium refill" $ADDRESSES
exit 0

```

4. Modify the settings of "ADDRESSES" and "INSTRUMENT" and the text in the example as you like
5. Click **File** → **Save as**
6. Go in the installation directory of **TOPSPIN** (default is /opt/topspin) into the subdirectory prog/bin and save the file as heliumwarn

Now the defined email addresses will get a notification if the actual helium level is the same or below the entry (in %) in the file

```
<TOPSPIN_HOME>/conf/instr/<spect>/heliumlimit
```

The notifications are sent only once but the procedure will be activated again if the helium level is again higher than the limit.

Note that this script is only an example, if you would like to add more recipients to the notification list you had to add the respective addresses, if you like to have less addresses you have to delete the respective entries.

## 16.2 Script for saving important files

TOPSPIN offers a script (*xwinnmr.save*) for saving important files that are not delivered by BRUKER. This script can save and extract user specific data from the TOPSPINHOME directory and puts them into a tar-file. Files that start with an "." are ignored. The script can be started as follows:

- Open a Linux shell and enter

```
<TOPSPIN_HOME>/prog/bin/scripts/ xwinnmr.save
```

Files and directories will be saved into the special tar-file *xwin\_backup.tar* which is, by default, in the directory *TOPSPINHOME/xwin\_backup*.

The *xwinnmr.save* script saves the same files and directories as the AU program *nmr\_save* (see chapter 6.16.6). Nevertheless there are differences between *xwinnmr.save* and *nmr\_save*. This is because of compatibility reasons with previous XWIN-NMR installations and corresponding individual scripts that calling the *xwinnmr.save* script. The differences between these two ways to create a backup are opposed in the following:

| Script name         | default directory  | name of tar-file       |
|---------------------|--------------------|------------------------|
| <i>xwinnmr.save</i> | <i>xwin_backup</i> | <i>xwin_backup.tar</i> |
| <i>nmr_save</i>     | <i>nmr_backup</i>  | <i>nmr_backup.tar</i>  |

The options for *xwinnmr.save* starting from a command prompt are:  
***xwinnmr.save*** [-options]

- ***-help*** additional information
- ***-path <x>*** absolute path where backup tar-file *xwin\_backup.tar* will be created
- ***-restore <x> <y>*** absolute path where backup tar-file *xwin\_backup.tar* can be found (<x>) and absolute path where the restoring should take place (<y>) (installation directory of TOPSPIN)
- ***-replace*** overwrite the old backup.tar file
- ***-silent*** do not display any messages
- ***-source*** absolute path of TopSpin/XWIN-NMR installation that should be saved

- **-spect <x>** spectrometer name (e.g. spect; default is the current name)
- **-verbose** print more information while backing up

For safety reasons the tar-file will not be replaced by executing the **xwin-nmr.save** script once again. So if you like to use the script several times, you have to move the `xwin_backup.tar` file from the directory `TOPSPIN-HOME/xwin_backup` on a backup medium (floppy, CD, DVD ...) and then start the script once again or you have to use the 'replace' option.

---

## 16.3 Where Do I Find What?

---

If you cannot find the required information in this manual, please check the manuals listed below.

Further useful information you can also find in the Release Letter for the respective TOPSPIN version.

- *TopSpin Users Guide*  
Window handling, displaying data, interactive data manipulation, user interface related items
- *Acquisition Reference*  
Description of all acquisition commands, parameters and files.
- *Processing Reference*  
Description of all processing commands, parameters and files.
- *AU Reference*  
How to write AU programs. Description of all AU macros, functions and variables. List of Bruker AU programs.

For further manuals please refer to:

- TOPSPIN → **Help**
- <http://www.bruker-biospin.de/NMR/nmrsoftw/passwd/docu/index.html>
- NMR-GUIDE



---

# Chapter 17

## Books and Web addresses on Linux

---

### 1. Books on SAMBA

- Using Samba - Robert Eckstein, David Collier-Brown, Peter Kelly  
this book is online available as part of the SWAT (samba configuration tool)

### 2. Web pages on Linux

- <http://www.tldp.org/>
- <http://www.hp.com>
- <http://www.redhat.com>
- <http://www.bruker-biospin.de/NMR/nmrsoftw/passwd/faq/index.html>





---

# Chapter 18

## Bruker addresses

---

### Germany

Bruker BioSpin GmbH  
Silberstreifen  
D-76287 Rheinstetten  
Tel: (++49) (721) 51 61 0  
Fax: (++49) (721) 51 71 01

Bruker BioSpin Software Department  
Silberstreifen  
D-76287 Rheinstetten  
Tel: (++49) (721) 5161 455  
Fax: (++49) (721) 5161 943

<http://www.bruker-biospin.de>

ftp server: <ftp.bruker.de>

Email: [nmr@bruker-biospin.de](mailto:nmr@bruker-biospin.de) (sales)  
[mbu@bruker.de](mailto:mbu@bruker.de) (service)  
[applik@bruker.de](mailto:applik@bruker.de) (application)  
[nmr-software-support@bruker.de](mailto:nmr-software-support@bruker.de) (software)  
[license@bruker.de](mailto:license@bruker.de) (licenses)

## USA

Bruker BioSpin Corporation  
15 Fortune Drive  
Manning Road  
Billerica, MA. 01821-3991

Tel: (+1)(978) 667 9580 5195 (center)  
          (+1)(978) 667 9580 5444 (application)  
Fax: (+1)(978) 667 6168 (center)  
          (+1)(978) 667 2955 (application)

<http://www.bruker-biospin.com>  
ftp server: <ftp.bruker.com>

Email: [applab@bruker.com](mailto:applab@bruker.com)  
          [center@bruker.com](mailto:center@bruker.com)  
          [software@bruker.com](mailto:software@bruker.com)

## Switzerland

Bruker BioSpin AG  
Industriestraße 26  
CH-8117 Fällanden

Tel: (+41)(1) 8 25 91 11  
Fax: (+41)(1) 8 25 96 96

web server: [www.bruker.ch](http://www.bruker.ch)  
E-Mail: [epweb@bruker.ch](mailto:epweb@bruker.ch)  
          [all\\_ap@bruker.ch](mailto:all_ap@bruker.ch)

## France

Bruker BioSpin S.A.  
34, rue de l'industrie  
F-67166 Wissembourg/Cedex

Tel: (+33)(3) 88 73 68 00  
Fax: (+33)(3) 88 73 68 79

web server: [www.bruker.fr](http://www.bruker.fr)  
E-Mail: [support-rmn@bruker.fr](mailto:support-rmn@bruker.fr) (customer support)

## England

Bruker BioSpin LTD.

Banner lane  
Coventry CV4 9GH  
Tel: (+44) (2476) 855200  
Fax: (+44) (2476) 465317  
web server: [www.bruker.co.uk](http://www.bruker.co.uk)  
Email: [service@bruker.co.uk](mailto:service@bruker.co.uk)  
[apps@bruker.co.uk](mailto:apps@bruker.co.uk)

Our webpage

<http://www.bruker-biospin.de/NMR/about/offices/index.html>

provides contact addresses of our facilities and offices worldwide



---

# Index

---

## Symbols

- .bashrc 191
- .cshrc 192
- .profile 192, 197
- .rhosts 214, 215
- /etc/bootparams 162, 169
- /etc/exports 161, 163
- /etc/hosts 156, 158, 162, 168, 169
- /etc/hosts.allow 215
- /etc/init.d/bfsd 164
- /etc/nsswitch.conf 158
- /etc/rc2.d/S95bfsd 164
- /etc/shadow 187
- /etc/startd 156
- /etc/sysconfig/network 159, 163
- /etc/sysconfig/network-scripts/ifcfg-eth1 163
- /etc/xinetd.conf/bootparamd 164
- /sbin/arp 168
- /sbin/chkconfig 169, 213
- /sbin/reboot 165
- /usr/diskless/bfsd 159
- /usr/etc/nfsd 158
- /var/log/messages 174

## Numerics

- 2D 39
- 3D 39
- 4D 39

## A

- acquisition
  - mode 27, 67, 74
- Acquisition Reference Manual 19, 75
- Acrobat Reader 17, 98
- Adobe 17
- Amix 51
- Amix-Viewer 35
- arp 122, 162, 168

- ASP\_ST2 162, 164
- AURELIA 35, 39
- Aurelia 35, 51
- auto command 176
- automatic
  - configuration of a datastation 25, 75

## B

- BASH 79
- bash 190, 191, 193, 194, 198
- bfsd 159, 168
- bfsd daemon 159, 164
- bind nameserver 158
- Boot Floppy 177
- bootfile 166
- bootmode 166, 170
- bootparamd 98
- bootparams file 162, 169
- Bruker addresses 241
- BSMS 79

## C

- CCU 23, 35, 39, 42, 43, 125, 155, 165, 167, 168
- cd 178
- cf 22, 25, 42, 62, 155
- cf makelist 76, 77
- chgrp 181
- chkconfig 159, 169
- chmod 179, 198, 199, 214
- chown 180
- clipboard 88
- command
  - cf 25, 33, 75
  - configuration 33
  - expinstall 33
- compress 184
- configuration 59
  - datastation 25, 75
  - default 25, 75

- directory 71
  - name 27, 66, 73
- conffout 79
- conventions 10
- copy and paste 198
- CPU 11
- cron 200
- crontab 200
- csh 192, 194
- cu 166
- CUPS 127
- curinst 156

## D

- DAEMON line 54, 55
- datastation
  - configuration 75
- dd 164
- default
  - printer 27, 67, 74
  - selection in expinstall 27, 66, 73
- Demo license 52, 54, 55, 56, 57
- desktop manager 18
- Developer license 52
- DHCP 107, 108
- Diskless 35, 118, 162
- diskless directory 161, 162, 164
- display resolution 11
- dns 158
- domainname 167, 215
- dump file 163

## E

- echo 199
- edhead 27, 69
- edlock 28, 70
- edprosol 28, 70
- edscon 28, 71
- edsolv 27, 69
- emacs 192
- environment variable 192, 193, 197
- Esc 172
- eth0 109, 112, 115, 162
- eth1 163
- ethernet cable 162

- Ethernet card 11, 53, 54, 56
- ethernet interface 162
- examId\_1H 25, 62
- executable 198, 199
- expinstall 22, 27, 42, 64, 71, 75
- expinstall command 25
- export 87, 191, 197
- export variable 194
- exportfs 161
- exports file 161, 163

## F

- FAQ 170
- FCU 13
- FEATURE line 54, 55, 56, 151, 152
- find command 171
- Firewall 96, 210
- firmware 27, 64, 77
- FLEXLM 35
- Flexlm license manager 53, 55, 56, 57
- Floating license 53, 54, 55
- fstab file 163
- ftp 185, 212
- ftp server 45
- Full license 52, 56, 57
- full name 190

## G

- Ghost 226
- Gimp 215
- glibc-devel 100
- glibc-kernheaders 100
- GLP 35, 39, 51
- GNU TOOLS 35
- Good Laboratory Practice 35, 39
- graphics card 11
- graphics file
  - export 87
  - import 88
  - png 87
  - tif 87
  - wmf 87
- grep 182, 200
- group 179, 189, 190, 191, 193, 202
- group id 202

gs 155  
gzip 184

## **H**

Hardware requirements 11  
hardware tests 176  
hardware\_list 76  
headcrash 75  
helium 233  
helium refill 104, 233  
heliumlimit 234  
heliumlog 233  
heliumwarn 235  
help 171, 176  
history 194  
history function 171  
home directory 190  
host 191  
hostid 54, 56, 151  
hostname 108, 121, 122, 160, 167, 168, 215  
hostname -f 108, 160  
hosts file 156, 158, 162, 168, 169  
hosts.allow 214  
HP xw4100 12

## **I**

ICON-NMR 35, 39, 51  
ifconfig 122  
Image 226  
install.net 161  
installation directory 25, 75  
Installation Guide 9  
installnmr 34  
iptables 157, 212

## **J**

Java 47

## **K**

KDE 18  
kill 184  
known bugs 170  
known problems 170

## **L**

license 51  
license.dat 53, 54, 55, 57, 151, 152  
licensing program 53  
LINUX 100  
LINUX function 197  
LINUX shell 198  
lists 199  
lmgr 53  
lmreread 57  
lmutil 53  
logoff 152  
loop structure 199  
LPRng 127  
ls 179, 193, 199, 202

## **M**

Mac address 56  
MAXENT 51  
minihub 162  
mkbootdisk 177  
mkdev 165  
Monitor Mode 166, 170  
mount 125  
mouse control 172  
mutt 100, 121

## **N**

nameserver 114  
nautilus 229  
netaddr 166  
netmask 166  
netstat 122  
network 10, 31, 55, 108, 122  
NFS 125, 158  
NFS Server 125  
nfsd daemon 158, 168  
nfs-utils 98, 100, 101  
nis 158  
NMR Administration Password 26, 33, 41, 62  
NMR SuperUser 33, 41  
nmr\_save 60, 80  
NMR-GUIDE 35, 39, 46, 47, 51  
NMR-SIM 35, 39, 51, 55  
nmr-software-support 205

nmrsu 22, 33  
Node address 56  
Node Locked license 54  
Norton Ghost 226  
Notification Tool 224  
nsswitch.conf file 158, 168

## O

only-processing-PC 61  
openmotif 100  
order number 55  
owner 179, 180, 193, 202

## P

paper format 27, 67, 74  
PARAVISION 51  
passwd 188  
password 188, 189, 190, 191  
patchlevel 45  
pdf 9  
performance 11  
PERL 35  
permissions 192, 193, 201, 202, 203  
physical address 54, 56  
ping 122  
png 87  
Powerpoint 88  
pp 199  
pre-scan-delay 27, 67, 74  
Presenter license 52  
Processing-Only license 52  
ps 183, 200  
pulse program 199  
pwd 198

## R

rcp 185, 214, 215  
RCU 13, 14  
rcu 176  
rcutest 176  
RED HAT 17, 171  
Red Hat Documentation 95  
Red Hat Enterprise Linux WS 3.0 12, 95  
Red Hat Network 217  
redhat-config-nfs 100

Register 219  
Release Letter 9, 21, 37  
remote files 122  
reset button 166, 167, 170  
reviveccu 158, 159, 161, 162, 172, 173  
rlogin 214, 215  
rm 199  
route 122  
rs232 165  
rsh 214  
rsh-server 212

## S

Samba 122, 212  
samba-swat 212  
scratch installation 10  
screenshot 215  
select  
    default printer 67, 74  
    spectrometer configuration name 27, 66, 73  
    the default printer 27  
server line 55  
service command 159, 169  
setenv 166, 170, 194  
sh 194, 199  
shadow 187  
shell script 198, 199  
shmm 173  
showconf 79  
Software requirements 17  
spect 156, 162  
spectrometer 118  
spectrometer-PC 25, 61  
startd process 156  
su - 164  
swapfile 162, 164, 169  
SWAT 122  
swat 122  
systemst 176

## T

tar 186, 187  
TCP/IP 107, 111  
tch 192, 194, 198  
TCU 13



Teaching license 52  
telnet spect 155, 156, 157, 176  
telnet-server 212

tif 87

## TOPSPIN

- CDROM 9, 53, 56
- commands 35
- configuration 24, 31
- Express Installation 19
- installation 23, 31, 36, 42, 43
- installation directory 25, 75
- license 31, 55
- Online Help 9
- owner of the program files 33
- programs 35
- Release letter 9
- start 24

topspin -client 152

TOPSPIN Plot Editor 35, 51, 55, 88, 127

touser 173

traceroute 122

transceiver 162

## U

umask 192

UNIX host 125

up2date 219

User Manager 89, 91

uucp 100, 101, 165

uxproc 154

## V

Very Important Files 80

VIF 80

## W

Web Browser 17

who 171, 172, 173

Windows 2000 123, 125

wmf 87

wobb 155

Word 88

WWW account 171

## X

XCDroast 229

Xconfigurator 197

xinetd.d 164

## Z

zg 155

