



# **AVANCE III Cabinet**

## **Avance III Wiring User Manual**

**Version 001**

**Bruker BioSpin**

The information in this manual may be altered without notice.

BRUKER BIOSPIN accepts no responsibility for actions taken as a result of use of this manual. BRUKER BIOSPIN accepts no liability for any mistakes contained in the manual, leading to coincidental damage, whether during installation or operation of the instrument. Unauthorized reproduction of manual contents, without written permission from the publishers, or translation into another language, either in full or in part, is forbidden.

This manual was written by

Stephane Kreiss

© June 18, 2007: Bruker Biospin GmbH

Rheinstetten, Germany

P/N: Z31812

DWG-Nr.: Z4D10273001

# Contents

	<b>Contents .....</b>	<b>3</b>
<b>1</b>	<b>Introduction .....</b>	<b>5</b>
1.1	Introduction .....	5
1.2	Disclaimer .....	5
1.3	Safety Issues .....	5
1.4	Warnings and Notes .....	6
1.5	Contact for Additional Technical Assistance .....	6
1.6	Declaration of Conformity .....	7
<b>2</b>	<b>Console Configuration .....</b>	<b>9</b>
2.1	AVANCE III MicroBay .....	9
2.2	AVANCE III OneBay .....	10
2.3	AVANCE III TowBay HR .....	11
2.4	Avance III TwoBay Solids .....	12
<b>3</b>	<b>Internal Wiring .....</b>	<b>13</b>
3.1	AVANCE III Wiring .....	13
<b>4</b>	<b>Main Power Wiring .....</b>	<b>25</b>
4.1	AVANCE III MicroBay & OneBay 230V 16A .....	25
4.2	AVANCE III TwoBay 230V 16A .....	26
4.3	AVANCE III TwoBay 230V 32A .....	27
4.4	AVANCE III TwoBay 380V 3x16A .....	28
4.5	AVANCE III OneBay Accesory Cabinet .....	29



# Introduction

# 1

---

## Introduction

1.1

The fundamentally superior precision and stability of fully digital signal generation and processing was introduced and established by the precedent-setting series of AVANCE™ NMR spectrometers. With its digital advantage, the Bruker AVANCE™ series set revolutionary standards for performance, long-term reliability and ease-of-use, whether for routine applications or the most demanding research.

The next-generation Avance III NMR spectrometer series features a second-generation digital receiver technology (2G-DR™), a new milestone in NMR detection fidelity. The 2G-DR provides significant benefits for the most demanding NMR applications, e.g., materials science, polymer analysis, trace analysis, LC-NMR, MR imaging and structural biology, particularly for measurements with Bruker Bio-Spin's ultrasensitive CryoProbes™.

This manual provides detailed information on the AVANCE III console configuration and wiring.

---

## Disclaimer

1.2

The AVANCE III consoles should only be used for their intended purpose as described in this manual. Use of the unit for any purpose other than that for which it is intended is taken only at the users own risk and invalidates any and all manufacturer warranties.

Service or maintenance work on the consoles must be carried out by qualified personnel.

Only those persons schooled in the operation of Bruker spectrometers should operate the unit.

Read this manual before operating the unit. Pay particular attention to any safety related information.

---

## Safety Issues

1.3

The spectrometer hardware is no more or less hazardous than any typical electronic or pneumatic hardware and should be treated accordingly. Do not remove any of the protective panels from the various units. They are fitted to protect you and should be opened by qualified service personnel only. The main panel at the rear of the console is designed to be removed using two quick release screws, but again, this should only be done by trained personnel. Please note that, unless disconnected, cooling fans on the rear panel will continue to run even with the panel removed.

There are two types of information notices used in this manual. These notices highlight important information or warn the user of a potentially dangerous situation. The following notices will have the same level of importance throughout this manual.



---

Note: Indicates important information or helpful hints

---



---

**WARNING:** Indicates the possibility of severe personal injury, loss of life or equipment damage if the instructions are not followed.

---

For further technical assistance, please do not hesitate to contact your nearest BRUKER dealer or contact us directly at:

BRUKER BioSpin GMBH  
am Silberstreifen  
D-76287 Rheinstetten  
Germany

Phone: + 49 721 5161 0  
FAX: + 49 721 5171 01  
Email: [nmr-software-support@bruker-biospin.de](mailto:nmr-software-support@bruker-biospin.de)  
Internet: [www.bruker-biospin.de](http://www.bruker-biospin.de)



## **DECLARATION OF CONFORMITY**

The undermentioned product

### **NMR Spectrometer AVANCE III Cabinet H03128T2**

conforms to the main requirements  
set by the commission for the  
Harmonization of Regulations of the EU Member States  
with regards to electromagnetic compatibility  
(EMI 89/336/EWG) and safety (Low Voltage Electrical  
Equipment: 72/23/EWG) regulations.

For the assessment the following norms were applied:

EMI: EN 61326:97+A1:98+A2:01+A3:03  
Safety: EN 61 010-1

Tested by: Nemko GmbH

Test report: FS-0603-63703

Manufacturer's Name: Bruker Elektronik D-76287 Rheinstetten  
Bruker S.A. F-67166 Wissembourg  
Bruker AG CH-8117 Fällanden

Declaration approved by:

Dr. Tonio Gianotti	Jean Yves Fraval	Werner Schittenhelm
Head of Development	Technical Manager	Direction

Rheinstetten: May 09, 2006







# Console Configuration


# 2


AVANCE III MicroBay

2.1









# *Internal Wiring*


# 3

*AVANCE III Wiring*

3.1

---

Figure 3.1. Avance III MicroBay 2 Channel with AQS Preamp Page 1/2



WIRING DIAGRAM


Part-No. for Drawing only  
AVANCE III

Drawn  
KST  
18.06.07

Modification



Figure 3.2. Avance III MicroBay 2 Channel with AQS Preamp Page 2/2




WIRING DIAGRAM

Part-No. for Drawing only  
AVANCE III

Drawn  
KST  
18.06.07

Modification



WIRING DIAGRAM

Part-No. for Drawing only  
AVANCE III


Drawn  
KST  
18.06.07

Modification





Figure 3.4. Avance III MicroBay 3 Channel with HPPR/2 Page 2/2



WIRING DIAGRAM


Part-No. for Drawing only  
AVANCE III

Drawn  
KST  
18.06.07

Modification



Figure 3.5. Avance III OneBay 3 Channel HR Wiring page 1/2



WIRING DIAGRAM


Part-No. for Drawing only  
AVANCE III

Drawn  
KST  
18.06.07

Modification



Figure 3.6. Avance III OneBay 3 Channel HR Wiring page 2/2



WIRING DIAGRAM


Part-No. for Drawing only  
AVANCE III

Drawn  
KST  
18.06.07

Modification



Figure 3.7. Avance III TwoBay 3 Channel HR Wiring page 1/2



WIRING DIAGRAM


Part-No. for Drawing only  
AVANCE III

Drawn  
KST  
18.06.07

Modification



Figure 3.8. Avance III TwoBay 3 Channel HR Wiring page 2/2



WIRING DIAGRAM


Part-No. for Drawing only  
AVANCE III

Drawn  
KST  
18.06.07

Modification



Figure 3.9. Avance III TwoBay 3 Channel Solids Wiring page 2/1



WIRING DIAGRAM


Part-No. for Drawing only  
AVANCE III

Drawn  
KST  
18.06.07

Modification



Figure 3.10. Avance III TwoBay 3 Channel Solids Wiring page 2/2



Part-No. for  
Drawing only  
AVANCE III

Drawn  
KST  
18.06.07

Modification

WIRING DIAGRAM








# Main Power Wiring

# 4

## AVANCE III MicroBay & OneBay 230V 16A


## 4.1





# Main Power Wiring

## AVANCE III TwoBay 230V 16A

4.2








# Main Power Wiring

## AVANCE III TwoBay 380V 3x16A

4.4





# Main Power Wiring



## Notes: