

CryoProbe Prodigy System

- Site Planning & Checklist

Version 001



Copyright © by Bruker Corporation

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form, or by any means without the prior consent of the publisher. Product names used are trademarks or registered trademarks of their respective holders.

© October 19, 2012: Bruker Corporation

Fällanden, Switzerland

P/N: Z33080

DWG-Nr.: Z4D12020

For further technical assistance on the CryoProbe Prodigy System, please do not hesitate to contact your nearest BRUKER dealer or contact us directly at:

BRUKER Biospin AG
Industriestrasse 26
CH-8117 Fällanden
Switzerland

Phone: + 41 825 97 97
FAX: + 41 825 94 04
E-mail: cryoprobe.service@bruker.ch
Internet: www.bruker.com


Prodigy System: Site Planning

The BRUKER CryoProbe™ Prodigy System (CPP-System) should be positioned in a lab such that it allows for convenient and safe operation. The customer needs access to the front of the Prodigy unit for operating the touch-screen, and to the LN2 dewar for periodic LN2 refills.

Example A (a suitable N₂ gas supply is available in the lab):

The Prodigy Platform components are shown in **blue**.

- Ideally, place the Prodigy unit to the right side, and the LN2 dewar to the left side of the magnet. A mirrored setup is also possible. This arrangement provides the user with simplicity and functionality when performing the regular LN2 dewar refill.
- The Prodigy unit must be placed outside the 5mT (50G) stray field of the magnet, but remain in its vicinity.
- The LN2 dewar is non-magnetic and will be placed directly adjacent to the magnet. There is only one length of the LN2 Transferline available (2m).
- An optional balance for the LN2 dewar is recommended (BH5450; with shielded magnets only). The balance is placed under the LN2 dewar.



The N₂ gas supply system is marked in **green**.

- The Prodigy unit needs a continuous supply of very dry N₂ gas (dew point <-60°C).
- The N₂ gas cylinder serves as backup for the aforementioned main N₂ gas source. It is optional, but strongly recommended when a sample safety system (BH0550 or BH0551) will be used.

The necessary safety infrastructure is marked in **orange**.

- The ventilation and air conditioning (A/C) system must provide fresh air. The use of totally recirculated air is not allowed.
- An additional supply of fresh air can be provided e.g. by an open window. Alternatively, the A/C can be switched to maximum fresh air supply while refilling the LN2 dewar.
- An oxygen monitoring and alarm system is required.
- The door of the lab must open outwards from the room.

Electricity supply:


- A single phase supply 100-120VAC or 220-240VAC is sufficient for the Prodigy unit and the optional balance. Depending on customer's needs a UPS is recommended.

The system requirements are summarised on p. 9.

Example B (if no suitable N₂ gas supply is available):

Example B corresponds to example A except for the N₂ gas supply (see **green** items in the adjacent figure).

- The required N₂ gas is generated from a compressed air supply with a nitrogen separator.
- If the dew point of the compressed air is >-25°C, an air dryer is required in addition.



Checklist for Basic Site Planning

| Lab space requirements | YES | NO | |
|------------------------|--|--------------------------|--------------------------|
| | <p>The free volume of the lab is at least 70 m³.</p> | <input type="checkbox"/> | <input type="checkbox"/> |
| | <p>It is possible to arrange the Prodigy Unit beside the magnet, but outside the 5mT (50G) stray field.</p> | <input type="checkbox"/> | <input type="checkbox"/> |
| | <p>It is possible to arrange the LN2 dewar beside the magnet.</p> | <input type="checkbox"/> | <input type="checkbox"/> |
| | <p>The LN2 dewar is freely accessible, and the site is adequately sized for the free movement of the LN2 transport vessel.</p> | <input type="checkbox"/> | <input type="checkbox"/> |

| N2 gas requirements (green items) | YES | NO | Customer Responsibility |
|--|--------------------------|--------------------------|--------------------------|
| <p>Case 1: An uninterrupted supply of nitrogen gas as specified will be available (at least from installation date onwards) at the NMR site (connector type Ø 8mm).</p> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| <p>Case 2: Compressed air with a dew point < -25 °C will be used → N₂-separator required.</p> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| <p>Case 3: Compressed air with a dew point > -25 °C will be used → Air dryer and N₂-separator required.</p> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

| Electrical requirements | YES | NO | Customer Responsibility |
|---|--------------------------|--------------------------|--------------------------|
| <p>Electric power for the Prodigy unit is available at an electrical power outlet. 100-120VAC/max 10A or 220-240VAC/max 5A, 50-60Hz, 0.6kW max. Length of power cable: 7.6m / 25ft (with adapter for Schuko and US Plug).</p> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| <p>In case of voltage fluctuations which exceed +/-10% of the nominal voltage value: Is there a voltage stabiliser available for the Prodigy system?</p> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| <p>Will an uninterrupted power supply (UPS) be available to connect the Prodigy system as well?</p> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |


| Electrical requirements | YES | NO | Customer Responsibility |
|--|--------------------------|--------------------------|--------------------------|
| Optional balance for the LN2 dewar: Is an additional electrical power socket for the electronic balance available? 100-120VAC or 220-240VAC / 50-60Hz Length of power cable: 7.6m / 25ft (with adapter for Schuko and US Plug). | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

| Safety (orange items) | YES | NO | Customer Responsibility |
|---|--------------------------|--------------------------|--------------------------|
| The exit doors of the magnet room open outwards. | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| An oxygen monitoring and alarm system is installed in the lab. | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| A fresh air supply system is available to provide adequate ventilation: normal operation @ fresh air supply rate $\geq 200\text{m}^3\text{h}^{-1}$ refill situation @ fresh air supply rate $\geq 400\text{m}^3\text{h}^{-1}$ (or open windows) | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| The LN2 dewar will be secured against tipping over (earthquake risk; similar to magnet). | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

| In case of an existing console: | YES | NO |
|---|--------------------------|--------------------------|
| For 400 MHz magnets: Is a BOSS 1, a BOSS 2 or BOSS 3 Shim system of PLUG type available? (if not, an upgrade has to be performed) | <input type="checkbox"/> | <input type="checkbox"/> |
| For 500 MHz magnets and above: Is a BOSS 3 Shim system of PLUG type available? if not, an upgrade has to be performed | <input type="checkbox"/> | <input type="checkbox"/> |
| Is a BVT3000 or a BVT3200 or a BVT3200A or a BSVT unit available? | <input type="checkbox"/> | <input type="checkbox"/> |
| Is a GREAT 1/10 or 3/10 or a GAB respectively a GAB/2 available? | <input type="checkbox"/> | <input type="checkbox"/> |
| Is a Nanobay400 'V3' (Z133003) with AQS preamplifier available? (*type label on the back plane) If yes, the system is ready for a CPP BBO. Notice: 19F capability on 1H is supported via manual rewiring at the probe for CPP BBO. | <input type="checkbox"/> | <input type="checkbox"/> |
| Is one of the following consoles with AQS preamplifier available? <ul style="list-style-type: none"> • Nanobay400 'V2' (Z119572) (*type label on the back plane) • Nanobay400 'V1' (Z108356) (*type label on the back plane) • Microbay400 (H03128MB) (*type label inside cabinet, top right corner) If yes, the console has to be upgraded. Notice: 19F capability on 1H is supported via manual rewiring at the probe for CPP BBO. | <input type="checkbox"/> | <input type="checkbox"/> |

| In case of an existing console: | YES | NO |
|--|--------------------------|--------------------------|
| Is a console (AV II or newer) with HPPR/2 preamplifier available? If yes, an upgrade to CryoProbe-compatible HPPR/2 modules can be necessary. For CPPs with H & F feature an upgrade to HPLNA 1H is recommended. | <input type="checkbox"/> | <input type="checkbox"/> |
| In case of an existing console, is Topspin ≥3.1 PL 4 available? if no, the Topspin software has to be upgraded | <input type="checkbox"/> | <input type="checkbox"/> |
| Take a copy or print out of the configuration information of the existing system: <topspin-home>/conf/instr/<instrument name>/uxnmr.info | <input type="checkbox"/> | |

*) Example of type label:



| Is a detailed Site Planning necessary? | YES | NO |
|---|--------------------------|--------------------------|
| A magnet from BRUKER, but not one of the following types: <ul style="list-style-type: none"> • 400 MHz: US+ or Ascend • ≥ 500 MHz: US, US+ or Ascend | <input type="checkbox"/> | <input type="checkbox"/> |
| A magnet from another company | <input type="checkbox"/> | <input type="checkbox"/> |
| A widebore magnet | <input type="checkbox"/> | <input type="checkbox"/> |
| The magnet is placed in a pit | <input type="checkbox"/> | <input type="checkbox"/> |
| The magnet is equipped with a sample changer system | <input type="checkbox"/> | <input type="checkbox"/> |
| Any other unclear situation | <input type="checkbox"/> | <input type="checkbox"/> |

If any of the questions in this section has been answered with YES, then contact cryoprobe.service@bruker.ch to clarify if an upgrade with a Prodigy system is possible.

Responsible partner at customer's site

| | |
|------------------------------------|--|
| Name | |
| Telephone | |
| Email | |
| Company / Institution | |
| Address | |
| Ordering No. of the Prodigy System | |
| Ordering No. of the spectrometer | |
| Date and Signature | |

Prodigy System: Requirements

| Equipment dimensions | |
|--|---|
| Prodigy Unit | 60 x 40 x 46 cm ³ (L x W x H); add 15cm for rear connections |
| LN2 dewar | 50 x 62 x 135 cm ³ (L x W x H; W includes dewar handles) |
| Electrical requirements | |
| Prodigy Unit | 100-120 VAC / 50-60 Hz / max 10 A / max 0.6 kW 220-240 VAC / 50-60 Hz / max 5 A / max 0.6 kW |
| Optional weight scale | 100-120 VAC or 220-240 VAC / 50-60 Hz |
| N ₂ gas requirements | |
| Case 1: In-house N ₂ gas supply used as main N ₂ source | <p>Pressurised N₂ Gas:</p> <ul style="list-style-type: none"> - N₂ content > 95% by volume - Dew point < -60°C @ 1 bar - Pressure 6 - 10 bar - Oil content: < 0.005 ppm (0.00425 mg/m³) - Solid impurities: Use 5 micron filters. Filters should retain 99.99% of the specified particles <p>Capacity:</p> <ul style="list-style-type: none"> - Prodigy Unit: 40 l/min (1.4 cfm) + VT gas: 11 l/min (0.4 cfm) - Sample protection option: Add 40 l/min (1.4 cfm) to above values |
| Case 2: Compressed air with dew point <-25 °C available | <p>Use an N₂-gas separator to generate on site the required N₂ gas (as specified in Case 1). <i>Please contact Bruker, if you need help choosing a suitable N2 gas separator.</i></p> <p>Requirements for compressed air:</p> <ul style="list-style-type: none"> - Pressure: 7-10 bar - Oil content: < 0.005 ppm (0.00425 mg/m³) - Solid impurities: Use 5 micron filters. Filters should retain 99.99% of the specified particles. - Capacity: >200 l/min (7 cfm) |
| Case 3: Compressed air with dew point >-25 °C available | <p>Use an additional Air dryer in order to obtain compressed air with a dew point < -25°C @ 1 bar (as specified in Case 2 above). <i>Please contact Bruker, if you need help choosing a suitable Air dryer.</i></p> |
| Backup N₂ gas supply (optional; in addition to any of the main gas supply variants above) | <p>Independent N₂ gas source, e.g. a N₂ gas cylinder fitted with a pressure reducing valve. Specifications: See Case 1 above.</p> |
| Ventilation requirements: | |
| Free room volume (exchangeable gas volume) $V_{\text{free}} = V_{\text{room}} - V_{\text{furniture}}$ | ≥ 70 m ³ |
| Fresh air supply rate | ≥ 200 m ³ h ⁻¹ During LN2 refill: ≥ 400 m ³ h ⁻¹ |
| Further local regulations may apply. | |

Safety in the Lab:

WARNING

Risk of suffocation.

Risk of injury due to very low temperature liquids & metal parts.

Contact with the skin may cause cold burns.

Contact with the eyes may cause blindness.

Therefore:

- The lab must have a free room volume of $\geq 70 \text{ m}^3$ with a steady fresh air supply rate of $\geq 200 \text{ m}^3 \text{ h}^{-1}$ for each CryoProbe Prodigy System.
- Before the refill process is started, the fresh air supply rate must be increased to $\geq 400 \text{ m}^3 \text{ h}^{-1}$.
- Windows and doors must be opened before starting the LN2 refill.
- The lab must be equipped with oxygen monitors to detect a possible drop in the oxygen level. One oxygen monitor must be above the magnet and one oxygen monitor approx. 30 cm off the floor of the magnet room in order to detect low oxygen levels. Where the magnet is located inside a pit, an additional oxygen monitor located approx. 30 cm from the bottom of the pit must be fitted. All detectors should be located outside the 0.5 mT (5 G) line.
- Persons must not accompany a liquid nitrogen transport vessel inside an elevator (observe local regulations).
- The fill or refill procedure of the LN2 dewar must be carried out by trained laboratory personnel or trained personnel from a nitrogen supply company.
- During the entire refill process protective gloves, goggles, apron and personal oxygen monitor must be worn.
- Never look directly into the openings of components without eye protection (e.g. the Transferline) because liquid nitrogen droplets may spill out.
- The transport vessel for dispensing liquid nitrogen must be equipped with a safety pressure release valve, be non-ferromagnetic and must be placed outside the 0.5 mT (5 G) line.
- If the LN2 dewar is placed on the scale (optional equipment): The scale must be positioned outside the 0.5 mT (5 G) line of a shielded magnet.



Contact

NMR Hotlines:

Bruker Corporation provides dedicated hotlines and service centers. Please select the NMR service center or hotline you wish to contact from our list available at:

http://www.bruker-biospin.com/hotlines_nmr.html

Contact our NMR service centers, so that our specialists can respond as quickly as possible to all your service requests, application questions, software or technical needs.

CryoProbe Prodigy Sales:

BRUKER BioSpin AG
Industriestrasse 26
CH-8117 Fällanden
Switzerland
phone: +41 44 825 91 11
fax: + 41 44 825 96 96
email: sales@bruker-biospin.ch
<http://www.bruker.ch>

BRUKER Instruments, Inc.
44 Manning Road
Billerica, MA 01821
U.S.A.
phone: +1-978-667-9580
fax: +1-978-667-0985
email: sales@nmr.bruker.com
<http://www.bruker.com>

CryoProbe Services:

BRUKER BioSpin AG
Service Department
Industriestrasse 26
CH-8117 Fällanden
Switzerland
phone: +41 44 825 97 97
fax: + 41 44 825 94 04
email: cryoprobe.service@bruker.ch
<http://www.bruker.ch>

BRUKER Instruments, Inc.
44 Manning Road
Billerica, MA 01821
U.S.A.
phone: +1-978-667-9580
fax: +1-978-667-0985
email: center@nmr.bruker.com
<http://www.bruker.com>

Your Local Contact:

http://www.bruker.ch/contact_bruker.html



Bruker Corporation

info@bruker.com
www.bruker.com