

Nitrogen Fill Movie

This document is recorded as a QuickTime Movie File. It shows how to refill a magnet's liquid nitrogen vessel and gives safety advices.

WARNING:

When this film was recorded some time ago, latex transfer hoses were in use. It has been reported that latex hoses may split during use, therefore it is **IMPERATIVE** that you **do not use latex** as the material for the transfer hose. Hoses with metal sleeving, teflon hoses or surrogated PFA hoses are acceptable for use.

SGI systems

To view the movie on an SGI system, type

```
movieplayer /CDROM/misc/n2fills.mov
```

in a UNIX shell.

WINDOWS

To start the movie under WINDOWS, just click on this line or on the picture below. Install QuickTime from www.apple.com/quicktime/download/ if required.

This movie was prepared by BRUKER INSTRUMENTS INC. Manning Park, Billerica, M.A. All rights are reserved.

Document Part Number: **Z31338**

Drawing Number: **1047001**